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LIGHT BOX AND OPTICAL SET

REF. QLG003

INDEX OF LANGUAGE- ÍNDICE DE IDIOMAS

ESPAÑOL.....	2
ENGLISH.....	42

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LIGHT BOX AND OPTICAL SET

LIGHT BOX AND OPTICAL SET	42
ITEMS INCLUDED	44
INTRODUCTION	44
SETTING UP THE LIGHT BOX:	45
Important points of consideration while performing experiments	45
Possible problems encountered while performing experiments and their remedies:.....	46
REFLECTION	46
Experiment 1: reflection- single ray	47
Experiment 2: reflection-multiple rays- convergent or divergent rays.	47
Experiment 3: reflection- lateral and vertical inversion.....	48
Experiment 5: reflection-multiple reflections, single mirror	50
Experiment 6: reflection- multiple reflections, multiple mirrors.....	51
Experiment 7: reflection- rotation of a plane mirror	52
Experiment 8: reflection- images formed in a plane mirror and parallax	53
Experiment 10: reflection- curved mirror surface- center/radius of curvature and focal length of a circular mirror.....	54
Experiment 11: reflection- concave mirror- circular aberration and caustic curve.....	55

Experiment 12: reflection-parabolic reflector	57
REFRACTION	60
Experiment 14: refraction- semi circular block (measurement of refractive index of the material of semi circular block)	60
Experiment 16: refraction- rectangular block (lateral displacement and apparent depth)	64
Experiment 17: refraction- total internal refraction, semi circular block	65
Experiment 18: refraction- total internal refraction, triangular prism ..	66
Experiment 19: refraction- double refraction, angle of minimum deviation.	68
Experiment 21: refraction- double refraction using colour dispersion .	70
Experiment 22: refraction- colour absorption by filters	70
Experiment 23: refraction- newton's experiment with spectra	71
Experiment 24: refraction- double refraction using bi-convex lens	72
Experiment 25: refraction- radius of curvature	73
Experiment 26: refraction- focal line or focal plane of bi-convex lens ..	73
Experiment 27: refraction- spherical aberration in bi-convex lens	74
Experiment 28: refraction-chromatic aberration in bi-convex lens.....	74
EXPERIMENTS WITH COLOURS	76
Experiment 30: color observations- colors of objects	76
Experiment 31: colour observations-addition of colors	77

Experiment 32: color observations- color shadows- umbra/penumbra77
Experiment 33: colour ofservations-colour shadow 78

ITEMS INCLUDED

- Light box
- Set of 8 color cards
- Set of 8 color filters
- 2 slit former plates
- 1 plane mirror on stand
- 1 semi-circular mirror- Double sided
- 1 parabolic mirror- double sided
- 1 rectangular slab
- 1 semi-circular slab
- 1 lens biconvex large
- 1 lens biconvex small
- 1 lens biconcave
- 1 prism 45 45 90
- 1 prism 60 60 60
- 1 prism 60 30 90
- Spare lamp

INTRODUCTION

The complete kit consists of a source of light rays and a set of various optical components that reflect and refract light and color blend light, and allows study of the phenomena of reflection and refraction and a variety of color experiments, simultaneously permitting measurements and observations, conveniently. The unit operates on 12V AC / DC power rated at 2 to 3 amperes.

The light source is located in a specially constructed light box, with the top view as shown in the

Figure below:


Figure 1

The light box comes completely assembled along with an adjustable collimating lens

(As shown in the Figure 1).

Study of Colors: The source-end of the box is fitted with a light bulb and has one front opening and two side openings. The side openings are fitted with two hinged doors, each with a mirror, to reflect the light emanating from these openings. All the three openings have provision to allow the mounting of color filters (provided in the kit) at their front. With all the three windows having color filters in place, rotating the hinged side window mirrors, the light beam reflected by them can be swung back and forth to overlap and blend with the fixed center-beam from the front opening. A colorful pattern as a result of color mixing of the three beams can be easily observed on a screen placed in front of the center-opening at a distance of about 15- 20 cm from the light-box.

Study of Reflection & Refraction: The other end of the box also has an opening with a provision for mounting a slit plate or a color filter. Between the source and the opening is a collimating lens, the position of which can be changed with respect to the light source by means of a knob at the top of the light box to obtain parallel, slightly converging or diverging beam of light, as desired. This beam emerging from the opening can be split into one narrow beam (more suitable for production of spectra) or alternatively one, two, three or four narrow slit rays by mounting the appropriate slit former into the groove provided alongside the window opening. These rays or beams may be further colored by placing a color filter into the wider groove provided on the outer side of the window opening. The rays emerging from the opening are used to perform various experiments in reflection and refraction. The effect of slightly converging or diverging beam gets more pronounced for clear observation with the use of a multiple slit former. For greater degree of convergence or divergence, use the combination of lenses as shown below.


Figure 2

SETTING UP THE LIGHT BOX:

- Position the light box on a horizontal work bench surface with the connecting lead facing at the top. The light box has three non-skid feet at the bottom that ensures proper level and stability during operation.
- Connect the lamp to a low voltage (12V) AC or DC power source through the captive connecting lead provided along with, terminating in a

pair of 4mm banana plugs. The power source should have rating of at least 3A. operating the light box at lower voltages decreases the intensity of light, whereas operating it at higher voltages than the rated ones significantly reduces the life of the lamp and hence is not recommended.

- Insert one of the multiple slit former into the groove provided in front of the window and adjust the position of collimating lens by sliding the knob at the top to obtain a set of parallel rays.
- Place a plain white sheet of paper in front of the windows, along the path of the rays and adjust its position and height in such a way that the path of rays is clearly visible on the paper.
- Place the accessory to be used in various positions, depending on the experiment as detailed later on, for performing the experiment.

Important points of consideration while performing experiments

- The results of an experiment can be observed best in semi or full darkness, on account of the clear and bright visibility of the rays in such condition.
- The optical devices (mirrors, lenses, prisms, slabs) should be handled by their finger-grips, so that no smears, smudges or scratches are inadvertently left on the active optical faces. The bases of all the refractive optical components are especially finished to provide a background for obscuring the path of light rays to cause the rays to reflect so the path of the rays through the optical components is clearly visible.
- For tracing the ray paths, the optical accessory being used should be placed on a plain sheet of paper on which its outline can be marked by

running a sharp pencil around its perimeter. Rays can be traced by marking two points on each ray (incident on the accessory and demanding from it) as far apart as possible. Make sure that the markings are done at the center of the ray, one close to the surface of the accessory and other as distant as possible. If the ray pattern is complicated, with the rays crossing each other, carefully number the dots representing each other so that they can be easily followed.

After removing the optical accessory, carefully join the respective numbered points by a straight line using a ruler to indicate the ray paths to and from the optical accessory. Thereafter the ray path, within the accessory, can be traced by joining the incident and emergent point by a straight line. The direction of propagation can be indicated by arrowheads on the ray paths. In case of any doubt in regard to the continuity of any ray or line, replace the accessory exactly in the same position (as indicated by the marked perimeter) and retrace the ray.

Possible problems encountered while performing experiments and their remedies:

- 1. Faint secondary rays emerging from the light box:** these are caused by reflections from the support wire of the filament and can be removed by rotating the bulb-holder on top by 180 so that the support-wire lies behind the filament.
- 2. Internally reflected rays inside the device:** faint, secondary reflected rays (in the case of refracting surfaces) are normal, as there is some light always reflected at such surfaces, and can be neglected.

3. This is more pronounced if the refracting optical component is placed very close to the light box, with the internally reflected rays becoming brighter and more evident inside the component.

4. When the optical components are placed very close to the light box, bright internally reflected rays may become more pronounced inside the component. These are on account of the rays entering the optical component through the top face and internally reflecting off the vertical faces. In addition, the rays passing over the top of the component may be evident some distance component. Both of these effects lead of erroneous observations and can be removed or minimized by.

- Blanking off the top of the slits to shorten the height of the rays (not very helpful).
- Moving the optical component further from the light-box.
- Placing the optical component along with the paper sheet at a level slightly higher than the surface on which the light-box rests.

REFLECTION

Reflection is a phenomenon where light falling on a surface bounces back, in accordance with the following laws:

1. The angle of reflection is equal to the angle of incidence.
2. The incident ray, normal ray and reflected ray, all lie in the same plane.

Besides these laws, there are other important phenomena associated with the reflection that we usually come across every day, i.e. lateral inversion, multiple reflections, periscope, shaving mirrors, rear-view mirrors etc.

Experiment 1: reflection- single ray

Adjust the collimating lens so that the resultant beam of light is parallel. Allow a single ray to pass out of the light-box (use the respective slit former plate) mark its position on the paper. Place the plane mirror halfway along this path, crossing it at an angle.

Figure 3

Mark the position of the front face of the mirror, the reflecting rear face of the mirror and that of the reflected ray. Remove the mirror and draw a normal (perpendicular) to the outline of the mirror at the point where the incident and reflected ray meet the outline. Measure the angles between the normal and incident ray and normal reflected rays to obtain the angles of incidence and reflection respectively. Check that within the limits of experimental errors, the first law of reflection is verified.

The angle of incidence and reflection are measured with respect to the normal because in later experiments, reflection of rays from the curved surfaces will be studied. Since it is almost impossible to measure the angles between the rays and the curved surface of the mirror, the possible way is to draw a normal to the curved surface and from this straight line, measure the angles of incidence and reflection.

Next, replace the slit plate by the multiple slit plate with three or four slits. As above, place the mirror in the path of the rays and measure the various angles of incidence and reflection.

Experiment 2: reflection-multiple rays- convergent or divergent rays

RAY	ANGLE OF INCIDENCE	ANGLE OF REFLECTION
A		
B		
C		

Figure 4

Place a 3-slit or 4-slit former plate in front of the window opening of the light box. Adjust the collimating lens to get a set of diverging rays along the sheet of paper and mark the ray paths. If higher degree of divergence is required, place a concave lens in the path of rays as shown in figure 4. Place a plane mirror so that the rays meet it at angles other than 90° . Mark the reflecting surface of the mirror and the path of reflected rays. Draw normal to the mirror surface at each point of reflection of the light rays. Measure the angle of incidence and reflection corresponding to each ray and tabulate the observations in the table below.

The same experiment can also be repeated with the converging rays using convex lens or combination of concave and convex lens.

On the basis of above observations, answer the following questions.

- Do these observations match with what you expected?
- Is the angle of incidence greater than, less than, or equal to the angle of reflection?
- You have come across one of the laws of reflection- what is it?
- What happens to the parallel, converging and diverging rays after they are reflected? Did the diverging rays, converging rays or parallel rays maintain their characteristic after reflection?

Experiment 3: reflection- lateral and vertical inversion

When you look at an object and then at its image in a mirror, the two are observed not to be identical.

What is different about them?

The above phenomenon is called lateral inversion.

Remove all slit-plates from the light-box and adjust the collimating lens for parallel beam of the light. Stand a pencil in the left or right corner of the opening and then look at the image in the mirror. Where does the pencil stand?

From the above experiment can you explain why lateral inversion occurs?

Set the light box to project two parallel rays of light as shown in figure 5. Hold the green color filter in front of the left hand beam as seen from the front end or slit end of the box and the yellow filter in front of the right

hand beam so that the light rays incident on the plane mirror are of different colors. Observe the imagen.


Figure 5

Face the mirror and look at the reflection of the rays

Is the incident ray from the left side of the light-box and the image on the left side the same color i.e. is the left reflected beam green or yellow?

Record the rays in color pencil showing the green and yellow beams.

What happens to an image on reflection in a plane mirror?

Look at yourself in a plane mirror.

Is the image you see of yourself in a mirror the same as the image your friends see of you.

Do you observe any vertical inversion i.e. if your face is reversed from left to right, why is it not reversed from top to bottom?.

Why does vertical inversion not occur?

If you turn your head sideways, so that it is horizontal and not vertical, your reflection will be reversed vertically. Try it and see.

What is meant by "LATERAL INVERSION REFLECTION2?"

Try another prediction. If you observe the image of the word MAXAYAXAM in the mirror, what do you expect to see?

What should the reflection of the following capital letter word look like if a plane mirror is placed vertically along the dotted line and the reflection is observed in the mirror from a position at the bottom of the page?

CARBON-DI-OXIDE

Write down the expected image before you actually try the experiment

Was our prediction correct?

Was this lateral inversion?

If you were told to hold the word CARBON-DI-OXIDE and observe its reflection in a mirror, how would you hold it?

Is this the way it was represented the mirror in the previous experiment?

If the word is written on transparent paper and presented to the mirror in the two manners described, how would it look to you, if you viewed the mirror through the paper? Try it and see.

Is the latest inversion due to the mirror, or due to the way the word is presented to the mirror?

Why do some of the letters show reversal, while others do not?

If both halves of an object or image are mirror images of each other about a central line or axis, they are called symmetrical about that axis.

Are the above bold words symmetrical? If so, about which axes?

Is your face symmetrical?

Place a large plane mirror so that it stands out vertically from your face, straight down the center of your nose, then look yourself in another mirror.

Do you look normal with a perfectly symmetrical mirror image type of your face?

Which half reflected do you like the best?

Move the mirror to the side of your nose.

Is having two noses an improvement?

Move the mirror to the other side of our nose.

Do you like yourself with no nose and two eyes close together?

Experiment 4: reflection- position of image in a plane mirror

When you look into the image of an object in a mirror, the reflected image of the object appears to be somewhat behind the mirror. On moving the object slightly backward or forward with respect to the mirror, the reflection of the object also moves farther or nearer respectively. This experiment will assist you in locating the position of the image formed in a plane mirror.

Project a set of converging rays across the sheet of paper and record their position and focal point. Use the lens combination as shown in figure 6 and move them relative to one another to adjust the focal length.


Figure 6

Place a plane mirror across the converging rays at an angle and record the path of the reflected rays. While looking in the mirror at the reflection of the converging rays, lift the mirror and observe the real converging rays. Repeat the same by removing and replacing the mirror vertically several times, noting the similarity between the real and reflected rays.

Would you say the point where the real rays meet is the reflection of the point where the reflected rays meet?

Mark these two points of convergence and the position of the reflecting surface of the mirror. Carefully draw a line joining the two points of convergence.

What angle does this line make with the mirror?

What is the distance of the real convergence point and the of the reflection convergence point from the mirror?

Repeat the above experiment on another paper for another set of converging rays with different focus and the plane mirror position closer to or further from the point of convergence. Locate the convergence point and the mirror position. Draw the lines and measure the angles as before. Mount a pin vertically exactly at the two convergent points on the paper, using a cardboard beneath the paper. This gives one pin in front of the mirror and another pin hidden behind the mirror.

Lift the mirror vertically until you can see the hidden pin, then replace it, repeating this several times.

Is the mirror vertically until you can see the hidden pin, then replace it, repeating this several times.

Is the hidden pin located at the position of the reflected image of the front pin?

Does any shift in your position of observation affects the location of the image position?

Leave the front pin and the mirror unaltered, try it and see.

Also, try and explain, why is the image of the reflected rays fainter than the actual rays observed by removing the mirror?

Experiment 5: reflection-multiple reflections, single mirror

This experiment requires knowledge of refraction of light.

Allow a single ray to fall on the plane mirror placed very close to the light-box at an angle to the light ray. Examine the reflected rays closely.

How many rays do you observe? Which one is the brightest and which one is the faintest?

Explain how the images occur?

HINT: use the theory of refraction and the fact that the reflecting surface of the mirror is at the back.

Look carefully down on the mirror from the above and then make an enlarged drawing to show how the three reflected rays occur

What happens to the reflected rays as the angle of incidence is.

Decreased from 45 to 0 ?

Zero?

Increased from 45 to 90?

Whether all the three reflections occur at different angles of incidence varying as described above?

Which reflected ray disappears and at which angle of incidence does it occur?

Why is this reflection the faintest of the three?

To study an angle of incidence of 90 place the mirror so that its side completely covers the slit and the mirror is parallel to the ray.

Is there something special in angle of incidence being 90?

What do you observe at the other side of the mirror?

Rotate the mirror keeping the slit covered by one and of the mirror at all times. What happens?

The above phenomenon is used in OPTICAL FIBERS where light follows a curved path within the fiber.

Shift the mirror forward and allow the ray to fall on the side of the mirror and observe the mirror from above.

How does the above phenomenon occur?

HINT: internal reflection.

Place the rectangular block in front of the mirror and repeat the experiment and observations.

Experiment 6: reflection- multiple reflections, multiple mirrors

Borrow a plane mirror from another light box set. Place the two mirrors at right angles to each other. Allow a single ray to fall on the one of the mirrors at an angle, at a point about 25-30mm from the corner where the two mirror meet. Observe the principal reflected ray (not the fainter secondary reflections) as it is reflected from both the mirrors. Record this position.

Figure 7

What do you notice about the directions of the original incident ray and the final emergent ray?

Hint: both the rays should be parallel to each other.

Does this result occur at whatever angle you send the ray into the right angled mirror?

Using geometry, explain why this should be so?

Look into the corner of the mirrors. What do you observe?

If, you shift your position, does the same thing occur?

If a large plane mirror is available, place it horizontally on the work bench and make the other two plane mirrors stand vertically on it so that all are mutually perpendicular to one-another. A set of reflectors like this was placed on the moon by the astronauts.

Move aside and up and down, looking into the triple corner from all directions.

What do you observe in the triple corner?

Aim a beam of light into the corner

What do you notice about the reflection?

Suggest, why the reflector placed on the moon was a triple right-angled reflector?

What type of light beam was aimed at it and with what result?


Since, scientist know this type of reflection occurs on earth, what have they gained by placing the reflector on the moon?

Carefully examine the reflectors at the back of a car or bicycle. What is the shape of the dimples in the glass or plastic?

What can you say about the nature of light?

Experiment 7: reflection- rotation of a plane mirror

Aim a single ray of light at the plane mirror and mark the incident ray, mirror position and the reflected ray. Let the mirror position at this instance be marked as M and that of reflected ray as R.

Rotate the mirror slightly around the incident point, so that the incident ray strikes the mirror at the same position, but at a slightly different angle. Mark the new mirror position as M, and the new reflected rays as R.

Using a protractor, measure the angle of rotation of the mirror M i.e. angle between M and M and the angle of rotation of reflected ray R, i.e. angle between R and R1. Repeat this for different angles of rotation of the mirror with respect to the incident ray. Ensure that the point of incidence remains fixed for all observation. Record the observation in the table below and compare the two angles for each set of measurement.

{1 No.	Angle of Rotation of Mirror (K'	Angle of Rotation of Reflected Ray (P'	Ratio $\left(\frac{\angle R}{\angle K}\right)$
1			
2			
3			
4			

Find out the correlation between the angle of rotation of mirror M and the angle of rotation of the reflected ray R.

This optical phenomenon finds immense practical application in various instruments for exaggerating or amplifying small movements such as meters, spot-galvanometers, sextant etc.

Experiment 8: reflection- images formed in a plane mirror and parallax

Place a cardboard sheet beneath the paper sheet and mount a pin vertically on top of it. Position the plane mirror about 40 mm behind the pin and observe the reflection of the pin in the mirror. The reflected image appears to be behind the mirror. In order to locate the position of the image of the pin, aim a single ray of light to hit the mirror and adjust it in such a way so as to make it reflect back to the pin.

Does the reflection of the ray seen in the mirror go to the reflection of the pin?

If look just above the pin and along the ray towards the mirror is the ray bent or is it straight?

Mark the mirror position and incident ray only. Now, move the light box to a different position, aiming the reflected ray at the pin as before. Again mark the incident ray.

Repeat this procedure, moving the light box to several different positions on both sides of the pin and mark the incident rays only. Remove the mirror and draw the various incident rays, continuing them until they meet.

Do they all meet at one point?

Now replace the mirror.

Do the lines you have drawn in front of and behind the mirror lead to where you see the reflection of the pin?

If you are in doubt, stand a tall pin or knitting needle vertically at the point, where the lines meet behind the mirror so that you can see the top of the pin or needle and the reflection of the pin at the same time.

Do they appear to be in line?

Move your head towards one side and then to the other.

Do the pin reflection and the needle remain in line? Is there any relative motion between the needle and pin reflection?

This is a phenomenon known as PARALLAX.

PARALLAX is the apparent sideways motion of a distant object, relative to a near object, in the same direction as that of the movement of the observer. Astronomers use this phenomenon of parallax to show which of the stars or heavenly bodies are nearer to the earth and which are distant ones. As the earth moves along its elliptical orbit path, distant stars or celestial objects appear to move, relative to the nearer ones, in the same directions as that the motion of the earth.

Join the pin position to its image position.

Does this line cut the mirror line at right angles?

Is the image position as far behind the mirror as the pin is in front of it?

Move the pin to another position, and predict where its image will be. Repeat the earlier part of this experiment to check if your prediction is correct.

Curved Concave Surface

Figure 9

Select the semi-circular curved mirror. Aim a set of parallel rays into the center of the concave surface inside curved surface of the mirror so that the rays are also parallel to the axis of symmetry of the mirror as shown in figure 9. Record the incident and the reflected ray and mark where the reflected rays meet. This point is called the focus of the mirror. The distance of focus from the center of the mirror is called focal length of the mirror.

If the local point appears blurred and broad, with too many rays overlapping through it, block the outer rays as they leave the light box and use only the central ones.

Experiment 10: reflection- curved mirror surface- center/radius of curvature and focal length of a circular mirror.

Set up the apparatus in previous experiment and on a plane white paper, trace the inside concave reflecting surface of the mirror. Carefully, move the mirror around the curve and continue tracing until you have a complete circle. Measure the diameter of this circle along different axis or directions and calculate its average diameter.

What is the radius of the circle?

Find the center of the circle i.e. the center of curvature of the curved surface of the mirror.

How does this radius compare with the focal length of the mirror determined in the previous experiment?

Another method of finding the center of curvature of the mirror is to aim a single ray at the side concave curved surface of the mirror so that it reflects straight back on itself. This happens if the ray meets the reflecting surface along its normal i.e. the ray must be perpendicular to the reflecting surface at the point of incidence and must be reflected back along its radius position through the center of curvature. Record this ray position and without moving the mirror, carefully move the light box to another position, where the ray again reflects back on itself. Record this new ray position as well. Repeat this procedure for the third time


Figure 10

The point where the rays meet is the center of curvature of the curved mirror and the distance from this point to the curved reflecting surface of the mirror is the radius of curvature.

The center of curvature lies at a distance of twice the focal length F from the reflecting surface of the mirror. The distance of $2F$ is also called the radius of curvature R of the mirror.

Place a concave lens of short focal length in front of the slits to produce diverging rays. Shift the mirror until the reflected rays retrace their paths. The center of curvature is the point where the incident rays appear to diverge from and the radius of curvature is the distance of the point from the mirror.

If in the above experiment, you are unable to obtain a position where the reflected rays retrace themselves, use

$$\frac{1}{F} = \frac{1}{u} + \frac{1}{v} = \frac{2}{R}$$

Where

U = distance between the mirror and the point from where incident rays appear to diverge

V = distance between the mirror and the point where reflected rays meet.

(Here, all distances are measured along the axis of symmetry)

Experiment 11: reflection- concave mirror- circular aberration and caustic curve

What happens when parallel rays are directed at the circular mirror, parallel to the central axis of symmetry but displaced sideways away from the axis of symmetry?

CIRCULAR ABERRATION is the phenomenon in which parallel rays far away from the center axis of symmetry of the mirror do not focus onto the same point as those near the center producing a distorted image.


Figure 11


Figure 12

To study this aberration, allow two parallel rays of light from the light box to fall on the concave mirror as shown in figure 12 and mark the point where the reflected rays meet. Repeat this procedure by shifting the light-box to different positions so as to obtain at least six positions of the focal point, keeping the incident pair of rays parallel to the axis of symmetry of the concave mirror.

Do the various focal points obtained as above lie on a curve or on a straight line?

When four rays, which are parallel to the axis of symmetry of the mirror, are focused do they really focus at a single point?

Using a four slit former, block the inner two rays and record the position of the outer two rays and their focal point. Also mark the mirror position. Now, block the outer two rays and record the inner two rays and their focal point.

Which rays focus nearer to the mirror- the inner rays or the outer rays?

Take another sheet of plain paper and repeat the above procedure of inner and outer ray tracing, with the four rays parallel to the axis of symmetry but displaced sideways.

Is the circular aberration i.e. failure of all rays to focus at a single point increased or decreased by shifting the rays to one side away from the axis of symmetry?

Remove the four slit former from the light box and adjust the collimating lens to project a diverging beam. Hold the light box slightly above the table and remove it about 60-70cm away.

Project this wide beam into the mirror at a slight downward angle and observe the brightly illuminated area. Outline this area.

This shape is called a CAUSTIC CURVE, which can also be obtained by freehand joining the foci of the parallel rays as obtained above. Alternately, light from a bulb falling directly on the concave mirror also produces caustic curve caustic means burning.

Why is magnifying glass often called a burning glass?

What causes the shape of this caustic curve?

Slide the single slit former sideways across the face of the light box so that the wide beam is gradually cut off. Determine to where each part of the beam is reflected.

A shape like the caustic curve, which is the outline of the limiting positions of a movable or flexible line, is sometimes called an envelope. These envelopes are presently outside the scope of discussion and will be discussed in greater detail in higher mathematics along with the equations that generate them.

On one of the records of the reflections obtained previously, draw lines from the center of curvature to the points of reflection. Each of these new lines is a radius and is a normal to the mirror curve.

For each reflection, mark the angle of incidence (between the incident ray and the normal) and the angle of reflection (between the reflected ray and the normal)

Does each and of incidence equal its respective angle of reflection?

This problem of improper focusing, also called circular aberration can be minimized by using a paraboloid instead of spherical or parabola instead of circular mirror.

Experiment 12: reflection-parabolic reflector

A parabola is an unusually shaped curve that can be obtained from one of the following methods.

- Plotting a graph of $Y=X^2$ or some other quadratic algebraic function.
- Record the flight of a projectile through the air, freely gliding under the effect of gravity only.
- Move a point P in such a way that at any position, its distance from a fixed point F also called Focus is same as that from a straight line AB shown in figure 13 below. This line is called the directrix, because it directs or decides which way the parabola will face.


Figure 13

Cut a cone in a plane parallel to the sloping side. This is called marking a conic section. Other conic sections are- circles, ellipses or hyperbolas.

How are these made?

Aim a set of parallel rays into a parabolic reflector along the path parallel to the axis of symmetry of the reflector as shown in figure 14


Figure 15

Gradually, move the light box side ways keeping the rays parallel to the axis of symmetry of the mirror. You will find the graph paper a useful tool in aligning the mirror and ray paths to ensure that the incident rays are parallel to the axis.

What do you notice about the position of the focal point?

Aim a broad parallel sided beam of light into the parabolic mirror and observe the effect.

- What shaped mirror would be used to produce sharp images of stars scattered in all directions over the field of view?
- What would happen if a point source of light a torch globe is placed at a focal point of the parabolic mirror? Observe after performing experiment.
- Why do radar antennae, radio telescopes, car head-lamp reflectors have a parabolic shape rather than a spherical or circular shape?

- In the examples given, where in relation to the reflector, would the receiving, transmitting or radiating device be placed? Also suggest the reasons for the same. You will observe that in a parabolic mirror or reflector, distant as well as near rays focus at the same point as shown in figure 16.


Figure 16

Perform the same experiment as for obtaining the caustic curve. Do you obtain another curve or only a point?

Experiment 13: reflection- convex mirror


Project a number of parallel rays to strike the outside convex surface of the semi-circular mirror, parallel to the axis of symmetry of the later (as shown in figure 17)

Mark the mirror position and the ray paths, clearly indicating the ray directions with arrow heads.

Where do the diverging rays appear to come from?

This virtual point of origin of these diverging rays can be located by backward extrapolating the rays through the mirror along the reflected rays as shown in figure 18. This point where the rays seem to originate is called the virtual focus of the convex mirror and the distance of this point from the mirror is called the virtual focal length.

How does this focal length compare with the focal length of the concave side of the mirror?

How does it compare with the radius of curvature of the mirror determined previously?

What is the variation in these two results?


Can this difference or variation be on account of the mirror construction?

Suggest why slightly convex mirrors are used as rear view mirrors in automobiles.

If a line is drawn from the center of curvature through the point where a ray strikes the mirror, the line is normal to the surface. Draw a number of these normal corresponding to different points of incidence on the mirror. For each reflected ray, measure the angle of incidence (between the incident ray and the normal) and the angle of reflection (between the reflected ray and the normal)

For each reflection, does the angle of incidence equal its respective angle of reflection i.e. the basic law of reflection?

Place a convex lens of short focal length in front of the slits and adjust the convex mirror beyond it until the reflected rays coincide with the incident rays as shown in figure 19. Draw back these rays to the point where they appear to meet. This point is the center of the curvature of the mirror and the distance of this point from the mirror is the radius of curvature R.

In case co-incidence is not obtained above use the relationship

$$\frac{1}{F} = \frac{1}{u} + \frac{1}{v} \quad ? \quad \frac{2}{R}$$

Where F and R are already defined

U= distance between mirror and point where incident rays appear to meet and

V= distance between mirror and point where reflected rays appear to meet.

REFRACTION

Refraction is the change in the direction of propagation of light or bending of light at the boundary of two different transparent media having different densities, which occurs on account of the light traveling in different media with different speeds.

The basic law of refraction also called Snell's law is

$$\mu_2 \sin i = \mu_1 \sin r$$

Where μ = refractive index of the medium into which light refracts with respect to the medium the light is incident from

i= angle of incidence and

r= angle of refraction

Experiment 14: refraction- semi circular block (measurement of refractive index of the material of semi circular block)

Allow a single ray to fall normally (at an angle of 90°) at the center of the flat side of the semi-circular block as shown in figure 20.

Figure 20

Mark the position of the block and carefully observe the ray path.

Do you notice any deflection in the ray path?

In case, you notice any bending or deviation the light ray path, the angle at which ray is striking block surface is not 90°. Adjust the slab to make its incident surface perpendicular to the incident ray (so that the emergent ray is not deflected and is along the same line as that of incident ray). Mark this point of incidence on the incident surface of the block- this is the center of curvature of semi- circular block.

As shown in figure 21, shift the light-box gradually by rotating it about its center of curvature, so that the single ray strikes the same point of incidence on the flat vertical surface of the block at an angle other than 90

Clearly mark the ray path incident into the block as well as emerging out the block so that even on removing the block, the complete ray path through it can be clearly seen.

Repeat this procedure by rotating the light box further as in previous step for different incident angles (say 20, 30...) and clearly mark all the incident rays as well as emergent rays and the name them individually so as to have their clear identity. Make sure that the incident ray must enter the block at the same point of incidence i.e. center of curvature of the block, as marked earlier.

From the above experiment, what observation do you make?

What happens to the ray of light-

When it enters the semi circular block from the air medium at an angle other than 90 ?

When it emerges out of the semi circular block and enter the air medium again?

Try and explain why no bending of ray is observed as the ray cross the circular face of the block, while this is not the case when it crosses the flat face of the block?

Using the procedure explained above, record the observation for at least four different rays. For this, you can use rays at different angle of incidence on either side of the normal to the flat surface of the block. Carefully draw all the rays meeting at the mid-point of the flat side and their subsequent paths.

As shown in figure 22 draw a circle having diameter bigger than that of the semi circular slab say about 100mm with its center at the center of curvature of the block or point of incidence of the rays. Carefully extended the incident rays and the refracted rays so as to cur across the perimeter of the circle. Also draw the normal N to the flat incident surface, passing through the point of incidence.


Figure 22

Also, draw the perpendicular from the normal to the point of intersection of each of the incident and refracted ray with the perimeter of the circle shown as dotted lines in the figure 22. These lines are called half-chords.

Complete the following observation table

P_w	Angle of Incidence (i)	Angle of refraction (r)	Difference (i-r)	Ratio $\frac{i}{p}$	Length of Half Chord i (hc-i)	Length of Half Chord r (hc-r)
A						
B						
C						
D						
E						
F						

Ray	sin(i)	sin (r)	Ratio $\frac{\sin i}{\sin r}$	Ratio $\frac{hc-i}{hc/r}$
A				
B				
C				
D				
E				
F				

From the above table of observations, try and answer the following questions.

Is the difference between angle of incidence and angle of refraction (i-r) always remains approximately same?

Is the ratio of angle of incidence and angle of refraction (i/r) always remains approximately the same?

Is the ratio sin I /sin R always remains aroximatly the same?

Is the ratio hc-1/hc-r always remains aroximately same?

You will observe that sin I /sin p is the same as hc-i/hc-p as both angles have the equal length hypotenuse.

This phenomenon was first discovered and explained by Snell, and thus is named after him.

From the table of obsertations, refractive index of material of the semi circular block with respect to the air (m = can be calculated)

Experiment 15: refraccion- rectangular block (measurement of refractive index of the material of rectangular block)

Allow a single ray to fall perpendicularly at an angle of 90 at the center of the long side of the rectangular block as shown in figure 23.


Figure 23

Do you observe any deviation in the light ray either entering into or emerging out of the block if not why not?

Here, the ray enters the other medium (material of rectangular block) from air mediums as it first strikes the incident surface of the block and vice versa while leaving the block. Therefore, what other conditions according to your observations apart from the change in medium in the ray path, must be met for the deviation to the ray refraction to take place?

Mark this position of black and the point of incidence of the ray on the block.

As shown in figure 24, carefully rotate the light box so that the incident ray strikes the long face of the rectangular block at the same point of incidence at an angle other than 90



Figure 24

Mark the incident and the emergent paths of the ray. Now, remove the block and complete the actual ray path by joining the point of incidence and point of emergence of the ray. Also draw the theoretical ray path which would have been followed in the absence of block. You will obtain a diagram similar to Figure 25.

- Observe, whether the emergent ray is parallel to the incident ray and its extension?

Draw the normal (perpendiculars, N) to the interfacing surfaces of the block, where the ray enters or emerges from the block. Mark the angle of incidence and refraction with respect to the normals at both the interfacing surfaces and measure them using a protractor. For the second interfacing surface, the refracted ray becomes incident ray and the emergent ray is the refracted ray.

Check whether the magnitude of refraction at the first interfacing surface gets reversed at the second interfacing surface?

As in previous experiment, draw circles around each point of the refraction and draw and measure the half chords to the normal at these points.

Using the equation $\sin i / \sin r = n_2 / n_1$ compute the refractive index from air to the material of block and also, the refractive index from block material to the air.

Try and answer the following questions

Which of the two refractive indices is greater than 1, and which one is smaller than 1?

What do you expect the product of the two refractive indices to be? What is it actually?

Is one refractive index inverse of the other?

Experiment 16: refraction- rectangular block (lateral displacement and apparent depth)

As shown in figure 26, corresponding to the observations of the previous experiment, “d” represents the lateral displacement of the emergent ray with respect to the incident ray. For a given angle of incidence, this displacement is directly proportional to the thickness, of the slab “t”.


Figure 26

Allow a single ray of light to fall on the shorter side of the slab. Plot the emergent ray on the opposite side and measure “d”. also measure the depth of the slab “t” along which the ray is traveling, and find the ratio d/t . now, allow the ray to strike the longer side of the slab at the same incident angle and measure d and t in this case, with the same slab and find their ratio. Use other slabs of the same material and do the same. Is d/t the same in all cases?

Real depth / apparent depth (this is strictly true only when i , is close to 90°) This phenomenon is observed in day-to-day life when an underwater object appears to be close to the surface than it actually is.

Experiment 17: refraction- total internal refraction, semi circular block

Position the semi-circular block in front of the light box so that a single ray is incident on the curved surface of the block at its normal mid-point and passes through the center of the flat side of the block.


Figure 27

Do you observe any refraction at either of the interfacing surfaces? If not why not?

Mark the position of the center of the straight face of the block, from which ray is emerging, and rotate the block slightly about this point until the ray inside the block meets the flat side at an angle.

Is the ray at this flat face refracted towards or away from the normal at the point of incidence?

Again, rotate the block further about the same center point of the flat side of block and marks and observe the emergent ray at this new position. Continue this rotation of the block gradually and you will observe that at a particular position of the block no refracted ray is emerging out of the flat face. Mark this position of the block and the incident ray entering the curved face block.

What is the angle of the incidence of the ray at the flat face at this position of block?

Rotate the block a little further and observe the ray.

What happens to the ray that was being previously refracted?

On rotating back the block to its previous position, you will again observe the ray to be refracted instead of internally reflected. Once again find the position of the block at which the refraction phenomenon just ceases, initiating the internal reflection.

Measure the angle of incidence at this position. This angle of incidence is called critical angle i_c for the material of the block. The critical angle for the material is defined as the angle of incidence beyond which a light ray entering the rarer medium from denser medium suffers total internal reflection in place of refraction.

You will notice that the further rotation of the block results in the increase in intensity of the reflected ray, while the refracted ray becomes weaker.

At $i = 90^\circ$, the refracted ray completely disappears while the reflected ray is very strong.

At this position, draw a circle around the block about the incident point on the flat face and mark the half chords corresponding to the critical angle.

Also find the value of $(\sin i_c)$.

Also compute $1/\mu$ from the value of μ obtained in experiment 14. Compare the three values.

The refractive index of the slab is also given by, $\mu = \sin r / \sin i = 1 / \sin i_c$

The critical angle for water to air is 49° . By what angle should a fish rotate its head to see the riverbank, if it initially sees the opposite river-bank?

Total phenomenon of internal reflection is used in the optical-fiber technology, a well-developed method of transporting light and other electromagnetic waves over a distance.

Experiment 18: refraction- total internal refraction, triangular prism

The phenomenon of total internal reflection can be observed by using any of the prisms. Allow a single ray of light to fall on the shortest side of the 30-60-90 prism as shown in figure 28 so that the refracted beam inside the prism strikes the hypotenuse.


Figure 28

Carefully adjust the light box and prism positions until incident ray gets totally internally reflected in place of refraction and the ray emerges through the third side. Mark the first position of the prism where this takes place and measure the angle of incidence with respect to the normal to the hypotenuse.

Do you find any similarity with the value of critical angle obtained in the previous experiment (since, material of the semi-circular block and the triangular prism is same)?

Using another prism 45-45-90 make the single ray fall perpendicularly on to the shorter side of the prism so as to be internally incident on the hypotenuse of the prism as shown in figure 29.

Figure 29

Does the ray get reflected?

Mark the position of the prism and the incident as well as emergent rays.

Again reposition this prism and to make the single ray strike the hypotenuse perpendicularly at about a quarter the distance from one end as shown in figure 30.


Figure 30

In this case, you will observe that the ray gets totally internally reflected two time and returns along the path parallel to the incident ray, but reversed in direction.

Trace the prism position and the ray path through it.

Allow two parallel rays to fall normally on the hypotenuse face of the 45° prism near one end at the center of the upper half of the face and trace the path of the each ray.

Are the ray positions reversed due to the double reflections?


Figure 31

This phenomenon is observed in 45-45-90 prism is used to invert images and even reverse the direction of light.

The other prisms can also be experimented with for studying the total internal reflection. Try to obtain at least five positions of the three prisms besides the above two, which result in internal reflection. Trace the ray diagrams, using different colour pencils for different rays.

Use a pair of 45-45-90 prisms and aim a set of two parallel rays to fall on one of the prism as in figure 31. Position the other prism with its hypotenuse facing the reflected rays from the first prism and adjust its position so as to make the rays fall at the center of upper half portion of hypotenuse as in the first prism.

Trace the path of each ray undergoing quadruple reflection two reflections in each prism.

Does the quadruple reflection of the light rays restore them to their original relationship, or are they still reversed? Explain , why?


Figure 32

Compare the practical results with the theoretical ones as shown in figure 32.

Pairs of prisms arranged in similar configuration are used quite often in many optical instruments, such as prismatic binoculars. This facilitates shortening the overall length of long telescope from say 50 to 75cm as in case of binoculars.

The phenomenon of total internal reflection in a 45-45-90 prism can also be used to obtain inversion with out deviation as shown in figure 33.


Figure 33

Setup the prism with three light rays incident on the shorter side, but parallel to the hypotenuse and observe the emergent rays coming out the opposite face after suffering reflections as well as total internal reflection.

Does your observation corroborate the result as is evident in figure 32?

Try to set up the three prisms in such a manner that three rays are reflected internally by all three prisms. (optical-fiber technology is all about repeated total internal reflection).

Many faint rays are seen while examining total internal reflection. What could be the reasons for these rays?

Mirages are also an effect of total internal reflection of light due to different layers of heated air having different densities. Can you explain how?

Experiment 19: refraction- double refraction, angle of minimum deviation.

Select equilateral triangular (60-60-60) prism and position it in front of the light box. Aim a single ray of light on to one face of the prism in such a way that it is almost parallel to the face adjacent to the incident face (as shown in figure 34)

Mark the rays and the prism position

Carefully, rotate the prism to make the incident ray strike the incident face at different angle of incidence with respect to the normal to the face, but striking the face at the same point of incidence.

What changes do you observe in the emergent ray coming out of the opposite face?

Repeating the same procedure as above, find the position of the prism corresponding which you get minimum deviation of the ray with respect to the direction of incident ray striking the outer face of the prism.

Repeat this experiment with 45-90-30 angle at position A using different prisms and in each case, find the angle of minimum deviation.

The angle A is called angle of the prism. From your observations, try and answer the following questions.

How does the angle of prism "A" affect the angle of minimum deviation?

What happens when 90° angle is used in position A?

Theoretically, if "D" is the angle of minimum deviation, then refractive index of the material of prism is given as

$$\mu = \frac{\sin\left(\frac{A+D}{2}\right)}{\sin\left(\frac{A}{2}\right)}$$


Figure 35

For a prism, measure the angle of incidences and angle of refractions as shown in figure 35.

Compute using the above formula and also by $\mu = \frac{\sin i}{\sin r} = \frac{\sin i_2}{\sin r_2}$ and compare the different values obtained.

At minimum deviation, whether i_1 equals r_2 and i_2 equals r_1 . Explain why?

Is "D" equal to $((i_1+r_2)-(i_2+r_1))$.

Experiment 20: refraction- double refraction, using parallel rays

As in previous experiment, allow a single ray to fall on the triangular prism as shown in figure 33. Let the side or face of the prism opposite angle A be denoted as the base of the prism. Carefully observe and answer the following.

Does the doubly refracted ray bend towards or away from the base?

Repeat this experiment with each angle of each triangular prism at angle A position in turn.

Whether the direction of bending of ray towards the base remains unchanged for any angle at position A?

Allow a pair of parallel beams to pass through each of the triangular prisms in turn and observe the emergent rays.

For each of the three prism angles of 30, 45 and 60 whether the parallel incident rays still remain parallel to each other after emerging out of the prisms?

Experiment 21: refraction- double refraction using colour dispersion

Using the wide slit former, allow a wide beam of light to be incident on one face of the equilateral prism. Adjust the prism so as to obtain maximum deflection of the ray. On a white card using it as a screen, observe the emergent spectrum.

Is the original beam from the light box white or coloured?

Whether the emergent beam contains any colour?

Which edge of the emergent beam is red and which is blue?

If more colours than red and blue are present, make a list of the order in which they appear. Such a spread of colours, caused by the dispersion of white light is called spectrum

By adjusting the prism also observe the emergent spectrum corresponding to different angles of deviation (from minimum deviation to near-internal reflection) on the white screen.

What difference is there between the spectra obtained at minimum deviation and near-internal reflection? Hint: how many colours are seen in each case? Why?

Would you say that refractive index of the material of prism is the same for all colours? If not, arrange the colours in increasing order of.

A prism is just a parallel-slab cut diagonally. Why, then, is there no dispersion of colours in a slab?

Spectrum is caused because white light is composed of many colours, and each colour bends slightly different amount during diffraction. This suggests that refractive index of a material varies for different colours. This slitting of light into different colours by the prism is called dispersion.

The shape of rain-drops resembles very much like prisms as shown in figure 36. Hence, they also disperse white light from the sun in its constituent colours resulting in the formation of colorful rainbows in the sky during rainy seasons.

Experiment 22: refraction- colour absorption by filters

In the setup as used in previous experiment, place a red filter in the groove at the front window of the light box, in front of the wide beam slit former plate.

Here the prism is positioned in the path of incident beam between the light source and the prism.

What change do you observe in the emergent spectrum?

Remove the filter from the window of the light box and position it just beyond the colour dispersed beam emerging from the prism, so that this dispersed beam passes through the filter before falling on the white card.

What happens to the colour which are not red?

Repeat this procedure using the other filters included in the kit, one by one in turn complete the following table of observations by putting mark in the appropriate box.

		Basic Constituents of White Light						
		Red	Orange	Yellow	Green	Blue	Indigo	Violet
RED Filter	Transmits							
	Absorbs							
ORANGE Filter	Transmits							
	Absorbs							
YELLOW Filter	Transmits							
	Absorbs							
GREEN Filter	Transmits							
	Absorbs							
CYAN Filter	Transmits							
	Absorbs							
BLUE Filter	Transmits							
	Absorbs							
MAGENTA Filter	Transmits							
	Absorbs							
VIOLET Filter	Transmits							
	Absorbs							

Experiment 23: refraction- newton’s experiment with spectra

Produce a clear spectrum using the colour dispersion caused by double refraction through the equilateral prism adjusted at the position of minimum deviation with respect to the incident beam from the light box. Replace the wide beam slit former with the single slit narrow beam former.

Perform Sir Issac Newton’s experiment by positioning another equilateral prism in the path of the colour dispersed beam emerging from the first prism (as shown in figure 37). Ensure that the faces of both the prisms are parallel and sufficiently close together.


Figure 37

This second prism again refracts the beam in the opposite direction and recombines the colour of the spectrum to make white light. The mixing of colours to produce white colour can also be observed from the newton’s color disc, which has coloured sectors of the basic colours and on spinning you see all the colours merging to form white colour.

Experiment 24: refraction- double refraction using bi-convex lens

In the experiment 20, you observed that a single prism can not bring parallel rays to meet at a focus since, the rays maintain their parallelism even after refraction. But, it is clearly evident that combination of sets of prisms can be used to focus the parallel incident rays so as to meet at a point.

Block the middle two slits of a four slit former using tape and mount it in the window of the light box to get to widely spaced parallel rays of light. Position the prisms at 60° and 30° angles as shown in the figure 38. Trace the path of both the emergent rays.


Figure 38

Now again rearrange the prisms as shown in the figure 39, using a pair of identical equilateral prisms. You will observe that this arrangement converge both the parallel rays to meet at a point.


Figure 39

If another pair of parallel to the existing ones are made to fall of the prism combination, do they all focus at the same spot?

You will observe that in order to focus all the parallel rays at a point, prisms with different angles are needed, or in other words a single prism is needed where the angle is continuously changing from a small angle to a large angle. Such an arrangement is provided by a lens.

Allow a set of four rays to fall on the face of the thinner lens included in the kit. Mark the position of the lens the incident and emergent rays paths and the point at which the rays converge to meet. This point to convergence is called focal point.


Figure 40

If the local length of a lens is the distance from the center of lens to the focal point what is the focal length of this lens?

Similarly, try to find the focal length of the thicker lens of similar shape included in the kit. Lenses of this shape are called bi-convex because their both sides are bulging outwards.

Experiment 25: refraction- radius of curvature

Trace the curve of the perimeter of one side of the thinner lens surface on a sheet of plane paper. Move the curved surface of the lens along this trace thereby extending this tracing a number of times until it forms a complete circle. Measure the diameter of the circle and compute its radius.

Repeat the procedure using the thicker biconvex lens.

Each of the radius of the circle is called the radius of curvature of the respective lens used to draw that circle.

From the value of focal length of the lens determined as in experiment 23, compare them with the respective radius of curvatures as obtained above.

How does the radius of curvature affect focal length of the lens?

If a lens with a focal length of 50mm is needed, what should be the radius of curvature of the lens?

You will observe that the radius of curvature of a lens need not be the same for its both curved surfaces. The lens maker equation, as given below, gives the relationship between the refractive index of the material of the lens, radius of curvature of the two curved surfaces of the lens and the focal length of the lens.

$$\frac{1}{f} = (n - 1) \left(\frac{1}{R_1} + \frac{1}{R_2} \right)$$

Experiment 26: refraction- focal line or focal plane of bi-convex lens

Position a biconvex lens in front of the light box and mark its outline.

Allow a pair of parallel rays to fall on it such that the incident rays are parallel to the axis of symmetry of the lens. Mark the incident and emergent rays and the focal point.

Adjust the light box so as to project the parallel rays at a small angle to the axis of symmetry of the lens but aimed at the center of the lens. Again mark the rays and focal point. Repeat the same procedure with rays angled at the center of lens from the other side of the axis of symmetry.

Again repeat the procedure by adjust the rays to strike the center of the lens at a larger angle from either side of the axis of symmetry.

Join all the focal points or foci obtained above with a line

What is the shape of this line straight or curved?

If we assume that the pairs of parallel rays are coming from the distant stars and passing through the lens, where should a photographic film be placed so that the star images are in focus?

The plane at which parallel rays incident to the lens at any angle meet is called the focal plane of the lens.

Some cameras have a shutter which open and closes near the lens. Others which have removable or interchangeable have a shutter like a roller blind which covers the film, except for the instant when the blind rolls up to make the exposure. Such a shutter is called focal plane shutter.

Experiment 27: refraction- spherical aberration in bi-convex lens

Allow four parallel rays to strike a convex lens (as shown in figure 41), parallel to its axis of symmetry. The inner two and the outer two rays meet at different foci. This defect is called spherical aberration.


Figure 41

In the figure the parallel rays striking the lens near the edge are bent to a nearer focus, B, than those passing through near the center, which meet at A. the magnitude of the defect, represented by AB in the diagram, is shown greatly exaggerated.

As in case of the concave mirror, plot the caustic curve for a convex lens. This curve is called a semi-cubical parabola.

Experiment 28: refraction-chromatic aberration in bi-convex lens


Pass the rays through the outer edges of the biconvex lens. You will observe that they have coloured edges after refraction.

Which edge do you observe to be red and which one blue?

Allow the parallel rays to pass through the outer edges of the lens, to make them to converge to a focus. Position a white card vertically near the point of focus and gradually move it towards and away from the lens.

At what position of the card from the lens, the image obtained on the card has reddish center and bluish edges?

Which colour has the shortest focal length? In this agreement with the dispersion obtained by a prism? Consider the lens as two prisms placed end-on-end.

This aberration is primarily due to the different refractive indices of the optical material for different colors (or wavelengths of light)

Suggest why cheap toy telescopes and binoculars have colour fringes around the images they produce. Makes of lenses for optical instruments such as, cameras, telescopes, etc, go to great trouble using combinations of lenses of different curvatures and made from glasses of different refractive indices so as to ensure that all colours come to focus at the same point. A lens combination which has elements designed to minimize the chromatic aberration is called an Achromatic lens.

For more information on aberrations, see appendix 1.

Select the lens which has two hollow faces curving in towards each other making it thinner at the centre than at the top or bottom. Such a lens is called bi-concave lens.

Allow a set of four parallel rays of light to fall on the lens parallel to its axis of symmetry. Mark the lens position and the incidental and emergent rays.

You will observe that after refracting through the lens, the rays are diverging away from each other. Remove the lens and extend the emergent rays back through the lens position towards the light box as shown in figure 43.


Figure 43

Do they appear to radiate from one point, or from several points?

A biconcave lens is a diverging lens and is said to have negative length. This is due to the reason that the point from which the rays appear to emerge is not beyond the lens, but between the light source and the lens, or in the other word, in the direction opposite to the direction of propagation of light ray.

What is the focal length of the lens used?

Such lenses find application in correcting myopia or short sightedness, wherein the vision of distant objects gets blurred. This is a condition of eye in which the image formed by the lens of the eye falls short of the light

sensitive retina. In combination with the converging lenses or convex lenses they are also used to increase their effective focal length.

EXPERIMENTS WITH COLOURS

For experiments involving colors, the other end source end of the light box is to be used. Different colored filters can be mounted in the three openings and the swinging mirrors can be adjusted in the desired position in order to achieve a rich mixing of colors on a screen placed in front of center-opening.

Experiment 30: color observations- colors of objects

In the refractions experiments, white light was split up into its constituent colours by means of dispersion through the prism to obtain the white light spectrum. In the experiment 22, various filters were placed in the path of light rays to observe which of the rainbow colours are transmitted by the filters and which are absorbed.

Close the side- openings use the mirrors and one by one, place the filters in the center- opening and observe the colour on the screen in front of it. This will familiarize you with the various colours.

Next place the colour cards, included in the kit, on the screen and observe their colours in lights of different colours obtained after filtering from the coloured filters in front of the light box.

The colour of an object is explained by stating that the white light is composed of many colours and that a white object reflects all the colours.

A red object absorbs all the colours except red, which it reflects, making the object appear red. If however, red is completely filtered from the incident light falling on the object, by using other coloured filters say cyan, the red object should absorb the incident colour transmitted by the filter cyan and should reflect nothing. Check your observations against this theory. Note your observations in the table given below.

Colour of light falling on Plate →	White	Red	Magenta	Orange	Yellow	Green	Cyan	Blue	Violet
Colour of Plate ↓									
RED									
MAGENTA									
ORANGE									
YELLOW									
GREEN									
CYAN									
VIOLET									
BLUE									

Next, illuminate the various color cards with white light and observe them through coloured filters by holding the filters close in front of your eyes.

Tabulate your results in a similar manner.

Colour of Filter Used →	White	Red	Magenta	Orange	Yellow	Green	Cyan	Blue	Violet
Colour of Plate ↓									
RED									
MAGENTA									
ORANGE									
YELLOW									
GREEN									
CYAN									
VIOLET									
BLUE									

After familiarizing yourself with the different colours, perform some experiments on your own observing a colour card illuminated by light of one colour and viewed in another color. Try to predict your observation before performing the experiment.

Experiment 31: colour observations-addition of colors

By using the side-openings windows additionally up to three colours can be projected onto the screen. Place the weakest color filter in the center aperture and the stronger ones on the sides so that intensity loss due to reflection from the side mirrors is compensated. In order to blank off any of the openings, use blank slit former the end without slits.

Initially, take any set of filters and note down the different colours formed by their overlapping.

Swing the mirror back and forth to obtain a good overlap of the beams and observe the result of addition of colours of various color combinations.

Next, try to find complementary colors, i.e. colour combinations that reproduce white light, or close to it. Note down the colours and their complement (the complement of colour A is the colour which on combining which on combining with it gives white light)

Add three colours to make white, then remove any one of them by moving its mirror. The color left on the screen is the complementary color to that removed.

Note: remove the color filter from the light box after use, since their prolonged exposure to the heat from the lamp may damage the filter.

Experiment 32: color observations- color shadows- umbra/penumbra

Remove the coloured filters and close the side-openings.

Allow the light emerging from the rear aperture to fall on the white screen about 30 to 45cm distant from the light box. Place a pencil in front of the

screen. A sharp shadow within a faint one will be observed. If the faint shadow is too faint, use a tracing-paper to diffuse the light coming from the light box.

The sharp shadow is called UMBRA and is the region where no light falls. The faint shadow enclosing the umbra is the PENUMBRA. This is the region where partial light falls.

The phenomenon of umbra-penumbra can also be observed in everyday-life in your own shadow.

Try to explain the phenomenon with a ray diagram.

Experiment 33: colour observations-colour shadow

Place a set of complementary colors in the openings and adjust the mirrors to obtain the white light on the screen.

Place a pencil about 8-10 cm in front of the screen and in a path illuminated by all three colour beams. Make the necessary adjustment, if needed.

Observe the shadows formed and record the colors of the shadows and the colors of their background. Explain these colored shadows.

Shift the pencil to patches illuminated by two colors and then by one color and observe the shadows. The number of shadows formed is the same as the colours illuminating the patch, while their colours are the same as a few of the patches.

The above phenomena occur with any set of colours (not necessarily complementary). Try to explain the number of shadows and their colours. A ray diagram will prove helpful here, too).

Appendix: the aberration of lenses

Under ordinary circumstances, a lens or a system of lenses will not produce a perfect image. The defects of the image arise in various ways from geometrical causes as well as from the physical properties of glass, and they become especially troublesome when the lens is used at a large relative aperture to cover a wide field of view. Even in unusual cases, lens designers are concerned with six types of aberrations, two of which are briefly described below.

It is impossible to correct a system of lenses for all aberrations at the same time. Sometimes the process of eliminating one defect intensifies another. Furthermore, a lens designed to give good images of objects at one distance will give imperfect image at others. Hence it is impossible to design a universally useful lens; each problem in lens design requires individual study. If an exceedingly sharp image is required, the system needs especially fine spherical and chromatic correction; and this can only be achieved by sacrificing the field of view to reduce the aberrations known as coma and astigmatism. A wide-angle photographic lens, on the other hand, calls for freedom from coma and astigmatism, and then the aperture of the lens must remain small to avoid spherical and chromatic errors. This explains why wide-angle lenses rarely have relative apertures (values of D/f) greater than $1/10$, while ordinary photographic lenses with smaller fields of view may have relative apertures as high as $1/15$.

Spherical aberration; this name is applied to the defect occurring in the image when the lens is large and the object on the axis of the lens, and light of only one colour is used. It is a property of lenses with spherical surfaces that rays from a point on the axis are not brought together at a single point focus unless the lens has a very small relative aperture.

The magnitude of spherical aberration is proportional to the square of the relative aperture for lenses of the same focal length, or to the focal length of the lenses of the same relative aperture. For all sorts of lenses the magnitude is proportional to the square of the diameter divided by the focal length. The magnitude also varies with the shape of the lens. For a simple lens of 10cm diameter and 100cm focal length (relative aperture 1/10), the magnitude is 1.67cm if the lens is symmetrically double-convex; 4.5cm if the lens is plano-convex with the plane side toward the distant object; and 1.17 cm if the plane side is turned toward the image.

It is possible to avoid, or at least reduce, this problem in three ways.

Reduce the diameter of the lens, or place a diaphragm in the beam near the lens so as to cut out all but the central rays. The amount of light transmitted by the lens is then greatly reduced.

This is undesirable.

Choose a lens of the best shape by combining lenses of different shapes. It is possible to combine a strongly convergent lens of a shape producing a small amount of aberration with a weaker diverging lens, which gives the same aberration in the opposite direction without at the same time completely annulling the convergence of the rays produced by the first lens.

Abandon the spherical surface and grind the lens to a slightly different shape. Aspherical lenses are sometimes used in spectacles, microscopes and projection lanterns.

Chromatic aberration: the bending of light by a lens varies with the colour. A simple lens, if it had no other defects, would bring the blue range of white light to a nearer focus than the red. The yellow rays, which affect the eyes

the most, come in between. A screen placed so as to catch the yellow focus would show an image with a purple border, the edges being bluish.

The correction of this defect is accomplished by combinations of lenses of different sorts of glass. Flint glass for instance, spreads the colours apart more widely than does the crown glass. A diverging lens of flint glass placed close to a converging lens can bring the spread coloured rays back to parallelism.