BANCO ÓPTICO
OPTICAL BENCH
REF.QLG004

INDEX OF LANGUAGE- ÍNDICE DE IDIOMAS

ESPAÑOL.....2
ENGLISH.....36

Banco óptico

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Optical bench

An optical bench made of well-painted square cross sectional aluminum pipe with closed ends, fitted on two moulded supports with levelling screws. A scale graduated in millimeter fitted on the top of the optical bench. There are specially designed, easy sliding tilt proof riders with index mark for mounting various accessories, height of accessories adjustable with the help of knobs fitted in the riders. We can perform the experiments related to the fundamental of optics, such as image formation, reflection/refraction through the optical elements.

Experiments using optical bench

1. Focal length of different mirrors
2. Focal length of a thin lens
3. Focal length of a combination of lenses
4. Cardinal points of a lens system
5. To construct a telescope and to determine the magnification
6. To construct a microscope and to determine the magnification
7. To verify inverse square law
8. Fresnel biprism experiment
9. Wavelength of He-Ne laser
10. Thickness of a thin paper
11. Diameter of lycopodium powder particle using corona rings
12. Comparison of luminous intensities of two sources of light

Determination of the focal length and strength of convex lens and concave lens

Aim: to determine the focal length of a convex lens.

Apparatus required: an optical bench with three uprights, convex lens, lens holder, two optical needles.

Theory


The focal length of the convex lens is given by the relation

$$\frac{1}{f} = \frac{1}{v} - \frac{1}{u}$$

Where f is the focal length of the convex lens

U is the distance of the object needle from the optical center of the lens

V is the distance of the image from the optical center of the lens

Procedure

Determination of the rough focal length

1. Mount the convex lens in the lens holder
2. Go out in the open and face the lens towards the distance building
3. Obtain the image of the building on a white painted wall or screen and move the lens forward and backward to get a sharp image on the wall.
4. Measure the distance between the lens and the wall. This will be equal to the rough focal length of the lens.

To set the lens

1. Clamp the lens holder in a fixed upright and keep the lens in the holder and keep the upright at about 60cm mark.
2. Adjust the lens such that its surface is vertical and perpendicular to the length of the optical bench.
3. Keep the upright in this position throughout the experiment.

To set the object needle

1. Mount a needle in the upright at the extreme left of the upright holding the lens
2. Move the object needle upright and keep it at a distance about 1.5 times the focal length of the lens from the upright holding the lens.
3. Adjust the height of the needle such that the tip lie on the horizontal line through the optical centre of the lens.
4. Note down the distance of the object upright on the optical bench.

To set the image needle

1. Close the left eye and see with the right eye from the other end of the optical bench. An inverted and enlarged image of the object needle will be seen. Tip of the image must lie in the middle of the lens.

2. Mount the second needle in the upright in the right of the optical lens
3. Adjust the height of the image needle so that its tip is seen in line with the tip of the image when seen with the right eye.
4. Move the eye towards right. The tips will get separated.
5. Remove the parallax and note the position of the image needle upright on the optical bench.
6. Repeat the experiment for different values of the object position.

Observations

Rough focal length of the convex lens =cm

Table for u and v

Sr. No.	Position of the lens (cm)	Position of the object from lens (cm)	Position of the image from lens (cm)	$\frac{1}{u}$ (cm ⁻¹)	$\frac{1}{v}$ (cm ⁻¹)
1.					
2.					
3.					
4.					
5.					

Calculations: calculate the focal length of the lens using the relation for each observation

$$\frac{1}{f} = \frac{1}{v} - \frac{1}{u}$$

And calculate the mean focal length

Precautions and sources of error

1. The tips of the object needle and the image needle should lie at the same height as the center of the optical lens.
2. Parallax should be removed.
3. The object needle should be placed at such a distance that only the real and inverted image of it should be formed.
4. Uprights should be vertical

Determination of the focal length and strength of convex lens and concave lens

Aim: to determinate the focal length of a concave lens

Apparatus required: an optical bench with three uprights, convex lens, lens holder, optical needles.

Theory: a concave lens always forms a virtual image. Its focal length cannot be found directly as that of the convex lens. For this purpose indirect method is used. An object needle is placed on one side of a convex lens and its real inverted image is located on the other side. The concave lens is placed between the convex lens and the image needle. The concave lens diverges the ray and the image is now formed. The focal length of the convex lens is given by the relation.

$$\frac{1}{f} = \frac{1}{v} - \frac{1}{u}$$

Where f is the focal length of the concave lens

U is the distance of the object needle from the optical center of the lens

V is the distance of the image from the optical center of the lens

Procedure:

Determination of the rough focal length:

1. Mount the concave lens in the lens holder
2. Go out in the open and face the lens towards the distant building
3. Obtain the image of the building on a white painted wall or screen and move the lens forward and backward to get a sharp image on the wall.
4. Measure the distance between the lens and the wall. This will be equal to the rough focal length of the lens.

To set the lens

1. Clamp the lens holder in a fixed upright and keep the lens in the holder and keep the upright at about 60 cm mark.
2. Adjust the lens such that its surface is vertical and perpendicular to the length of the optical bench.
3. Keep the upright in this position throughout the experiment

To set the object needle:

1. Mount a needle in the upright at the extreme left of the upright holding the lens.
2. Move the object needle upright and keep it at a distance about 1.5 times the focal length of the lens from the upright holding the lens.
3. Adjust the height of the needle such that the tip on the horizontal line through the optical center of the lens.
4. Note down the distance of the object upright on the optical bench.

To set the image needle:

1. Close the left eye and see with the right eye from the other end of the optical bench. An inverted and enlarged image of the object needle will be seen. Tip of the image must lie in the middle of the lens.
2. Mount the second needle in the upright in the right of the optical lens.
3. Adjust the height of the image needle so that its tip is seen in line with the tip of the image when seen with the right eye.
4. Move the eye towards right. The tips will get separated.
5. Remove the parallax and note the position of the image needle upright on the optical bench.
6. Repeat the experiment for different values of the object position.

Observations

Rough focal length of the concave lens= ...cm

Table for u and v

Sr. No.	Position of the concave lens (cm)	Position of the object from concave lens (cm)	Position of the image from concave lens (cm)	$\frac{1}{u}$ (cm ⁻¹)	$\frac{1}{v}$ (cm ⁻¹)
1.					
2.					
3.					
4.					
5.					

Calculations: calculate the focal length of the lens using the relation of each observation.

And calculate the mean focal length

Precautions and sources of error

1. Focal length of the convex lens must be less than that of the concave lens so that the combination is convex.
2. The lenses must be clean.
3. The tip of the object needle and the image needle should lie at the same height as the center of the optical lens.
4. Parallax should be removed.
5. The object needle should be placed at such a distance that only the real and inverted image of it should be formed.
6. Uprights should be vertical.

Studying newton rings for the determination of wavelength of light

Aim: to determine the wavelength of sodium light using newton’s rings method.


Apparatus required: a travelling microscope, sodium lamp arrangement, newton ring apparatus consisting of optically plane glass plate and convex lens of about 100cm focal length placed in an arrangement having an optically plane glass plate inclined at 45, a spherometer, convex lens of short focal length.

Sodium lamp arrangement: it consist of a U shaped glass tube with two tungsten electrodes in the form of spiral coated with BaO. This whole of the tube is enclosed in a vacuum jacket. The U tube is filled with Neon gas at 100mm Hg pressure and some specks of metallic sodium are deposited on inner surface of tube. To start the discharge, a potential of 400V is applied through a step up transformer. These lamps are available at 35W, 80W and 140W. the discharge first passes through the Neon gas and produces pink colour.

The heat is generated which vaporized the metallic sodium specks, which expands to fill whole of the tube volume. The passage of discharge through the sodium vapour produces persistent yellow colour light. This light source emits two close lying wavelengths of 5890A and 5896A and for practical purposes considered to be monochromatic in nature. The sodium lamp is connected to the mains through a special AUTOLEAK transformer, which has high leakage reactante resulting drooping current-voltage characteristics.

Theory: circular interference fringes produced by enclosing a very thin film of air of varying thickness between the surface of a convex lens of large radius of curvature and optically plane glass plate are known as Newton's rings. In order to produce these fringes, light from an extended monochromatic source is rendered parallel by a convex lens. It falls on the glass plate G inclined at an angle of 45 to the vertical and is reflected

normally on the lens N. an air film of gradually increasing thickness is enclosed between lower surface of this lens and glass plate P. the light reflected from upper and lower surface of the air film produces interference fringes. At the center the lens is in contact and thickness of air film is zero. The center will be dark as a phase change of is introduced due to reflection at the lower surface of the air film. As we proceed outwards from the center, the thickness of the air film increases gradually being same all along the circle with center at the point of contact. Hence the fringes are concentric circles and are localized in air film. These fringes are viewed by means of low power microscope.

For maximum intensity or bright ring the condition is

$$2\mu t = \frac{(2l - 1)\lambda}{0}$$

The thickness of the film can be expressed in terms of radius of convex surface of lens R and radius r of the circular fringes.

$$t = \frac{r^2}{2P}$$

Hence diameter of bright ring, taking $\mu = 1$ for air film, is

$$D_n^2 = 2(2n / \mu) \nu P$$

For minimum intensity or dark ring the condition is

$$2\mu t = (2n + 1/2)\lambda$$

Hence diameter of nth dark ring, taking $\mu = 1$ for air film, is

$$D_n^2 = 2n\nu P$$

The wavelength of light can be determined by measuring the diameter of different or bright rings by the relation given as

$$\lambda = \frac{D_m^2 - D_n^2}{4(m/n)P}$$

Procedure

1. Level the microscope table and set the microscope tube in the vertical position. Find the Vernier constant of the horizontal scale.

2. Clean the surface of the glass plate P, the lens N and the glass plate G1 and fix them in their respective positions. Place the arrangement in front of the sodium lamp such that the center of sodium lamp opening coincides with the center of convex lens (acting as a collimator) and that of the glass plate G2. The inclination of the glass plate is adjusted so that light falls at 45° angle of incidence on it after its collimation from the convex lens.
3. Adjust the position of microscope so that it is vertically above the center of lens N. focus the microscope to observe alternate dark and bright rings are clearly visible.
4. Adjust the position of microscope till the point of intersection of two cross wires coincides with the center of ring system and one of the cross wires is perpendicular to the horizontal scale.
5. Slide the microscope to the left till the cross wire lies tangentially at the center of 4th dark ring. Note the reading of the Vernier of the microscope. Slide the microscope backward by slow motion and note position of 6, 8, 10 and 12th dark rings respectively. Keep on sliding to the right and note the readings when the cross wire again lies tangential to the center of 4, 6, 8, 10 and 12th dark rings.
6. Remove the lens N in contact with the glass plate P and find its radius of curvature using spherometer.

Observations

- a) Determination of diameter of Newton rings
Vernier constant of linear scale

Ring No.	Microscope reading = LSR (in mm) + VC x VSR		Diameter of ring
	Left	Right	
4			
6			
8			
10			
12			

LSR= linear scale reading

b) determination of focal length of plano-convex lens

least count of spherometer= mm mean distance between two legs= L=

Q.L.m	Spherometer reading on		h = RCS - PCS
	Convex surface (RCS)	Plane surface (PCS)	
1			
2			

Radius of curvature of lens

$$\frac{l^2}{6h} + \frac{h}{2}$$

Precautions

1. The lens and glass plates should be cleaned
2. A lens of larger radius of curvature preferably plane-convex type should be used. This leads to a wedge shaped film of air of very small thickness between lens and glass plate.
3. The point of intersection of cross wires should coincide with center of rings system
4. The micrometer screw should always be moved in same direction to avoid backlash error.
5. The radius of curvature of the surface of lens in contact with the glass plate should be measured accurately.
6. The amount of light from the source should be adjusted to maximum visibility. Too much light causes general illumination and decreases the contrast between dark and bright fringes.

Verification of the inverse square law for light radiation and determination of the absorption coefficient of light using photocell.

Aim

To determine Planck's constant and verify inverse square law using photocell.

Apparatus required: photocell attached with optical bench, light source, regulated power supply, colour filters.

Theory

In the phenomenon of photoelectric effect when light of frequency greater than the binding energy of the electrons of the metal is bombarded on the metal, the electrons are emitted from the metal. The emission of photoelectrons is possible only when the frequency of the incident light is more than the critical frequency, known as the threshold frequency. The kinetic energy of the emitted electron depends upon the energy of the incident light but is independent of the intensity of the incident light. However, the number of photoelectrons emitted depends on the intensity and is independent of the frequency of the incident light.

A photoelectric cell is a device, which converts light energy into electrical energy. The electromagnetic field is quantized i.e. light is made of quantized photons, each having energy $h\nu$ where h is Planck's constant ($h = 6.63 \times 10^{-34}$ joule-sec) and ν is the frequency of the incident light. The electrons inside the metal are bound with an energy e_0 where ν_0 is known as photoelectric work function. It is given by the relation

$$e_0 = h\nu_0 \text{ where } \nu_0 \text{ is the threshold frequency}$$

when the light of frequency (energy $h\nu$) greater than the threshold frequency ν_0 is incident on the metal, a part of this incident energy is used by the electron to come out of the metal surface and the rest is restored in it as kinetic energy

$$h\nu = \frac{1}{2} mv^2 + e_0$$

$$h\nu = eV + h\nu_0$$

$$V = \frac{h}{e}\nu - \frac{h}{e}\nu_0$$

A graph is plotted between stopping potential V and frequency of the incident light, which comes out to be a straight line. The slope of the straight line gives

$$\tan\theta = \frac{h}{e}$$

Planck constant $h = e \tan\theta$

Verification of inverse square law using photocell: if E is the intensity of illumination and L is the luminous intensity of an electric lamp at a point distant r from it, then according to the inverse square law

$$E = \frac{L}{r^2}$$

Where r is the distance between the light source and the photocell.

A plot between $1/r^2$ and E comes out to be a straight line, which verifies inverse square law of radiation.

Procedure to calculate Planck's law

- Place one filter in the space provided on the apparatus, set the suitable intensity of light and switch on the power supply
- Adjust the de-accelerating voltage to zero volt and then increase the de-accelerating voltage so that the photocurrent becomes zero. Now note this value of the de-accelerating voltage in the table.
- Repeat the same procedure for all the filters and note the corresponding de-accelerating voltages.

Sr. No.	Filter	Wavelength (nm)	Frequency (Hz x 10 ¹⁴)	Stopping voltage (Volt)
1.	Red	635		
2.	Yellow	585		
3.	Green	500		
4.	Blue	460		

- A graph is plotted between stopping potential V and frequency of the incident light, which comes out to be a straight line. Now calculate the value of Planck's constant using the relation

$$h = e \frac{V_2 - V_1}{\nu_2 - \nu_1}$$

Procedure for verification of inverse square law

- Set the lamp and the photocell enclosed in a box on an optical bench. Connect the power supply and the galvanometer to read the deflection or microammeter to read the small current.
- Place red filter (635nm) in front of the lamp.
- Keep the position of the photocell fixed and the voltage constant. Now increase the distance of the illuminating lamp from the photocell in small steps and note down the galvanometer deflection in each step and also the position of the lamp.

Table

Sr. No.	Distance between photocell and lamp r (cm)	$\left(\frac{1}{r^2}\right) \text{ cm}^{-2}$	Current I (μA)
1.			
2.			
3.			
4.			
5.			
6.			
7.			
8.			
9.			
10.			

- Draw a plot between 1/r² and comes out to be a straight line, which verifies inverse square law radiation.

Precautions and sources of error: 1 the experiment should be performed in the dark room to avoid any stray light to photocell.

2. the observations should be taken by altering the anode potential in small steps of 0.05volts.
3. corresponding to the zero anode potential, the deflection of the light spot on the scale should be adjusted at its maximum value
4. the smoth graphs should be plotted
5. the stopping potentials should be read carefully

Diffraction gratings experiment for the determination of wavelength of light


Aim: to determine the wavelength of sodium light using diffraction grating

Apparatus required: a spectrometer, spirit level, diffraction grating with clamping arrangement sodium lamp arrangement.

Sodium lamp arrangement: it consists of a U shaped glass tube with two tungsten electrodes in the form of spiral coated with BaO. This whole of the tube is enclosed in a vacuum jacket. The U-tube is filled with Neon gas at 100mm hg pressure and some specks of metallic sodium are deposited on inner surface of tube. To start the discharge, a potential of 400V is applied through a step up transformer. These lamps are available at 35W, 80W and 140W. the discharge first passes through the Neon gas and produces pink colour. The heat is generated which vaporizes the metallic sodium specks, which expands to fill whole of the tube volume. The passage of discharge through the sodium vapour produces persistent yellow colour light. This light source emits two close lying wavelengths of 5890A and 5896 and for practical purposes considered to e monochromatic in nature. The sodium lamp is connected to the mains through a special autoleak

transformer, which has high leakage reactance resulting drooping current-voltage characteristics.

Optical spectrometer: the spectrometer is compact apparatus for obtaining pure spectrum. This spectrometer and consists of following parts: i) collimator; its purpose is to produce a parallel beam of light. It consists of horizontal tube on the arm of a spectrometer, which has converging achromatic lens at one end and a vertical slit at the other on a slid tube. The focal length of the converging lens is equal to the length of the tube. The distance between the lens and slit can be varied using rack and pinion arrangement or by sliding the inner tube to obtain parallel rays. The tube rests on two screws C1 and C2 which can be slightly titled up and down for horizontal alignment. The slit is in the form of two sharp edges of which one is fixed while other can be move parallel for changing its width and length. ii) telescope: it is an astronomical telescope with achromatic objective and Ramsden's eye piece. It is mounted on another arm fixed rigidly to the circular scale graduated in 0.5° . the telescope along with the circular scale can be turned around the vertical axis passing through the center of the spectrometer. It can be fixed to any position through a screw T1 and can then be given slow rotation using tangent screw T2.

The position of telescope can be read through two Vernier scale V1 and V2 180° fixed to the prism table. The telescope rests on two screws T3 and T4 for horizontal alignment. The rack and pinion arrangement is provided on the side of tube for focusing. iii) prism table: it consists of upper plate a and a lower plate B separated by three springs, through which leveling screws P, Q and R are passed. A set of parallel equidistant lines are engraved (say between screws P and Q). the prism is always placed with one of the reflecting faces perpendicular to these lines. A series of circles concentric with vertical axis of rotation are ruled on this plate, which help in placement of the prism correctly. The height of the prism table can be adjusted by clamping the screw A1, which fixes the prism table to verniers V1 and V2. The table can rotate about the vertical axis and may be fixed at any position by means of screw A2. After fixing by A2. The table can rotate about the vertical axis and may be fixed any position by means of screw A2. After fixing by A2, it can be moved slowly by tangent screw A3 at the base of the spectrometer. The position of the table can be read on two Vernier scales.

Diffraction grating: A plane diffraction grating consists of an optically plane glass plate, on which are ruled number of equidistant parallel straight lines. These lines divide the glass plate into opacities (b =grooves created by the lines) and transparencies (a = thickness of glass plate between two opacities) the thickness of which is of the order of wavelength of light. The quantity $(a + b)$ is called grating element. The number of lines in a plane transmission grating is of the order of 6000/cm. in laboratories, actual grating is not used as it is very costly. The original grating is usually ruled on a plate of speculum metal by a sharp diamond point fitted to a dividing engine. For ordinary use, replicas of this grating are obtained by depositing a thin film of gelatin on it. This film is removed when dry and pasted on optically plane glass plate.


ABC...H in the figure shown below represents the section of plane transmission grating supposed perpendicular to the plane of paper. Let the width of clear space be a and that of ruling be b . the distance $(a+b)$ is called grating element. The points A and C separated by grating element are called corresponding points. When a parallel beam of monochromatic light is incident normally on the grating, it suffers diffraction. In this figure ABC...H represents the incident plane wavefront while AMN...K represents the wavefront after diffraction through grating at an angle. There is no path difference in reaching the plane ABC...H. The only phase change that occurs is due to difference in path traversed by different rays in reaching from ABC...H to AMN...K

The path difference between the corresponding points A and C is $(a+b)\sin \theta$ the same is true any two corresponding points. The transmitted light after falling on the lens is focused at point P. the point P will be bright or dark according as the rays reinforce with each other. They would reinforce and give brightness or a maximum if.

$$(a + b) \sin\theta_n = n\lambda$$

They would interfere and produce darkness or minima if

$$(a + b) \sin\theta_n = (2n+1)\lambda/2$$

When $n=0, \theta_n=0$, it gives central bright fringe or the first orderspectrum ($n=1$)

$$(a + b) \sin\theta_1 = \lambda$$

For the second order spectrum $n=2$

$$(a + b) \sin\theta_2 = 2\lambda$$

$$\sin\theta_2 = 2\sin\theta_1$$

The above relation is helpful in ascertaining the second order spectrum after the first as in certain cases, the second order spectrum is missing and

first is followed by the third order spectrum. If the light is not monochromatic, then each wavelength will give its own primary maxima. Hence each order will consist of a spectrum.

The number of lines N are marked on the grating. The value of the grating element is given by

$$(a+b) = 2.54/N \text{ cm}$$

Absent spectra: for diffraction through a grating, the condition for the maxima is

$$(a + b) \sin\theta = n\lambda$$

If the path difference between rays from the extreme ends of one single slit is, then intensity of light is zero. It is because in such a case a slit can be divided into two half slits and the path difference between the corresponding points on two half slits is $\lambda/2$. The condition for zero intensity is $a \sin\theta = \lambda/2$. When above given two conditions are simultaneously obeyed the beams from all the slits destructively interfere to produce zero intensity and the corresponding order of the spectrum is absent from the diffraction pattern. Thus for absence for n th order spectrum we have the condition that $a+b/a=n$

Procedure

- a) settings
- iii) adjust the position of telescope eyepiece so that the cross-wires are clearly visible. Focus the telescope on a distance object and set it for

parallel rays. Level the spectrometer by the leveling screws and then the prism table with help of spirit level.

ii) fix the grating stand on the circular table with two screws in the holes drilled on one of the lines parallel to the line joining two of the screws meant for the purpose say P and Q. the face of the stand to which the clamps are attached should be at the center of the table. Take out the grating from the box, holding it from the edge and without touching its surface, fix it carefully to the frame with the ruled surface towards the telescope.

b) Optical leveling of the grating table

- i) Rotate the table so that the plane of grating is approximately inclined at an angle of 45° to the collimator axis. Rotate the telescope to receive the light from the grating surface. If the image is not symmetrical w.r.t. horizontal cross-wire, adjust with the help of third screw R. in this position, the plane of the grating will be vertical.
- ii) Rotate the table carrying the grating so that the plane of grating is perpendicular to the axis of collimator. Look for the first order spectrum on one side of the direct image of the slit. Turn the telescope so that the first order diffracted image coincides with the vertical cross-wire. If this image is not symmetrical w.r.t. the horizontal cross-wire, adjust the level of prism table with the help of screw P and Q. in this position, grating lines are parallel to the axis of the spectrometer. Now turn the telescope to the other side to that vertical cross-wire again coincides the first order diffracted image. If careful adjustments are done, then

diffracted image of the slit will be symmetrical to the horizontal cross-wire at all positions.

c) Setting the grating normal to the incident light

- i) Place the telescope in line with the collimator so the vertical cross-wire falls exactly at the center of image of the slit. Note the scale reading. Add 90° to the reading and set the telescope at this reading to place it perpendicular to the axis of the collimator. Clamp it in this position.
- ii) Rotate the grating table till the plane face of the grating is facing both the collimator and the telescope. Turn the slowly till the center of the slit fall exactly on the vertical cross-wire. In this position the plane of grating is an angle of 45° to the incident light. Note the reading. Turn the table through 45° from this position so that the plane of the grating is normal to the incident light with its plane facing towards the collimator. The grating is now set normal to the incident light with its ruled surface away from the collimator. Clamp the table in this position.
- iii) Place the eye in the front of collimator and move it gradually toward the telescope till the first order diffraction image is visible. Bring the telescope in this position and observe the image through it. Clamp the telescope in this position. If the resolution of telescope is sufficiently high two distinct narrow lines corresponding to 589nm and 589.6nm will lie side by side in the field of view. Ordinarily two lines will appear as one in the first order spectrum. Turn the tangent screw of the telescope till the cross-wire coincides the center of the image of the slit. Note the reading on both the Vernier scales. Similarly

observe the first order diffracted image and note the readings on both the Verniers.

- iv) Similarly observe the second order diffracted image on either side of the image and note the reading on both the Verniers.
- v) Note the number of lines marked on the grating and replace it carefully with its ruled surface upwards in the box.

Observations

- i) Vernier constant of the spectrometer
- ii) Number of lines per inch on the grating
- iii) Grating element
- iv) Direct reading of the telescope

Order of Spectrum	Telescope Reading			Angle of diffraction		
	Left	Direct	Right	Left	Right	Mean
I st order						$\theta_1 =$
II nd Order						$\theta_2 =$

Wavelength of sodium light as calculated from first order spectrum = $\lambda = (a+b)\sin\theta$

Wavelength of sodium light as calculated from second order spectrum = $\lambda = (a+b)\sin\theta/2$

Mean wavelength

Precautions

1. Grating should be held from the edges and the ruled surface should not be touched.
2. The telescope should be focused on the brightest image of the slit while seeing the reflected image.
3. The ruled surface should face away from the collimator.
4. The light should fall on whole of the grating surface
5. The adjustment is incorrect if the angles of diffraction on the left and the right of the direct position are not equal.
6. In the second order spectrum with a good grating the two D lines generally appear to be separated. In such a case, the angle of diffraction should be found separately for each line and wavelengths calculated corresponding to each.

Studying specific rotation using polarimeter


Aim: to determine specific rotation of the sugar solution using a Laurent's half shade polarimeter.

Apparatus required: Laurent's half shade, sodium lamp arrangement, sugar, weighing balance, graduated cylinder, two beakers, funnel, a pipette and a glass rod.

Sodium lamp arrangement: it consists of a U shaped glass tube with two tungsten electrodes in the form of spiral coated with BaO. This whole of the tube is enclosed in a vacuum jacket. The U-tube is filled with neon gas at 100mm Hg pressure and some specks of metallic sodium are deposited on inner surface of tube. To start the discharge, a potential of 400V is

applied through a step up transformer. These lamps are available at 35W, 80W and 140W. the discharge first passes through the neon gas and produces pink colour. The heat is generated which vaporizes the metallic sodium specks, which expands and fills whole of the tube. The passage of discharge through the sodium vapour produces persistent yellow color light. This light source emits two close lying wavelengths of 5890Å and 5896Å and for practical purposes considered to be monochromatic in nature. The sodium lamp is connected to the mains through a special autoleak transformer, which has high leakage reactance resulting drooping current-voltage characteristics.


Laurent's half shade polarimeter: it is an instrument used to measure specific rotation of certain solutions. The optical parts of the spectrometer

are shown in the figure below.

Nicol prism: it is an optical device made of calcite and is used in many instruments for producing and analyzing plane polarized light. A calcite

crystal of length three times its width is taken. The principal section is shown below.


Its end faces are cut down as so to reduce the angles at B and D from 71 to 68 with the principal section. The crystal is now cut apart along AC perpendicular to both principal plane and end faces such that AC makes an angle of 90 with ends AB and CD. The two cut faces are ground, polished optically flat and cemented together with Canada balsam whose refractive index lies between that of calcite for ordinary and extraordinary rays. The sides of the prism are blackened to absorb totally reflected ray.

A ray incident on the face of prism is split up into the ordinary and extraordinary rays. The ordinary ray is totally reflected at interface coating whereas the extraordinary passes through giving rise to plane polarized light having vibrations in principal plane. The nicol prism can be used both as polarizer as well as analyzer. When the principal section of two prisms is parallel, the extraordinary ray is transmitted beam decreases and finally no light is transmitted becomes the ordinary ray for the second prism and it suffers total internal reflection.

Laurent's Half shade device: it consists of a semicircular plate ADB of glass and ACB of quartz. The quartz plate is cut with its optic axis parallel to the line of separation AOB. The thickness of quartz plate is such that it introduces a phase difference of between ordinary and extraordinary vibrations. The thickness of the glass plate is such that it absorbs same amount of light as by the quartz plate.

Suppose that light after passing through the polarizer P is incident normally on the half shade plate and has vibrations along OP. On passing through the glass half, the vibrations will remain along OP but on passing through the quartz half, it will split into ordinary and extraordinary components. The vibrations of O-component will be along OD and those of E component along OA. A phase difference of is introduced between two components on passing through the quartz plate. The O-vibrations will lead in phase by and will take place along OC instead of OD. The resulting vibrations on emergence from the quartz plate will be along OQ so that $\angle POA = \angle QOA$.

If the analyzing nicol prism is fixed with its principal axis parallel to OP, the plane polarized light will pass through glass half while it will be partially obstructed on passage through quartz half. Hence glass half will appear brighter. When the principal plane of analyzing Nicol is parallel to OQ, the quartz half will be brighter than glass half due to same reason. When principal plane of analyzing nicol is parallel to AOB, the two halves will be equally bright since vibrations emerging from two halves make equal inclination with principal plane. When principal plane of analyzing nicol is perpendicular to AOB, the components OP and OQ will have equal intensity but intensity of light emerging is very small and two halves will be equally dark. The eye can easily detect slight change in the intensity, when two halves are equally dark. The readings are therefore taken for this position.

Working of polarimeter: the light from the monochromatic source after passing through a narrow slit is rendered into a parallel beam by lens L. The light is made polarized by the nicol prism (polarizer) and after passing through half shade device, glass tube (polarimeter tube) containing the solution is made to fall the analyzer nicol prism. The light is viewed through Gaussian telescope. The analyzing nicol prism can be rotated about the axis of tube and its rotation can be measured on the graduated circular scale divided in degrees, with the help of a Vernier.

Theory

Unpolarised and polarized light: the light is a transverse wave in which the oscillations of electric and magnetic fields are mutually perpendicular to each other and to the direction of propagation. The oscillating magnetic field is very weak and plays no significant role in visual effects. Hence an ordinary light can be supposed to consist of oscillations of electric field

in all possible directions transverse to the direction of propagation. The oscillating electric vectors in various directions can be resolved into rectangular components lying in the transverse plane. The phenomenon of restricting the vibrations of light in a single direction in a plane perpendicular to the direction of propagation is called polarization of light. Such a light lacks symmetry and confirms the transverse nature of light. The plane in which the vibration takes place in the polarized light is called plane of vibration while the plane perpendicular to the plane of vibration is called plane of polarization. The polarized lights are of the following types

- a) Plane polarized light: the polarized light in which vibration of electric vector is restricted in a single direction in plane perpendicular to the direction of propagation is called polarized light.
- b) Circularly polarized light: in such polarized light, the direction of the oscillation of the plane rotates periodically but the magnitude remains constant.
- c) Elliptically polarized light: in such a light the direction of oscillations of electric vector rotates periodically but magnitude varies within certain maximum and minimum limits.

Methods of polarization: the polarized light can be produced by the following different methods.

Polarizer crystals: certain organic crystals like tourmaline can be used to produce plane polarized light. When ordinary light is passed through such crystals, only the vibrations along the optic axis of the crystal are transmitted, thereby producing plane polarized light. If this polarized light is allowed to pass through another similar crystal, the intensity of transmitted light varies as a function of angle between optic axes of two crystals. Intensity is maximum when axes of two crystals are parallel while

it becomes zero, when two axes are perpendicular. The intensity variation is expressed as Malus' law given as

$$I(\theta) = I_0 \cos^2 \theta$$

Polarization by reflection: if an ordinary beam of light is incident on a transparent medium, then its reflected component is found to be partially polarized. The degree of polarization increases with the angle of incidence and at a particular value of angle of incidence, called polarizing angle, the reflected beam is completely polarized. At polarizing angle of incidence the polarized reflected light is perpendicular to the refracted unpolarized light and the refractive index of the transparent medium is related to the polarizing angle.

$$\mu = \tan p.$$

This is referred to as Brewster's law.

Polarization by scattering, when an ordinary light passes through the medium, it suffers scattering from the particles of the medium. If the size of particles is comparable to the wavelength of light, the light scattered in the direction perpendicular to the direction of incidence gets plane polarized.

Double refraction: when the ordinary light passes through a calcite crystal, the refracted ray splits up into two components which are

- a) Ordinary ray: this component obeys the law of refraction and has vibrations perpendicular to the principal section of the crystal. The velocity of ordinary rays is the same in all directions. These rays get

polarized such the direction of oscillation of electric vector is perpendicular to the optic axis.

- b) Extra-ordinary ray: this component does not obey the laws of refraction n has vibrations in the principal section of the crystal. This component has different velocities in different directions. These rays get polarized in such a way the direction of polarization is parallel to the optic axis.

This phenomenon is called double refraction. These crystals are said to be positive in which $v > v_0$ e and those in which the opposite is negative. The double refracting crystals have property of transmitting orthogonally polarized component of light with different velocities and this property is also called Birefringence.

Optical activity: it is the phenomenon of rotating the plane of polarized light when it is allowed to pass through a material. If a material rotates the plane of polarized light towards left. It is called leavo-rotatory while those rotating the plane of polarized light towards right is called Dextrorotatory. If polarized light passes through a column of an optically active material of length L and concentration c , then the angle of plane of polarized light is

$$\theta = sLc \quad s = \frac{\sigma}{Ja}$$

This is called Biot's Law and "s" is called specific rotation, which is dependent on temperature of materia and wavelength of incident light.

Procedure:

- a) Preparation of 20% solution: take a clean dry beaker. Add about 20g of sugar and calculate the volume of solution to have a 20% strength as

$$\theta = sLc \quad s = \frac{\sigma}{Ja}$$

Add nearly half of the volume of water in the beaker and stir well till whole of the sugar is dissolved. Add more water if necessary taking care that volume of solution is less than the calculated value of the volume. Transfer the solution to the graduated cylinder. Complete the volume by adding water in little amount by using pipette. Filter the solution in a clean beaker and cover it.

- b) Setting: find the Vernier constant of the circular scale. Place the polarimeter so that its aperture is in front of sodium lamp. Look through the telescope and adjust the position of the eye piece so that two halves of half shade divide are equally dark and are clearly focused.
- c) Remove the brass caps of the polarimeter tube. Clean the tube as well as glass windows. Now replace one of the caps in position taking that no strain is exerted on the glass window. Hold the tube in a vertical position and fill it with water. Slip the second glass window gently on the tube taking care that no ir bubbles are left underneath. Screw the cap gently so that no strain is exerted on the glass window.
- d) Taking observation: place the tube in position in the polarimeter and cover it. Rotate the analyzing prism by rotating the circule scale till halves of the shade are equally dark. Note the reading on the scale. These are the measurements for 0% sugar solution and

Formula used: the thickness of a thin paper can be determined using the relation.

$$t = \frac{\lambda l}{2\beta}$$

Where λ = wavelength of light used

β = distance between two successive bright (or dark) bands and l = horizontal length of glass plate i.e. distance between edge of air wedge and tip of paper.

Procedure:

1. Find the Vernier constant of the horizontal scale of the microscope.
2. Level the base of the microscope
3. Make the microscope tube vertical
4. Arrange the apparatus with the edge of the thin plate inserted under one end of the microscope slide which is placed on the glass plate.
5. Adjust the source and the microscope tube until parallel interference fringes are formed in the air wedge and are seen in the field of view of the microscope till the horizontal cross= wire coincides with one of the centers of the straight fringes. Take the reading on the scale.
6. Slide the microscope with the help of slow motion screw and fix the cross wire on the 5, 10 and 15 fringes and take the reading on the scale.

7. Take 6 reading in all and from these determine the distance traversed for 15 fringes. Divide it by 15 so as to calculate the fringe width.
8. With the help of microscope measure the distance between the points M and N.

Observations

Vernier constant of microscope=cm

Wavelength of sodium light used $\lambda = 5893 \times 10^{-8}$

(I) Table for fringe width β

S.No.	Order of fringe	Microscope reading(cm)	Traverse for 15 fringes(cm)	Mean traverse for 15 fringes(cm)	Fringe width β (cm)
1.	n				
2.	n+5				
3.	n+10				
4.	n+15				
5.	n+20				
6.	n+25				

iii) Table for measurement of l

S.No.	Microscope reading on		l (cm)
	Thin plate N (cm)	Line of contact M (cm)	
1.			
2.			
3.			

Mean value of $l = \dots$ cm

Calculations:

Thickness of plate

$T = \lambda \cdot l / 2 \cdot \beta = \dots$ cm

Precautions and sources of error

1. The glass plate and the travelling microscope slide should be clean
2. The glass plate should be perfectly horizontal
3. The edges of the paper should be sharply cut
4. The field of view for the measurement of fringe width should be in the central of the air film.

Determination of the thickness of thin metallic wire using laser beam

Aim: to determine the wavelength of He-Ne laser using diffraction grating.

Apparatus required: He-Ne laser source, an arrangement of two slits each of width 0.075mm and separated by distance 0.0541mm a high power microscope to measure the slit width, meter scale.

Formula used

$$\beta = \frac{\lambda D}{d}$$

Where λ = wavelength of laser light obtained from He-Ne (6328Å)

D = distance of screen from the slit (S)

d = separation between the centers of two slits.

And b ? fringe width

Procedure:

1. Place the pair of slit across the diameter of the He-Ne laser and close to the output of the He-Ne laser beam source
2. Adjust the position of the laser tube horizontal.
3. Rotate the frame of slits about horizontal axes and in the vertical plane by means of the tangential screw fitted in the frame of the slits to make the slits vertical.
4. Switch on the power supply of the He-Ne laser source
5. Rotate the frame about a vertical axis so that the frame of the slits is perpendicular to the direction of travel of laser beam. The bright

and fringes of small height but equal width are observed on the screen.

6. Note the position of the centers of successive bright fringes cm and meter scale place on the screen and perpendicular to the length of the slits with the help of micrometer eye piece.
7. Find the width of each slit with the help of a microscope.
8. Find the distance between the centers of two slits with the help of the microscope.

Observations

Least count of the microscope =cm

Least count of meter rod=....cm

Table 1

Mean value of beta=cm

Reading of microscope for slit widths

S.No.	Crosswire on edge	Reading of microscope (mm)		
		Main scale	Vernier scale	Total
1.	A			
2.	B			
3.	C			
4.	D			

Width of slit $S_1 = A \sim B = \dots \text{cm}$

Width of slit $S_2 = C \sim D = \dots \text{cm}$

Distance $AD = A \sim D = \dots \text{mm}$

Distance $BC = B \sim C = \dots \text{mm}$

$$D = \frac{AD + BC}{2} \dots \text{mm} = \dots \text{cm}$$

Distance between the slit frame and screen= $D = \dots \text{cm}$

Calculation

$$\beta = \dots \text{cm}$$

$$D = \dots \text{cm}$$

$$d = \dots \text{cm}$$

Wavelength of He-Ne laser beam

$$\lambda = \frac{\beta d}{D} \dots \text{cm} = \dots \text{\AA}$$

Result

Theoretical value of wavelength of He-Ne laser beam=6328Å

Experimental value of wavelength of He-Ne laser beam=Å

Precautions and sources of error

1. Do not stare at the laser source as it may cause the damage to the eye
2. The graph used as the screen should be vertically straight
3. the diffraction grating should be very near to the laser source

aim: to determine the high resistance by method of leakage using a Ballistic galvanometer.

Apparatus Required: Ballistic Galvanometer, accumulator, morse key, two way key, standard condenser (capacity of the order of 1.0 or 0.5 μF), given resistor, stop watch and connection wires.

Formula used

The high resistance R is calculated using the relation

$$R = \frac{t}{2.3026C \log_{10} \frac{\theta_2}{\theta_1}}$$

Where

T is the time period of the leakage of condenser through the resistance

C is the capacity of the standard condenser

θ_0 is the first throw of spot of light when initially the condenser is discharged through ballistic galvanometer.

θ_1 is the first throw of spot of light when the condenser is discharged through the ballistic galvanometer after a leakage of charge for time t through R.

Procedure

- 1) make the electrical connections as in the figure
- 2) close K1 to ii and press the morse key , i.e. charge the condenser for 40.
- 3) Release the morse key K2 so that the condenser is discharged through the galvanometer. Note down the first throw θ_0
- 4) Repeat the procedure of the points ii and iii several times i.e. time charge condenser and then discharge through Ballistic Galvanometer. Obtain mean value of θ_0
- 5) Closing K1 ii and pressing morse key K2, charge the condenser for the same time. Keeping morse key pressed, open K1 ii and close K1 I . start the stop watch
- 6) After a measured time t seconds, (say 5 or 10), release morse key and note down the first throw θ_1 in the galvanometer
- 7) Repeat procedure 5 and 6 for different values of time

Observations

Capacity of condenser=..... μF

Sr. No.	galvanometer θ_0	Mean θ_0	Leakage time t (sec)	Throw in ballistic galvanometer θ_1	$\frac{\theta_2}{\theta_1}$	$\log_{10} \frac{\theta_2}{\theta_1}$
1.
2.
3.
4.
5.

Calculations: plot a graph between leakage time t on x-axis and

$$\log_{10} \frac{\theta_2}{\theta_1}$$

on Y-axis

The slope of the curve

$$= \frac{\log_{10} \frac{\theta_2}{\theta_1}}{t}$$

The high resistance R is calculated using the relation

$$R = \frac{t}{2.3026C \log_{10} \frac{\theta_2}{\theta_1}}$$

Result: resistance of the given resistor is....ohms.

Precautions and sources of error

1. The galvanometer coil should be made properly free
2. Tapping key should be used across the galvanometer
3. Condenser should be free from dielectric loss
4. After observing θ_0 , the galvanometer coil should be at rest for observing the value of θ_1 .

Aim: to determine the diameter of a thin wire using a He-Ne laser

Apparatus required: laser light source (He-Ne or semiconductor laser), an optical bench with four uprights, micrometer eyepiece, adjustable fine slit width, piece of experimental wire, clamp for holding thin wire, screen, tracing paper and pencil.

Formula used: the diameter of the thin wire is given by the relation

$$d = \frac{\lambda D}{\beta}$$

Where

d= diameter of the wire

lambda= wavelength of the laser light used

D distance between the screen and the wire

Beta= fringe width, i.e. separation between the centers of the two consecutive maxima.

Procedure

- Adjust the screen and the optical bench perpendicular to each other. Up the uprights on the bench.
- Hold the straight thin wire tightly in its clamp and mount it on an upright on the optical bench. The wire should be vertical and tightly fixed.
- Place the screen on the upright at one of the extreme end position of the optical bench.
- Put the laser source on its upright about 2 to 3 meter away at another end of the optical bench.
- Switch on the power supply of the laser source and illuminate the wire with the laser light. The upright having the wire should be near the laser source.
- Adjust the position of the wire upright and the laser source so that a well defined diffraction pattern on both sides of central maxima

is observed on the screen. The pattern may also be obtained on the wall.

- Place a long strip of the trace paper on the screen and with sharp pencil mark the positions of the maxima on this paper.
- Record the positions of the wire and screen uprights to calculate the value of D.
- Record the distance between central maxima and centers of five consecutive maxima with the help of micrometer reader to find the value of beta.
- Repeat the experiment at least to times for different values of D and calculate the value of the diameter of the wire.
- The same experiment can be repeated for wires of different thickness.

Observations:

The value of wavelength of the laser light is=.....m

Observation table for the value of D and beta

Sr. No.	Type of wire	Readings for value of D			Readings for value of β			Value of $\frac{D}{\beta}$	Mean value of $\frac{D}{\beta}$
		Position of wire upright a (cm)	Position of screen upright b (cm)	D = b-a	Position of (cm)				
					5 th Max. on L.H.S. x_2 (cm)	Central Max. x_1 (cm)	5 th Max. on R.H.S. x_2 (cm)	$\beta = \frac{x_2 - x_1}{5}$ (cm)	
1. 2.	A								
1. 2.	B								

- Laser source should be handled carefully and placed on the upright specially provided for it.
- The experimental wire should be thin and must be placed vertically
- For the observation of the beta, distance from central maxima must be measured.

Determination of the refractive index of different glass prisms as a function of wavelength

Aim: to find the refractive index of a prism using a spectrometer.

Apparatus required: sodium lamp arrangement, glass prism, spectrometer, reading lens

Note: for different wires, separate values of D/beta should be determined

Calculations: the value of diameter of wire $d = D/\beta$


Result

The measured value of the diameter of the wire =cm


Percentage difference =

Precautions

- The laser source should remain switched on only either during the setting of the experiment or recording the observations
- One should never directly into the laser beam with the naked eye.

The glass prism: it is a piece of transparent glass medium having triangular cross-section and is bounded by two plane refracting surfaces meeting along a common edge calling refracting edge. The angle subtended between these refracting surfaces is called angle of prism. Other binding surfaces of the prism are grounded and hence opaque. The light falling on one of the refracting surfaces suffers deviation in its path of propagation on passing through prism.


The various angles shown in the above diagram are related as


$$i + e = A + \delta$$

$$A = r_1 + r_2$$

The variation of angle of deviation with angle of incidence is shown below

δ_m is the angle of minimum deviation. The prism is set in minimum deviation position for most of the measurements. In this position we have:


Hence refractive index (μ) of the prism is defined as

$$\mu = \frac{\sin\{(A + \delta_m)/2\}}{\sin(A/2)}$$

Hence the determination of prism requires to measurements of angle of prism and angle of minimum deviation. Since depends upon the wavelength of light by the cauchy's reaction, so monochromatic source of sodium light is used.

Procedure:

- a) Setting the spectrometer
- i) Similar process of horizontal alignment is done for the prism table by adjustment of screws P, Q and R.
 - ii) The width of the slit is adjusted so as to obtain a narrow slit. This slit is viewed through telescope and distance of eye piece from the objective is varied so as to obtain sharp and bright image of narrow slit.
 - iii) Note the least count or Vernier constants of verniers V1 and V2.
- b) Determination of angle of prism by rotation of telescope method
- i) Place the prism such that its refracting edge is coincident with center of prism table and opaque face opposite to this edge is parallel to the line joining screws P and Q.
 - ii) Rotate the telescope towards left of prism to observe the shape bright image of slit after reflection of incident light from face AB of prism. Position the telescope carefully through fine adjustment so that slit is coincident with vertical cross-wire. Note the reading of both the Vernier scales. Reading of two verniers must have difference of 180° .
 - iii) The angle of rotation 1 and 2, shown as difference in readings of each Vernier at two positions, are double the angle of prism. Hence angle of prism is found to be


Vernier constant

Vernier V1			Vernier V2		
1 st Face	2 nd Face	Difference (α)	1 st Face	2 nd Face	Difference (α_0)

$$? : (\alpha + \alpha_0)/4$$

Determination of angle of minimum deviation

- i) Place the prism with its center coinciding with the center of prism table and face opposite to the refracting edge facing you
- ii) Look through the face AC and locate the position of refracted image of slit with the naked eye. Bring the telescope in this position and adjust its position till the image coincides the vertical cross-wire of the telescope.
- iii) Now turn the prism table in the direction such that telescope has to be moved towards the line of collimator axis to keep the image in field of view.
- iv) When the prism table is rotated in this way, the angle of deviation decreases. At the position of minimum deviation, the image becomes stationary. Further rotation of prism table moves the slit image in opposite direction. The position of prism table is adjusted at position that the slit just becomes stationary. Clamp the table and telescope. The position of minimum deviation is correct if the image remains stationary even on giving one rotation to tangent screw either way.

- v) Keep the table fixed and remove the prism gently. Turn the telescope in line with the collimator such the slit image is coincident with vertical cross-wire. Clamp the telescope and note readings on verniers V1 and V2.

Telescope Reading		Difference $\phi_k = \alpha - \beta$	Vernier V2		Difference $\phi_k = \alpha - \beta$
Min Dev. (α)	Direct (β)		Min Dev. (α)	Direct (β)	

$$\delta_k = (\phi_k + \phi_o)/2$$

The refractive index (μ) of the prism is given by the expression

$$\mu = \frac{\sin \{(A + \phi_k)/2\}}{\sin (A/2)}$$

Precautions and sources of error

- i) The axis of telescope, the collimator and plane of the prism table must be horizontal.
- ii) The position of the eye piece of the telescope should be so adjusted that cross-wires are visible without any strain on the eye.
- iii) The telescope should be focused for infinity and collimator should give parallel beam of light.
- iv) The slit should be narrow


- v) The prism should be placed with its refracting edge at the center of the prism table when finding the angle of prism
- vi) The prism table should be leveled so that the refracted image of the slit lies at same height from the horizontal cross-wire when observed from both the refracting faces of the prism.
- vii) Both the verniers should be read to avoid error due to the non-coincidence of the center of the circular scale with the axis of rotation of the telescope or the prism table.
- viii) To find the angle of minimum deviation, the prism table should coincide with center of prism.
- ix) The position of minimum deviation, the prism table should coincide with center of prism.
- x) The position of minimum deviation must be located carefully by fine movement of tangent
- xi) The position of minimum deviation must be located carefully by fine movement of tangent screw associated with rotating of prism table.

Studing specific rotation using polarimeter


Aim: to determine specific rotation of the sugar solution using a biquartz polarimeter

Apparatus required: biquartz polarimeter, mercury lamp arrangement, sugar, weighing balance, graduated cylinder, two beakers, funnel, a pipette and a glass rod.

- a) Biquartz polarimeter: it is an instrument used to measure specific rotation of certain solutions. the optical parts of the spectrometer are shown in the figure below.

- l) Nicol prism: it is an optical device made of calcite and is used in many instruments for producing and analyzing plane polarized light. A calcite crystal about three times its length as its width is taken. The principal section is shown below.


Its end faces are cut down so as to reduce the angles at B and D from 71 to 68 with the principal section. The crystal is now cut apart along AC perpendicular to both principal plane and end faces such that AC makes an angle of 90 with ends AB and CD. The two cut faces are ground, polished optically flat and cemented together with Canada balsam whose refractive index lies between that of calcite for ordinary and extraordinary rays. The sides of the prism are blackened to absorb totally reflected ray.

A ray incident on the face of prism is split up into the ordinary and extraordinary ray. The ordinary ray is totally reflected at interface coating (of Canada balsam) whereas the extraordinary passes through

giving rise to plane polarized light having vibrations in principal plane. The nicol prism can be used both as polarimeter as well as analyzer.

When the principal section of two prisms is parallel the extraordinary ray is transmitted through analyzer. If one of prisms is rotated the intensity of the transmitted beam decreases and finally no light is transmitted (when principal sections of two prisms are perpendicular). In this extraordinary ray of the first prism (polarizer) becomes the ordinary ray for the second prism (analyzer) and it suffers total internal reflection.

Biquartz: a biquartz is simple and sensitive device which can be used instead of half shade arrangement in a polarimeter. It consists of two semicircular quartz plates ACB and ADB each 3.75mm thick, one of left handed (L) and the other right handed (R) quartz cut with their optic axis perpendicular to their refracting faces. These are cemented together along AB to form a circular plate. When plane polarized white light is incident normally on the plate, each half of this plate will rotate each of the colour in the opposite directions. For yellow light, the rotation is about 90

If POQ is the direction of incident vibration, then for yellow light, it will OQ' in the right handed half and OP' in the left handed half. Hence vibrations of yellow light will along P'OQ' perpendicular to the direction POQ. Thus if the principal plane of analyzing nicol is parallel to POQ yellow light will be quenched in both the halves, while other colours will be present in same proportions thereby giving same tint of grayish violet colour to both the halves. This is known as the tint of passage. If the analyzer prism is rotated to one side one half of the field of view appears blue while other appears red. If rotation to analyzer is given in opposite direction, the colours of two halves are interchanged.

As the transition from blue to red is rapid. The zero position can be adjusted accurately.

1) Working of polarimeter: the light from the monochromatic source after passing through a narrow slit is rendered into a parallel beam by lens L. the light is made polarized by the nicol prism polarizer and after passing through biquartz plate, glass tube (polarimeter tube) containing the solution is to fall on the analyzer nicol prism. The light is viewed through Gaussian telescope.

The analyzing nicol prism. The light is viewed through gaussian telescope. The analyzing nicol prism can be rotated about the axis of tube and its rotation can be measured on the graduated circular scale divided in degrees, with the help of a Vernier.

Theory

Unpolarised and polarized light: the light is a transverse wave in which the oscillations of electric and magnetic fields are mutually perpendicular to

each other and to the direction of propagation. The oscillating magnetic field is very weak and plays no significant role in visual effects.

Hence an ordinary light can be supposed to consist of oscillations of electric field in all possible directions transverse to the direction of propagation. The oscillating electric vectors in various directions can be resolved into rectangular components lying in the transverse plane. The phenomenon of restricting the vibrations of light in a single direction in a plane perpendicular to the direction of propagation is called polarization of light. Such a light lacks symmetry and confirms the transverse nature of light. The plane in which the vibration takes place in the polarized light is called plane of vibration while the plane perpendicular to the plane of vibration is called plane of polarization. The polarized light are of the following types:

- a) Plane polarized light: the polarized light in which vibration of electric vector is restricted in a single direction in a plane perpendicular to the direction of propagation is called plane polarized light.
- b) Circularly polarized light: in such polarized light, the direction of the oscillation of the plane rotates periodically, but the magnitude remains constant.
- c) Elliptically polarized light: in such a light, the direction of oscillations of electric vector rotates periodically but magnitude varies within certain maximum and minimum limits.

Methods of polarization: the polarized light can be produced by the following methods

Polarizer crystals: certain organic crystals like tourmaline can be used to produce plane polarized light. When ordinary light is passed through such

crystals, only the vibrations along the optic axis of the crystal are transmitted, thereby producing plane polarized light. If this polarized light is allowed to pass through another similar crystals. Intensity is maximum when axes of two crystals are parallel while it becomes zero, when two axes perpendicular. The intensity variation is expressed as Malus law.

Polarization by reflection: if an ordinary beam of light is incident on transparent medium, then its reflected component is found to be partially polarized. The degree of polarization increases with the angle of incidence and at a particular value of angle of incidence, called polarizing angle, the reflected beam is completely polarized. At polarizing angle of incidence the polarized reflected light is perpendicular to the refracted unpolarized light and the refractive index of the transparent medium is related to the polarizing angle. This is referred to as Brewster's law.

Polarization by scattering: when an ordinary light passes through the medium, it suffers scattering from the particles of the medium. If the size of particles is comparable to the wavelength of light, the light scattered in the direction perpendicular to the direction of incidence gets plane polarized.

Double refraction: when the ordinary light passes through a calcite crystal, the refracted ray splits up into two components which are:

- a) Ordinary ray: this component obeys the law of refraction and has vibrations perpendicular to the principal section of the crystal. The velocity of ordinary rays is same in all directions. These rays get polarized such the direction of oscillation of electric vector is perpendicular to the optic axis.
- b) Extra-ordinary ray: this component does not obey the laws of refraction and has vibrations in the principal section of the crystal.

This component has different velocities, in different directions. These rays get polarized in such a way the direction of polarization is parallel to the optic axis. This phenomenon is called double refraction. These crystals are said to be positive in which $v_o > v_e$ and the opposite negative. The crystals exhibiting double refraction have property of transmitting orthogonally polarized component of light with different velocities and this property is also called birefringence.

Optical activity: is the phenomenon of rotating the plane of polarized light when it is allowed to pass through a material. If a material rotates the plane of polarized light towards left, it is called leavo-rotatory while thowse rotating the plane of polarized light towards right is called dextrorotatory. If polarizes light passes through a column of an optically active material of length L and concentration c, then angle of rotation of plane of polarized light is

$$\theta = sLc \quad s = \frac{\sigma}{Ja}$$

This is called Biots Law and s is called specific rotation, which is dependent on temperature of material and wavelength of incident light.

Procedure

- a) Preparation of 20% solution: take a clean dry beaker. Add about 20g of sugar and calculate the volume of solution to have a 20% strength as

$$\text{Volume} = (m \times 100/20)\text{cc} \quad m = \text{mass of sugar}$$

Add nearly half of the volume of water in the beaker and stir well till whole of the sugar is dissolved. Add more water if necessary taking care that volume of solution is less that the calculated value of the volume. Transfer the solution to the graduated cylinder. Rinse the beaker with small quantity of water and add to the cylinder. Complete the volume by adding water in little amount by using pipette. Filter the solution in a clean beaker and cover it.

- Setting: find the Vernier constant of the circular scale. Place the polarimeter so that its aperture is in front of mercury lamp. Look through the telescope and adjust the position of the eye piece so that two halves of the biquartz plate are grey-violet in colour and are clearly focused.
- Remove the brass caps of the polarimeter tube. Clean the tube as well as glass windows. Now replace on of the caps in position taking that no strain is exerted on the glass window. Hold the tube in a vertical position and fill it with water. Slip the second glass window gently on the tube taking care that no air bubbles are left underneath. Screw the cap gently so that no strain is exerted on the glass window.
- Place the tube in position in the polarimeter and cover it. Rotate the analyzing prism by rotating the circular scale till halves of the biquartz are equally intense grey-violet. Note the reading on the scale. Turn the analyzer through 180 again for grey-violet. Note the reading on the scale. Turn the analyzer through 180 again for grey-violet colour position and note the reading on the circular scale. These are the measurements for 0% sugar solution and rotation caused by a sugar solution is measured taking this reading as reference.

- d) Now fill the tube with 20% sugar solution place it in position and note the scale readings when two halves of biquartz plate are equally grey-violet. Repeat by turning through 180. The amount of rotation introduced in plane of polarization is given by subtracting corresponding reading for pure water.
- e) Take 30cc of 20% sugar solution and add same amount of water to reduce its strength of solution to 10%. Repeat the same procedure as in e
- f) Similarly repeat for 5% and 2.5% solutions
- g) Draw graph between θ and C and calculate the slope of graph.

Observations

Vernier constant of circular scale=

Length of polarimeter tube (in decimeters)=L=

Temperature of water=

Precautions and sources of error:

1. The tube and glass windows should be clean
2. There should be no air bubble in the glass tube when solution is filled in it.
3. The caps should be so screwed that there is no leakage and glass window are stain free. before filling the tube, it should be rinsed with the same solution
4. The reading should be taken in the equally grey violet position and not in equally blue or red positions

Concentration of Solution	Scale readings through solution		Rotation		Mean rotation $\sigma; (\theta_1 + \theta_2)/2$	$\sigma \cdot A$
	First position (α)	Second position (β)	First position $\theta_1 = (\alpha - \alpha_0)$	Second position $\theta_2 = (\beta - \beta_0)$		
0%	$\alpha_0 =$	$\beta_0 =$				
20%						
10%						
5%						
2.5%						

Mean/C =

Specific rotation = $s = \theta / LC =$

Value of specific rotation from the tables at given temperature =

Percentage error =

Verification of Malus's law

Aim: to verify Malus's law for the plane polarized light


Apparatus required: optical bench, diode laser, analyzer and polarizer with graduated circular scales, photodetector, uprights for holding the accessories.

Theory: Malus's law describes how the intensity of light transmitted through a linear polarizer varies as a function of the angle between the polarizer transmission axis and the plane of polarization of the incident light. Malus's law is

$$I(\phi) = I_0 \cos^2 \phi$$

Where I_0 is the transmitted intensity when the polarizer transmission axis and the plane of polarization of the incident light are parallel.

To verify Malu's law, you will use two polarizing sheets and a photodetector as shown in the figure. By convention, we will call the first polarizing sheet the polarizer and the second the analyzer. The polarizer ensures that the light reaching the analyzer is linearly polarized. Since the laser is already linearly polarized, the polarizer can also be used to adjust the amount of laser light used in the experiment so that the photodetector is not saturated.

Only one of the polarizing sheets is mounted in a rotation mount; use that one as the analyzer. The polarizer can simply be held in a filter holder. Measure the power transmitted through the analyzer as a function of the rotation angle of the analyzer. Measurements at 5° and 10° intervals should be sufficient, except near 0° and 90° where more points will be needed. It is useful to first find the minimum and maximum points so that

I_0 can be properly determined. Plot your results and compare them with the theory.

Procedure

1. Level the optical bench with the help of spirit level and the screws provided for the purpose.
2. Arrange the diodelaser, the polarizer, analyzer and the photodetector on the optical bench.
3. Adjust the height of the uprights so that the light from the laser after passing through the polarizer and analyzer enters the photodetector.
4. Keep the angle of the polarizer at any fixed value by adjusting the circular scale provided with it. Bring the angle of the analyzer to 0° on the circular scale of the analyzer and note the reading in the digital meter of the photodetector, which is a measure of the intensity I of the light coming from through the analyzer.
5. Keep the angle of the polarizer same, rotate the analyzer, changing its angle from 0° to 360° in steps of 10° and note the corresponding reading (I) in the digital meter.
6. Find out the first value of θ for which I is maximum. For this position of the analyzer, the planes of the transmission of the polarizer and the analyzer are parallel, i.e. the angle between the planes of the transmission of polarizer and analyzer is zero. Find the angle for the other values of θ .
7. Plot a curve between along-x-axis and I along y-axis.
8. Plot another curve between

$$\cos^2 \phi$$

along x-axis and I along y –axis. It should be a straight

line

Obsevatios

Least count of the circular scale of the polarizer

Least count of the circular scale of the analyzer

Angle of polarizer

Precautions and sources of error

1. The source of light, polarizer and analyzer should be at the same height.
2. The experiment should be performed in the dark room so that no other light reaches the photodetector
3. The voltage applied to the source of light should be constant throughout the experiment.
4. The whole of the reflected light should pass through the polarizer.
5. The whole of light emerging from the polarizer should be incident on the photodetector.


Diffraction grating experiment for the determination of wavelength of light

Aim: to determine the wavelength of sodium light using diffraction grating.

Sr. No.	Angle of analyser σ (degree)	Intensity I	Angle between the analyser and polariser ϕ (degree)	$\cos^2 \phi$
1.	0°	I_{\max}	0	1
2.	10°			
3.	20°			
4.	30°			
-	-			
-	-			
-	-			
-	-			
-	360°			

Apparatus required: a spectrometer, spirit level, diffraction grating with clamping arrangement, sodium lamp arrangement. Glass tube filled with neon gas sodium lamp arrangement: it consists of a U shaped glass tube with two tungsten electrodes in the form of spiral coated with BaO. This whole of the tube is enclosed in a vacuum jacket. The U-tube is filled with Neon gas at 100mmHg pressure and some specks of metallic sodium are deposited on inner surface of tube. To start the discharge, a potential of 400V is applied through a step up transformer.

These lamps are available at 35W, 80W and 140W. the discharge first passes through the neon gas and produces pink colour. The heat is generated which vaporizes the metallic sodium specks, which expands to fill whole of the tube volume. The passage of discharge through the sodium vapor produces persistent yellow colour light. This light source emits two close lying wavelengths of $\lambda = 5890 \text{ \AA}$ and 5896 \AA and for practical purposes considered to be monochromatic in nature. The sodium lamp is connected to the mains through a special AUTOLEAK transformer, which has high leakage reactance resulting drooping current voltage characteristics.


Optical spectrometer: the spectrometer is compact apparatus for obtaining pure spectrum. This spectrometer consists of following parts:

i) collimator: its purpose is to produce a parallel beam of light. It consists of horizontal tube on the arm of a spectrometer, which has converging achromatic lens at one end and a vertical slit at the other on a sliding tube. The focal length of the converging lens is equal to the length of the tube. The distance between the lens and slit can be varied using rack and pinion arrangement or by sliding the inner tube to obtain parallel rays. The tube rests on two screws C1 and C2, which can be slightly tilted up and down for horizontal alignment. The slit is in the form of two sharp edges of which one is flexed while other can be move parallel for changing its width and length.

ii) telescope : it is an astronomical telescope with achromatic objective and Ramsden's eye piece. It is mounted on another arm fixed rigidly to the circular scale graduated in 0.5° . the telescope along with the circular scale can be turned around the vertical axis passing through the center of the spectrometer. It can be fixed to any position through a screw T1 and can then be given slow rotation using tangent screw T2. The position of telescope can be read through two Vernier scale V1 and V2, 180° .


The telescope rests on two screws T3 and T4 for horizontal alignment. The rack and pinion arrangement is provided on the side of tube for focusing.

Prism table: it consists of upper plate a and lower plate B separated by three springs, through which leveling screws P, Q and R passed. A set of parallel equidistant lines are engraved (say between screws P and Q). the prism is always placed with one of the reflecting faces perpendicular to these lines. A series of circles concentric with vertical axis of rotation are ruled on this plate, which help in placement of the prism correctly. The height of the prism table can be adjusted by clamping the screw A1, which fixes the prism table to verniers V1 and V2. The table can rotate about the vertical axis and may be fixed at any position by means of screw A2. After fixing by A2, it can be removed slowly by tangent screw A3 at the base of the spectrometer. The position of the table can be read on two Vernier scales.

Diffraction grating: a plane diffraction grating consists of an optically plane glass plate, on which are ruled number of equidistant parallel straight lines.

These lines divide the glass plate into opacities (b = grooves created by the lines) and transparencies (a = thickness of glass plate between two opacities), the thickness of which is of the order of wavelength of light. The quantity ($a+b$) is called grating element. The number of lines in a plane transmission grating is of the order of 6000/cm. in laboratories, actual grating is not used as it is very costly. The original grating is usually ruled on a plate of speculum metal by a sharp grating are obtained by depositing a thin film of gelatin on it. This film is removed when dry and pasted on optically plane glass plate

ABC...H in the figure shown below represents the section of plane transmission grating supposed perpendicular to the plane of paper. Let the width of clear space be "a" and that of ruling be "b". the distance ($a+b$) is called grating element. The point A and C separated by grating element are called corresponding points. When a parallel beam of monochromatic light is incident normally on the grating, it suffers diffraction. In this figure ABC... H represents the incident plane wavefront while AMN...K represents the wavefront after diffraction through grating at an angle θ . There is no path

difference in reaching the plane ABC...H the only phase change that occurs is due to difference in path transversed by different rays in reaching from ABC ...H to AMN....K. the path difference between the corresponding points A and C is $(a+b)\sin \theta$. the same is true any two corresponding points. The transmitted light after falling on the lens is focused at the point P. the point P will be bright or dark according as the rays reinforce with each other. They would reinforce the give brightness or maximum if

$$(a + b) \sin \theta_n = n\lambda$$

They would interfere and produce darkness or minima if

$$(a + b) \sin \theta_n = (2n+1)\lambda/2$$

When $n= 0$ it gives central bright fringe or the first order spectrum ($n=1$)

$$(a + b) \sin \theta_1 = \lambda$$

For the second order spectrum $n= 2$

$$(a + b) \sin \theta_2 = 2\lambda$$

$$\sin \theta_2 = 2\sin \theta_1$$

The above relation is helpful in ascertaining the second order spectrum after the first as in certain cases, the second order spectrum is missing and first is followed by the third order spectrum. If the light is not monochromatic, then each wavelength will give its own primary maxima. Hence each order will consist of spectrum.

The number of lines N is marked on the grating. The value of the grating element is given by

$$(a + b) \sin \theta = n\lambda$$

If the path difference between rays from the extreme ends of one single slit is λ , then the intensity of light is zero. It is because in such a case a slit can be divided into two half slits and the path difference between the corresponding points on two half slits is $\lambda/2$. The condition for zero intensity is $a \sin \theta = \lambda$. When above given two conditions are simultaneously obeyed the beams from all the slits destructively interfere to produce zero intensity and the corresponding order of the spectrum is absent from the diffraction pattern. Thus the n th order spectrum we have the condition that $a + b/a = \lambda$

Procedure

a) Settings

- I) Adjust the position of telescope eyepiece so that the cross-wires are clearly visible. Focus the telescope on a distance object and set it for parallel rays. Level the spectrometer by the leveling screws and then the prism table with help of spirit level.
- II) Fix the grating stand on the circular table with two screws in the holes drilled on one of the lines parallel to the line joining two of the screws meant of the purpose say P and Q . The face of the stand which the clamps are attached should be at the centre of the table. Take out the grating from the box, holding it from the edge and without touching its surface, fix it carefully to the frame with the ruled surface towards the telescope.

b) Optical leveling of the grating table

- I. Rotate the table so that the plane of grating is approximately inclined at an angle of 45° to the collimator axis. Rotate the telescope to receive the light from the grating surface. If the image is not symmetrical w.r.t. the horizontal cross-wire, adjust with the help of third screw R . In the position, the plane of the grating will be vertical.
 - II. ii) Rotate the table carrying the grating so that the plane of grating is perpendicular to the axis of collimator. Look for the first order spectrum on one side of the direct image of the slit. Turn the telescope so that the first order diffracted image coincides with the vertical cross-wire. If this image is not symmetrical w.r.t. the horizontal cross-wire adjust the level of prism table with the help of screws P and Q in this position, grating lines are parallel to the axis of the spectrometer. Now turn the telescope to the other side so that vertical cross-wire again coincides the first order diffracted image. If careful adjustments are done, then diffracted image of the slit will be symmetrical to the horizontal cross-wire at all positions.
- #### c) Setting the grating normal to the incident light
- I) Place the telescope in line with the collimator so the vertical cross-wire falls exactly at the center of image of the slit. Note the scale reading. Add 90° to the reading and set the telescope at this reading to place it perpendicular to the axis of the collimator. Clamp in this position.
 - II) Rotate the grating table till the plane face of the grating is facing both the collimator and the telescope. Turn the slowly till the center of the slit fall exactly on the vertical cross-wire. In this position the plane of grating is at angle of 45° to the incident light. Note the reading. Turn the table through 45° from this

position so that the plane of the grating is normal to the incident light with its plane facing towards the collimator. The grating is now set normal to the incident light with its ruled surface away from the collimator. Clamp the table in this position.

- III) Place the eye in the front of collimator and move it gradually towards the telescope till the first order diffraction image is visible. Bring the telescope in this position and observe the image through it. Clamp the telescope in this position. If the resolution of telescope is sufficiently high, two distinct narrow lines corresponding to 589nm and 589.6nm will lie side by side in the field of view. Ordinarily two lines will appear as one in the first order spectrum. Turn the tangent screw of the telescope till the cross-wire coincides the center of the image of the slit. Note the readings on both the Vernier scales. Similarly observe the first order diffracted image and note the readings on both the verniers.
- IV) Similarly observe the second order diffracted image on either side of the image and note the readings on both the verniers.
- V) Note the number of lines marked on the grating and replace it carefully with its ruled surface upwards in the box.
- Vernier constant of spectrometer
 - Number of lines per inch on grating =N=
 - Grating element (a+b)=2.54/N in cm
 - Direct reading of the telescope

Wavelength of sodium light as calculated from first order spectrum
 $=\lambda_1 = (a+b)\sin\theta_1$

Order of Spectrum	Telescope Reading			Angle of diffraction		
	Left	Direct	Right	Left	Right	Mean
I st order						$\theta_1 =$
II nd Order						$\theta_2 =$

Wavelength of sodium light as calculated from second order spectrum
 $\lambda_2 = ((a+b)\sin\theta_2)/2$

Mean wavelength $(\lambda_1 + \lambda_2)/2$

Precautions

- Grating should be held from the edges and the ruled surface should not be touched
- The telescope should be focused on the brightest image of the slit while seeing the reflected image.
- The ruled surface should face away from the collimator.
- The light should fall on whole of the grating surface.
- The adjustment is incorrect if the angles of diffraction on the left and the right of the direct position are not equal.

In the second order spectrum with a good grating the two D lines generally appear to be separated. In such case, the angle of diffraction should be found separately for each line and wavelengths calculated corresponding to each.

Determination of the thickness of thin metallic wire using laser beam.

Aim: to determine the diameter of a thin wire using a He-Ne laser.

Apparatus required: laser light source (He-Ne or semiconductor laser), an optical bench with four uprights, micrometer eyepiece, adjustable fine slit width, piece of experimental wire, clamp for holding thin wire, screen, tracing paper and pencil.

Formula used: the diameter of the thin wire is given by the relation

$$d = \frac{\lambda D}{\beta}$$

Where d = diameter of the wire D distance between the screen and the wire. Beta fringe width i.e. separation between the centers of the two consecutive maxima.

Procedure:

Adjust the screen and the optical bench perpendicular to each other. Put the uprights on the bench.

Hold the straight thin wire tightly in its clamp and mount it on an upright on the optical bench.

Place the screen on the upright at one of the extreme end position of the optical bench.

Put a laser source on its upright about 2 to 3 meter away at another end of the optical bench.

Switch on the power supply of the laser source and illuminate the wire with the laser light. The upright having the wire should be near the laser source.

Adjust the position of the wire upright and the laser source so that a well defined diffraction pattern on both sides of central maxima is observed on the screen. The pattern may also be obtained on the wall.

Sr. No.	Type of wire	Readings for value of D			Readings for value of β			Value of $\frac{D}{\beta}$	Mean value of $\frac{D}{\beta}$
		Position of wire upright a (cm)	Position of screen upright b (cm)	D = b-a	Position of (cm)				
					5 th Max. on L.H.S. x_2 (cm)	Central Max. x_1 (cm)	5 th Max. on R.H.S. x_3 (cm)	$\beta = \frac{x_3 - x_2}{5}$ (cm)	
1.	A								
2.									
1.	B								
2.									

Place a long strip of the trace paper on the screen and with sharp pencil mark the position of the maxima on this paper.

Record the position of wire and screen uprights to calculate the value of D. record the distance between central maxima and centers of five

consecutive maxima with the help of micrometer reader to find the value of beta.

Repeat the experiment at least two times for different values of D and calculate the value of the diameter of the wire.

The same experiment can be repeated for wires of different thickness.

Observations

The value of wavelength of the laser light is =.....m

Observation table for the value of D and beta.

Note: for different wires, separate values of D/beta should be determined.

Calculations the value of diameter of wire $d = D/\beta$

Result:

The measured value of the diameter of the wire=cm

The actual value of the diameter of the wire=.....cm

Percentage difference=.....

Precautions and sources of error

1. The laser source remain switched on only either during the setting of the experiment or recording the observations.
2. One should never directly into the laser beam with the naked eye.

3. Laser source should be handled carefully and placed on the upright specially provided for it.
4. The experimental wire should be thin and must be placed vertically.
5. For the observation of the beta, distance from central maxima must be measured.