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

## Instruction Manual Optical Bench, code QLG004


### **OPTICAL BENCH**

An optical bench made of well-painted square cross sectional aluminium pipe with closed ends, fitted on two moulded supports with levelling screws. A scale graduated in millimetre fitted on the top of the optical bench. There are specially designed, easy sliding tilt proof riders with index mark for mounting various accessories, height of accessories adjustable with the help of knobs fitted in the riders. We can perform the experiments related to the fundamental of optics, such as image formation, reflection/refraction through the optical elements.

### **Contents:**

- Optical bench with holder
- Lamp holder with lamp
- Optical needle
- Cross slit
- Multiple needle with stand
- Single needle with stand
- Lens Holder
- Diaphragm
- Jolly Photometer
- Mirror

- Ground glass screen
- Screen
- Object Screen
- Grating holder


**Optical bench:**

An optical table is a vibration control platform that is used to support systems used for laser- and optics-related experiments, engineering and manufacturing. The surfaces of these tables are designed to be very rigid with minimum deflection so that the alignment of optical elements remains stable over time.

**Experiment 1:** To find the focal length of the convex lens using light source.

**Theory:**

The focal point of a lens is found by allowing a bundle of mutually parallel rays to enter the lens (i.e., from an object infinitely far from the lens). The lens alters the direction of these rays, making them emerge as a convergent or divergent bundle. The point to which they converge (or from which they diverge) is called the focal point. The diagrams illustrate this; F labels the focal points.

- A lens which converges a bundle of parallel rays is called a converging lens, or positive lens (its focal length is taken as positive.) The converging lens is thicker at its center than at its edge.
- A lens which diverges a bundle of parallel rays is called a diverging lens, or a negative lens (its focal length is taken as negative.) The diverging lens is thicker at its edge than at its center.
- Light rays from a point source (object) passing through a lens emerge convergent to a point or divergent from a point. In either case, that point is called the image of the source.
- When the emergent rays converge to a point, the image is called real.

When the emergent rays diverge from a point, the image is called virtual. Such images can be seen only by looking through the lens, toward the light source. By our convention, with the rays passing through the lens from left to right, you must have your eye to the right of the lens and look through the lens to see the image which is to the left of the lens.

The relation between  $u$ ,  $v$  and  $f$  for a convex lens is

$$\frac{1}{f} = \frac{1}{v} - \frac{1}{u}$$


Where  $f$ =focal length of convex lens

$u$  =distance of object needle from optical center of the lens (-)

$v$  = distance of image needle from optical from optical center of lens (+)

Then

$$f = \frac{uv}{u - v}$$

### Components Required:

- Convex Lens (Not Included)
- Optical Bench
- Object Rider – 3 Nos.
- Screen
- Light source with stand
- Object screen (wire gauge at center)
- Lens Holder

### Procedure:

1. Place the optical bench on an appropriate place.
2. Place 3 object rider on the optical bench.
3. Insert the light source into light source holder.
4. Situate the object screen on the first upright near to the light source.

5. Locate the receiving screen near to the object screen.
6. Now place the convex lens in holder whose focal length has to determine in between the receiving screen and object screen.
7. Blurred image of the wire gauge can be seen on the receiving screen.
8. Adjust the distance between the receiving screen and the lens to get the clear image of the wire gauge.

**Observation and Table:**

Distance between lens and object screen	Distance between receiving screen and lens (v) cm	$f = \frac{uv}{u-v} \text{ cm}$

**Experiment 2:** To find the focal length of a convex lens by displacement method.

**Components Required:**

- Optical bench with four uprights
- Single Optical needle
- Multiple Optical needle
- Convex lens (Not Included)
- Lens Holder

**Theory:**

the relation between u, v and f for a convex lens is

$$\frac{1}{f} = \frac{1}{v} - \frac{1}{u}$$

$$f = \frac{uv}{u-v}$$


where,

$f$  = focal length of convex lens

$u$  = distance of object needle from optical centre of the lens

$v$  = distance of image needle from optical centre of the lens.

**Note:** According to sign-convention,  $u$  has negative value and  $v$  has positive value. Hence,  $f$  comes positive.

**Fig.** Focal length of convex lens.

**Procedure:**

**To determine rough focal length:**

1. Mount the concave mirror in mirror holder.
2. Go out in the open and face the mirror towards distant tree or building.
3. Obtain the image of the tree or the building on a white painted wall (screen) and move the mirror forward and backward to get a sharp image on the wall.
4. Measure the distance between the mirror and the wall (screen). This will be equal to the rough focal length of the mirror.

To set the lens

5. Clamp the holder with lens in a fixed upright and keep the upright at 50 cm mark.
6. Adjust the lens such that its surface is vertical and perpendicular to the length of the optical bench.
7. Keep the upright fixed in this position throughout.

**To set the object needle:**

1. Take the thick optical needle as object needle (O). Mount it in outer laterally moveable upright near zero end.

2. Move the object needle upright and clamp it at a distance nearly 1.5 times the obtained rough focal length of the lens.
3. Adjust height of the object needle to make its tip lie on horizontal line through the optical centre of the lens.
4. Note the position of the index mark on the base of the object needle upright.

**To set the image needle:**

1. With left eye closed, see with the right open eye from the other end of the optical bench. An inverted and enlarged image of the object needle will be seen. Tip of the image must lie in the middle of the lens.
2. Mount the other optical needle (single needle) in the fourth upright near the other end of the optical bench.
3. Adjust the height of the image needle so that its tip is seen in line with the tip of the image when seen with right open eye.
4. Move the eye towards right. The tips will get separated. The image tip and the image needle tip have parallax.
5. Remove the parallax tip to tip.
6. Note the position of the index mark on base of the image needle upright.
7. Record the position of the upright of the lens, the object needle and the image needle in the table.

**Observations:**

Rough focal length of the given convex lens = .....cm

Observed distance between		$f = \frac{uv}{u - v}$
Object Needle and lens u (cm)	Image Needle and lens v (cm)	

**Experiment 3:** To find the focal length of the concave lens using convex lens by contact method.

**Theory:**

Concave lens

A concave lens is thinner at the center than at the edges. So the light beams passing through the lens are spread out or diverged. Therefore, the concave lens is called a diverging lens. The image formed by a concave lens is virtual and diminished.


Since a concave lens will not produce a real image, a convex lens is used to measure its focal length. There are two methods to find the focal length of a concave lens.

**Lens in contact method**

When a concave lens of focal length  $f_2$  is placed coaxially in contact with a convex lens of focal length  $f_1$  and if  $F$  is the focal length of the combination, then,

$$\frac{1}{F} = \frac{1}{f_1} + \frac{1}{f_2}$$
$$f_2 = \frac{Ff_1}{f_1 - F}$$

Therefore, the focal length of the concave lens,

### Components Required:

- Optical Bench
- Object Riders
- Lens Holder
- Screen
- light source
- Light source Holder
- Object Screen
- Convex Lens (Not Included)
- Concave Lens (Not Included)

### Procedure:

1. Place the optical bench with object riders on an appropriate place.
2. Insert the light source in light source holder.
3. Now situate the object screen near to the light source.
4. Insert the convex lens and concave lens in contact with each other in one holder (convex lens towards light source).
5. It is placed between the object screen and the receiving screen at a fixed distance away from the light source. It is taken as 'u'.
6. The position of the receiving screen is adjusted to get a clear image of the wire gauge.

7. Measure the distance between the combination of lenses and the receiving screen. It is taken as 'v'.
8. Calculate the focal length of the combination as mentioned in table.
9. Repeat the experiment by placing the combination of lenses at different distances away from the hole slit.
10. Then calculate the mean value of F.
11. Using the values of  $f_1$  and F, the focal length of the given concave lens can be calculated.

**Observations:**

S.No.	Distance between the combination of lenses and		Focal Length $F = \frac{uv}{u - v}$
	Object u (cm)	Image v (cm)	
1			
2			
3			

- Focal length of the combination lens, Mean F = .....cm
- Focal length of convex lens,  $f_1 = \dots\dots\dots$ cm
- Therefore, the focal length of the given concave lens,

$$f_2 = \frac{Ff_1}{f_1 - F}$$

**Experiment 4:** To demonstrate the additive effect of different color filters.

**Primary Colour:**

Red, Green, Blue are primary colour. When red, green and blue color are mixed or added together with the proper intensity, white (W) color is obtained.


The mixing together (or addition) of two or three of these three primary colors with varying degrees of intensity can produce a wide range of other colors.

**Secondary Color:**

- Red color and green color add together to produce yellow (Y) color
- Red color and blue color add together to produce magenta (M) color
- Green color and blue color add together to produce cyan (C) color

Yellow, Magenta and Cyan are called secondary colors. They can be produced by the addition of equal intensities of two primary colors.

### Additive and subtractive effect of different Colours:

### Complementary Colors:

Any two colors that when mixed together in equal intensities produce white are said to be **complementary colors** of each other.

### Complementary Colors:

- Red and Cyan
- Green and Magenta
- Blue and Yellow

### Components Required:

- Light source
- Receiving Screen
- Optical Bench
- Object Rider – 3
- Slit/grating Holder
- Diaphragm
- Color Filter Set (Not Included)

### Procedure:


1. Place the optical bench on an appropriate place.
2. Insert the object rider into the optical bench.

3. Place the light source with stand at the starting end of the optical bench.
4. Now situate the diaphragm at the first object rider in front of the light source.
5. Locate the slit holder next to the diaphragm.
6. Insert the two-color filter in the slit holder.
7. Switch on the light source and open the diaphragm hole.
8. The resultant color of the two different color filter obtained on the receiving screen as describe in the table below.

Color slits in slit Holder		Color on Screen
Red	Blue	
Red	Green	
Green	Blue	

**Inverse square law:**

The inverse square law describes the intensity of light at different distances from a light source. Every light source is different, but the intensity changes in the same way. The intensity of light is inversely proportional to the square of the distance. This means that as the distance from a light source increases, the intensity of light is equal to a value multiplied by  $1/d^2$ . The proportional symbol,  $\propto$ , is used to show how these relate. The relationship between the intensity of light at different distances from the same light source can be found by dividing one from the other. Visible light is part of the electromagnetic spectrum, and the inverse square law is true for any other waves or rays on that spectrum, for example, radio waves, microwaves, infrared and ultraviolet light, x rays, and gamma rays. The intensity of visible light is measured in *candela* units, while the intensity of other waves is measured in Watts per meter squared ( $W/m^2$ ).

### Experiment 5: Differentiate between the intensities of two light source.

#### Component required

- Optical bench
- Light Source
- Light source (Candle or bulb) (Not included)
- Jolly's photometer

#### Procedure:

1. Place the optical bench in an appropriate place.
2. Mount the light source on the first end of the optical bench.
3. Now situate the jolly's photometer with stand on the optical bench. (Light source should face the transparent sheet of the photometer)
4. Choose another source of light (Candle /bulb).
5. Place the both source at the equal distance from the photometer and observe the light intensity obtained on the screen from the eye piece.

Move the photometer back and forth between the two light sources until the two polythene blocks are equally bright.

If we know the intensity of one light source, the other can be obtained using the formula

$$I_2 = I_1 \frac{d_2^2}{d_1^2}$$

#### Precaution:

1. Each side is illuminated by only one of the sources. The tin foil prevents light from passing between the blocks.
2. The light sources and the photometer should be at the same height.

3. The photometer should be positioned on the straight line between the two sources.
4. The two light sources should be the *only* sources of light in the room. No bright sunlight!
5. To ensure uniformity of illumination at the photometer, the distance from the photometer to the nearest light source should be large compared to the size of the photometer.

**Note:** This optical bench can also be used in numerous other experiments like wavelength of the laser light, diffraction of the light, interference using the laser source with this accessory.