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17. ACCESSORIES

- Tubing
- Reservoir bottle (2 l)
- Tubing connection with NS 29/32 for connection of distillate drain tube and reservoir bottle
- Stand for funnel
- Funnel (diameter 75 mm)

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans, but a non adequate disposal would damage the environment.


The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

After cleaning

Let the apparatus run in distilling mode for one hour. Discard the distillate.

To clean the cooling spiral (15)

After the apparatus has cooled down completely, drain the water out of the evaporator through the outlet valve (9). The outlet valve (9) must remain open.
 Hang the ends of the outlet valve (9) and the overflow tubing (8) into the vessel for disposal of the acid. For safety reasons, place the vessel in a receive tank
 Screw the screw cap (5) tightly onto the pressure compensation opening of the level tube (6).
 Disconnect the water supply tubing from the tap and insert the funnel into the free end of this tubing. Check that the funnel fits tightly.
 Support the funnel into the funnel stand (see fig. on this page). The bottom of the funnel cone must be level with the degassing outlet (13) (see fig. on this page)
 Slowly pour in approximately 300 ml of acetic acid, 10 - 50 %, using the graduated jug. Be carefully and avoid spilling.

CAUTION!
 Acetic acid flows out outlet valve and overflow tubing into the vessel.
 - As soon as the deposits have dissolved pour 300 ml of tap water into the funnel
 - After the liquid has flowed through hang the outlet valve and the overflow tubing back into the sink.

- To rinse the cooling spiral
- Slowly pour approximately one litre of tap water into the funnel.
 - As soon as the water has flowed through close outlet valve (9)
 - Remove the funnel, and reconnect the water supply tubing to the water tap.
 - Unscrew the screw cap (5), and keep it on the base of the apparatus, behind the level tube (6).

After cleaning

Let the apparatus run in distilling mode for one hour. Discard the distillate.

15. RETURN OF THE APPARATUS FOR THE REPAIR

Carefully clean the instrument, and send it to your supplier. Please describe your complaint in detail.

16. DISASSEMBLY AND REASSEMBLY OF THE GLASS PART

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Spanish 2-10
 English 11-19
 French 20-28

INDEX OF CONTENTS

1. Description of the apparatus.....11
 2. Technical data12
 3. General Instructions and safety regulations.....12
 4. Connection to the Electrical Mains/Water Supply.....12
 5. Picture of the apparatus.....13
 6. Unpacking the apparatus13
 7. Setting up the apparatus.....14
 8. Installing of the Reservoir Bottle14
 9. Protective cap.....14
 10. Connecting of the Tubing.....14
 Connecting of the cooling water supply14
 Connecting of the water drainage14
 Connecting of the outlet valve.....14
 11. Putting into Operation.....15
 Preparing to operation15
 To adjust the water supply15
 12. Interruption of water supply during working.....16
 13. To switch off.....16
 14. Cleaning16
 To clean the evaporator.....17
 To clean the cooling spiral18
 15. Return of the apparatus for the Repair.....18
 16. Disassembly and reassembly of the glass part18
 17. Accessories19

1. DESCRIPTION OF THE APPARATUS

Distillation apparatus for producing of distilled water TYPE DP 4000 is identified for preparing of smaller amounts of distilled water.
 By way of this apparatus, it is possible to prepare approximately 4 liter of distilled water at one hour.
 Apparatus is put together from two main parts - glass part and electrical part. Both parts are connected with plastic flange
 In upper part of glass part is cooling coil. This upper part is partitioned off the lower part with glass partition. Boiling of tap water proceeds in the lower part of the glass part. Vapour comes through the tube which is connected to the centre of partition to upper part where is cooled and distilled water is drained of the distillation apparatus through the nozzle which is situated above glass partition.

2. TECHNICAL DATA

Power consumption (W)	3100
Rated current (A)	15
Electrical protection – tube fuse 5 x 20 F 8 A (A)	1
Mains voltage (V/Hz)	230/50
Distilled output	
Conductivity of distillate (l/h)	~4
Referred to 25 °C (µS/cm)	1,5 - 2,5
Distilled temperature (°C)	95
Cooling water consumption (l/h)	~50
Dimensions of the apparatus (mm)	190x190x600

3. GENERAL INSTRUCTIONS AND SAFETY REGULATIONS

- To avoid possible hazards pay good attention to the general instructions and safety regulations.
- When you operate with this apparatus wear protective clothing, goggles, and protective gloves.
- Read and follow Operating manual carefully step by step.
- Every user must be acquainted with this Operating manual and keep it within reach at all times.
- The area in which the apparatus is to be used must conform to the appropriate electrical regulations.
- Connect the apparatus only to the earthen plugs and sockets.
- Set up the apparatus on solid, flat surface only.
- Use the apparatus in a dry environment only; do not use it in a potentially explosive atmosphere.
- The available voltage and type of current must correspond to the data given on the type plate.
- The electrical connection must have front end protection with fuse of at least 15 A.
- Connections to the electrical mains and water supply must be made in such way as to avoid creating dripping hazards
- Repairs within the base of apparatus (electric circuits, heating) must be made only by authorised electricians
- Never use force.
- Use original manufacturer's accessories and original manufacturer's spare parts only.
- Use hoses delivered by the manufacturer only. Do not use old hoses.
- Never send the apparatus to your supplier for repair without having cleaned it first.

CAUTION!

DO NOT OPERATE THE DISTILLATION APPARATUS UNLESS IT IS PROPERLY CONNECTED TO AN EARTHEN ELECTRICAL INSTALLATION WHICH IS IN PERFECT CONDITION. MAKE SURE THAT EARTH CONNECTION COMPLIES WITH VALID ELECTROTECHNIC REGULATIONS. OTHERWISE, SERIOUS HAZARD OF PERSONS AND FIRE ENDANGER CANNOT BE EXCLUDED.

4. CONNECTION TO THE ELECTRICAL MAINS/WATER SUPPLY

Electrical mains	230 volts 50 Hz
Required front end protection	min. 15 A máx. 16 A
Water supply	minimal pressure 200 kPa maximal pressure 600 kPa
Fluctuation of pressure	10 % of water pressure max. In case of a higher rate of fluctuations, we recommend using a commercially available pressure regulator.

WARNING!

Disconnect apparatus from the mains by pulling plug out of the socket before cleaning. There is risk of corrosive burns. Wear protecting clothing, gloves and goggles during cleaning. Dispose of the acetic acid in accordance with approved procedures. For your personal safety we recommend disconnecting the mains supply before cleaning. Keep the reservoir bottle (2) connected. Allow the water still to cool down completely. We are recommending mark the apparatus with table for instance: "Do not switch on! Cleaning with acid".

To clean the evaporator


- Drain the water out of the evaporator through the outlet valve (9) after the apparatus has cooled down completely
- Close the valve again
- Hang tubing of the outlet valve (9) into the vessel for the disposal of the acetic acid. For safety reasons, place the vessel in a receiving tank.
- Screw the screw cap (5) tightly onto the pressure compensation opening of the level tube. Insert the funnel into the free end of the tubing leading from the overflow tubing (8). If necessary, dip the end of the tubing into hot water (approximately 80°C) for about 3 minutes before fitting. Check that funnel fits tightly.
- Support the funnel in the stand. The bottom of the funnel cone must be level with the distillate drain tube (4) (view picture on this page)
- Slowly pour in approximately 1 liter of acetic acid 10 - 50 %, using the graduated jug. Fill the funnel up to no more than one third of its capacity, to avoid spilling. Fill the evaporator up to no higher than the lower edge of the of the distillate drain tube (4). To regulate the filling rate: Place the funnel into a higher position to expedite the filling of the evaporator; return the funnel to a lower position to control the rate of addition of acid as it approaches the lower edge of the distillate drain tube, which point must not be exceeded.

NOTE!

Take care that no acid drips into the reservoir bottle. If acid drip into the reservoir bottle, carefully rinse the bottle and the bottle top.

- As soon as the deposits have dissolved drain the acetic acid out of the evaporator through the outlet valve (9) to the vessel.
- Close the outlet valve (9).
- Hang outlet valve tubing back into the sink.

To rinse the evaporator

- Slowly pour approximately one liter of tap water into the funnel
- Carefully fill the evaporator up to the lower edge of the distillate drain tube (4)
- Drain it through the outlet valve (9) and close the outlet valve again
- Repeat the rinsing three times
- Remove the funnel, and hang the overflow tubing into the sink
- Unscrew the cap (5), and keep it on the base of the apparatus, behind the level tube

12. INTERRUPTION OF WATER SUPPLY

If there is insufficient supplying of water or if the water flow is interrupted, steam will start escaping from the degassing outlet (13). The water level in the evaporator will drop. After approximately 5 minutes the heater will be automatically switch off.

Caution

- Close the water tap.
- Switch off the apparatus (mains switch), and allow it to cool down for at least 15 minutes.
- Clear the defect in the water supply.
- Switch off the temperature limiter: take a piece of non-conducting material (length minimal 50 and maximal outside diameter 3 mm), insert it to the hole (14) on the back side of the base of the apparatus and slightly press on it. Remove used piece of non-conducting material from the hole (14).
- To continue distilling repeat the procedure starting at XI. Putting into operation.

13. TO SWITCH OFF

- Press the main switch. The green lamp will go out.
- Close the water tap.
- Pull out mains supply plug out of the mains socket.

CAUTION!

The water still must have cooled down completely before the water can be drained through the outlet valve (9). Close outlet valve afterwards.

14. CLEANING

Deposits (e.g. lime) build up during the distillation process on the interior parts of evaporator (11). This gradually reduces the quality of the distillate.

When there are considerable deposits in the evaporator, or when the quality of the distillate no longer meets the requirements for the intended application, the evaporator must be cleaned.


When there are heavy deposits in the cooling spiral, they must be cleaned up latter; these deposits however do not affect the quality of the distillate.

Producer cannot accept any warranty claims if the heater is damaged as a result of overheating because the apparatus has not been sufficiently cleaned.

The following items are required for cleaning:

- Protective clothing, gloves and goggles
- A suitable funnel
- A funnel stand, minimal length 70 cm
- A graduated jug of minimal capacity 2 liters
- A suitable vessel for the disposal of the acetic acid
- Acetic acid, 10 - 50 %
- Tap water for rinsing

5. PICTURE OF THE APPARATUS

FRONT VIEW

BACK VIEW

- 1 Glass stopcock of the reservoir bottle
- 2 Reservoir bottle
- 3 Tubing connection
- 4 Distillate nozzle
- 5 Protective cap
- 6 Level tube
- 7 Tubing connection for cooling water supply
- 8 Overflow tubing connection

- 9 Outlet valve
- 10 Mains switch
- 11 Evaporator
- 12 Green lamp mains
- 13 Degassing outlet
- 14 Hole for switch off the temperature limiter
- 15 Cooling coil

ATTENTION!

Parts 1,2 and 3 are not a matter of standard accessories. There is possible to order them separately as a special accessories.

6. UNPACKING THE APPARATUS

- Open carton box
- Carefully remove upper part of interior cover made from polyurethane foam and lift up the apparatus from the box.
- Check the apparatus immediately for glass breakage and other damage which may have occurred during transportation.
- To each package of complete unit is attached Allen key with wrench widths 3 mm and hoses see below:

Material	Outside diameter [mm]	Inside diameter [mm]	Wall thickness [mm]	Length [m]
PVC	19	14	2,5	3,5
Silicon	16	12	2	0,4

- When you want to bring apparatus to other place never hold it by side tube (6), hold only body of condenser and the base.

7. SETTING UP THE APPARATUS

- The apparatus set up always on a solid flat surface, have a respect to the fact that under the apparatus can temperature reach value 50°C
- The apparatus set in a dry environment.
- The lower edge of the apparatus must be level with, above edge of the sink.
- The minimum safety distance from the supply/water drain should be 1 m, to avoid the risk of splashing.
- Never set up the apparatus closer than this.

8. INSTALLING OF THE RESERVOIR BOTTLE

- In order to achieve conductivity of distillate water mentioned in „Technical data“ apparatus and the Reservoir Bottle must form closed system. Carbon dioxide from the air dissolves in water and would cause an increase in conductivity. A suitable reservoir bottle and bottle top are available as an original accessory (see XV.Ordering data for accessories)
 - Fit glass stopcock (1) to the hole near the bottom of the reservoir bottle (2).
 - Place reservoir bottle to the right of the apparatus.
 - Outlet of the stopcock must be attainable for taking of distilled water.
- Fit the tubing connection (3) to the head of the bottle.
- Fit the supplied silicone tubing to the tubing connection (3)
- Connect free end of silicone tubing to the tubing connection of the distillate drain tube (4).

9. PROTECTIVE CAP

Unscrew the protective cap (5) and place it for storage on the base of the unit behind the level tube (6).

10. CONNECTING OF THE TUBING

Use the PVC tubing supplied for the connections.

⚠ ATTENTION!
DIP THE END OF THE PVC TUBING INTO HOT WATER (APPROXIMATELY 80 °C) FOR ABOUT 3 MINUTES BEFORE USE. THIS WILL MAKE THE TUBING FLEXIBLE SO THAT IT CAN EASILY FITTED ON THE GLASS TUBING CONNECTION.

Connecting of the cooling water supply

- Fit the end of the tubing supplied (outer diameter 19 mm) to the bent tubing connection (7) of the still body.
- Lead the tubing to the water tap and in a case of necessity cut it off to suitable length.
- Fit the second end of the tubing to the water tap.

Connecting of the water drainage

Water from the cooler (15) is led from the cooler to the tubing (6), which is also supplier of the water to the evaporator. Surplus water is led by the bent tubing connection (8) to the sink. The water level in the evaporator is given by position of the bent tubing connection (8) on the tubing (6).

- Fit the end of the PVC tubing (outer diameter 19 mm) supplied to overflow tubing connection (8)

- The second end lead to the sink and in a case of necessity cut it off to suitable length.
- Insert the second end to the sink and safe the tube against falling out off the sink
- Don't allow the tubing sag.

Attention: When you are fitting the end of PVC tubing on the overflow tubing connection (8) support the level tube (6) with your hand to avoid glass breakage.

To connect the outlet valve (9)

- Fit one end of PVC tubing (outer diameter 14 mm) to the outlet valve (9) cut the second end off at suitable length and insert this end to the sink and safe it against falling out.

11. PUTTING INTO OPERATION

Preparing to the operation

Make sure

- Make sure that main switch (10) is in the „off“ position (green lamp „mains“ do not lighting).
- The screw cap (5) is unscrewed.
- Stopcock is closed
- The stopcock of the reservoir bottle and the outlet valve (9) are closed.
- Plug the mains supply plug into the mains socket

To adjust the water supply

- Slowly turn on water tap to allow just a small amount of water to flow through the overflow (rate of water flow approximately 0, 5 to 1 l/min).
- As soon as water in the evaporator (11) almost covers the heating spiral: press the mains switch (10). The green lamp „Mains“ (12) will light up.
- As soon as the water in the evaporator starts boiling, slowly reduce the water supply at the water tap until steam starts to escape from the degassing outlet (13) at the top of still body. Do not close water tap.
- Slowly increase the water supply until it is just sufficient to stop the steam from escaping
- The water supply is thus optimally adjusted

Note:

When operating the apparatus for the first time, let it run distilling in distilling mode for one hour. Discard the distillate.

ATTENTION!

1. DO NOT TOUCH THE STILL BODY. DANGER OF BURNS.

2. RISK OF ELECTRIC SHOCK

3. KEEP OF ACCIDENTAL STROKE OF THE HARD MATERIAL TO GLASS PART OF DISTILLING APPARATUS. THERE IS HAZARD OF BREAKAGE AND IN THE CASE WHEN IT OCCURS IN TIME OF OPERATING OF THE APPARATUS OR SHORTLY AFTER IT'S SWITCHING OF, THERE IS RISK OF SCALDING TOO.