ESTERILIZADOR INFRARROJO INFRARED STERILIZER STÉRILISATEUR INFRAROUGE

Ref. BRG016, BRG025

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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TABLE OF CONTENTS

1. Introduction	4
2. Technical specification	4
3. Operation	5

1. INTRODUCTION

Nahita infrared sterilizers are designed for the rapid and effective sterilization of small laboratory instruments by means of infrared radiation. They reach high temperatures within a few seconds, eliminating microorganisms without the need for chemical agents or prolonged exposure times.

2. TECHNICAL SPECIFICATIONS

Code	BRG016	BRG025
Maximum internal temperature	825 °C	
Heating temperature	Two options: 400 y 825°C	
Maximum object diameter	14 mm	35 mm
Heating zone length	150 mm	100 mm
Preheating time	10-20 min	
Weight	1.3 kg	1.4 kg
Power	150 W	240 W
Dimensions	230x130x185 mm	190x130x208 mm
Voltage	AC220V 50Hz	

3. OPERATION

■ Operating procedure:

- a. Place the sterilizer on a stable and level work surface and ensure that the power switch is in the OFF position.
- b. Connect the unit to the main power supply.
- c. Turn on the power switch, selecting the required operating temperature mode: 400 °C or 825 °C.

■ Operating environment:

- Indoor use only
- Altitude: ≤ 2000 m
- Operating temperature range: +5 °C a +40 °C
- Relative humidity: $\leq 80\%$
- Free from vibrations and air currents that may affect performance.
- No conductive dust, explosive gases, or corrosive vapours in the surrounding air.
- Maintain a minimum safety clearance of 10 cm around the instrument during operation

■ Safety precautions:

- Do not connect or disconnect the power plug, or operate the power switch, with wet hands.
- Do not unplug the power cord while the instrument is in operation.
- Do not perform maintenance or cleaning while the instrument is switched on.
- Do not install the instrument on an uneven, unstable, or vibrating work surface.