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Instructions on environment protection

- At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances.
- It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.
- The materials are recyclable as mentioned in its marking.
- By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.
- Please inquire at the community administration for the authorized disposal location.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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4. Description

Simple construction permits the use of water or even air for demonstration of the concepts involved with the hydraulic press. Comprises a pair of graduated glass syringes of capacities 50 and 20 ml, both having finely ground pistons and provide a cross sectional area ratio 3:1 respectively. Both the syringes mounted on a stable, non-skid sheet metal base with loading platform at their top and their nozzles linked to each other through a 3-position valve, by means of which they may be opened to the atmosphere, isolated, or interconnected as desired. Both the pistons provided with limit stop to prevent their complete ejection from the barrels. Thick wall of syringes makes it capable of withstanding pressures involved. The graduated syringe barrels also permit simple Boyle's Law experiments to be performed.


Fig. 1


Fig. 2


Fig. 3

5. Procedure

- (1) Adjust the three-way valve as shown in fig. (3) and press both the syringes slowly to remove out whole air from the syringes.
- (2) Connect a rubber tube in the nozzle fitted with the three-way valve.
- (3) Deep another end of the rubber tube inside the beaker filled with colored water and turn the knob as shown in fig. (1).
- (4) Pull the small syringe piston to fill the syringe with water and turn the knob as in fig. (2).
- (5) Press the small syringe slowly. The water from the small syringe transmits to the large syringe via conduit aluminum block. The piston of large syringe will forcefully lift up even on placing a weight on the base plate.
- (6) Repeat the steps 3-5 again.

6. Precautions

- (1) Never remove the pistons more than the limited heights from the syringe.
- (2) Always apply slowly pressure in the syringes.

Mathematically, if force F1 is applied to a piston of surface area A1, then increase in the pressure that is transmitted to the other piston through the enclosed fluid is

$$P = \frac{F1}{A1}$$

This pressure applies an upward force to the other piston having surface area A2 , whose magnitude is given by P x A2 Therefore, under ideal conditions (assuming there is no loss of efficiency due to leakage of fluid) the piston is capable of supporting force F2, as given by

$$F2 = P \times A2 = F1 \frac{A2}{A1}$$

As is evident from the above equation, the applied force is increased by a factor of inverse ratio of their respective surface areas and is called the mechanical advantage of the system. Such an arrangement finds application in hydraulic jacks used to lift motor vehicles.

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

1. Purpose

It is used to demonstrate the transmissibility of the fluid pressure and how a mechanical advantage may be obtained by using interconnected cylinders of different diameters.

2. Other items required

- (1) Rubber tube
- (2) Colored water
- (3) Beaker

3. Theory

The pressure is inversely proportional to the surface area on which it is acting. Thus, if some force is applied to the smaller pistons, then it exerts higher magnitude of force on the bigger piston and is capable of lifting higher load than as compared to the loading force on smaller piston.