DESTILADOR DE VIDRIO 4 L/H, QUARTZ
GLASS WATER STILL 4 L/H, QUARTZ
DISTILLATEUR EN VERRE 4 L/H, QUARTZ

The logo for Nahita Blue, featuring a stylized square icon to the left of the brand name 'nahita' in a bold, lowercase font, with 'blue' in a smaller, lowercase font directly below it.

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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Technical Assistance Department.

Cleaning

- ◆ You should clean the distiller regularly. So, in this way, the results will be better and you can get a distiller more durable.
- ◆ The cleaning periodicity depends on the salts which water supply contains.
- ◆ You should clean the filters of the tubes whenever you clean the still. Steps to follow:
 - Put acetic acid or hydrochloric acid into the boiler through cleaning mouth (2.5). The water should be hot.
 - Wait for the attack to the residues.
 - Rinse with clean water once or twice.
 - Fill with water until you see the water level is on the heater.

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6. ACCESSORIES

There are accessories available for substitute or replacement components and for complement the device applications with the following references:

<i>Code</i>	<i>Description</i>
99002110	Boiler
99002111	Condenser
99002112	Level
99002113	Heater
99002114	Flow limiter

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non-adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Water distiller Quartz. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Water distiller Quartz, so it should be available for all the users of this equipment.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ You should carefully handle the water distiller avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the water distiller to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the water distiller.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the water distiller warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the water distiller warranty.
- ◆ Neither fuses (1A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.

Manufacturer reserves the right to modify or improve the manual or equipment.

water goes out through the outlet (2.4)

- ◆ If water steam went out through steam hole (2.2) you would have to increase water supply.

Security

- ◆ The water distiller must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the water distiller in a horizontal plane stable table, having a safety area of at least 30 cm per side.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ Do not place the water distiller near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the water distiller for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the water distiller by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair the water distiller by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the water distiller (burns, hurts...) or damages in nearby equipment.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita

- ◆ Please put the water distiller on top of a horizontal, plane and stable table making a free space at least at 30 cm per side.
- ◆ Follow all required applicable laws about safety and health (electrical security devices, inflammable, toxic and/or pathological products, etc.).
- ◆ Do not put the water distiller near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth wire and to the other end to the water distiller connector.
- ◆ Electrical wiring must be, at least, of 2,5 mm² section.
- ◆ Power point must let intensity of 11 A pass.
- ◆ The tap water should have screw mouth of $\frac{3}{4}$ ".

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220V 50Hz $\pm 10\%$.

If you are not using the water distiller for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Assembly

- ◆ Be sure that heater elements (washer of silicone rubber, guide fitting and PBTP nut) are in the order and position that the fig.1 shows.
- ◆ Insert the heater through the screw mouth (see fig.2): hold the heater of the fitting (1.3) and turn the PBTP nut (1.2).
- ◆ Place the condenser onto the vertical tube. Be sure that the O-rings are on the surface of the tube and connect the silicon tubes.
- ◆ Put the water inlet tube on the inlet still (2.6).
- ◆ Be sure of putting the filter, of tightening and connecting the other end to the tap.
- ◆ Once you have assembled check the silicone tubes hold suitable position to let water pass and to avoid water pressures.

Switch on

- ◆ Before switch on connect the end wiring of the heater to the socket, which is on the metallic base (2.12) and turn the tap on. Turn off the tap when you see that the water level is on the heater.
- ◆ Switch ON/OFF (2.10). This will light showing that the distiller is connecting to the power point. The warning light will also light (2.9) showing that the water supply is not enough.
- ◆ Turn the tap on slowly until the pilot (2.9) turns off.
- ◆ Before few minutes, the heater will be getting red and boiling will start. The stilled

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish	2-8
English.....	9-15
French.....	16-22

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT.....	10
2. DESCRIPTION	11
3. TECHNICAL SPECIFICATIONS	12
4. INSTALLATION / SETTING UP.....	12
5. MAINTENANCE AND CLEANING.....	14
6. ACCESSORIES.....	15

1. USES OF THE INSTRUMENT

This is the Nahita 4 l/h water distiller with automatic and continuous production that provides a high quality distillation. The results are high purity, free of metal ions and pyrogenic substances. Quartz resistance manufactured in high quality borosilicate glass allowing continued use of the equipment and ensures its durability.


It provides safety mechanism that automatically switches off the quartz resistance when the cooling water pressure is reduced and reconnect when recover.

Include with an adapted bracket and small dimensions that allow easy handling and cleaning


The device is designed for clinical medicine, hospitals, researching laboratories, agro-alimentary or pharmaceutical industries and schools or universities.

2. DESCRIPTION

- 1.1 Junction
- 1.2 PBTP nut
- 1.3 Ferrule
- 1.4 Guide fitting
- 1.5 End wiring of the heater

- 2.1 PBTP nut
- 2.2 Steam hole
- 2.3 Condenser
- 2.4 Water outlet
- 2.5 Cleaning mouth
- 2.6 Water inlet
- 2.7 Supporting flange
- 2.8 Boiler
- 2.9 Water level warning light
- 2.10 Switch ON/OFF
- 2.11 Metal base
- 2.12 Socket for the end wiring of the heater

3. TECHNICAL SPECIFICATIONS

<i>Code</i>	59002114
<i>Output</i>	4 l/h
<i>Conductivity at 20 °C</i>	2 µS/cm
<i>Water consumption</i>	1,0 liters/min.
<i>Consumption</i>	1 x 2500 W
<i>Power supply</i>	230 V / 50-60Hz / 11 A
<i>Dimensions</i>	55 x 15 x 50 cm
<i>Weight</i>	4 Kg

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the water stiller, take off the involving plastic and take off the poliespan protection in which it comes fitted.
- ◆ Please make sure that it does not present any damage because of the shipment. If it does please communicate to your transport agent or suministrator immediately, so that they can make the claims in the correct time limit.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original package and with all accessories and documents supplied

Please check that all the accessories are enclosed with the equipment:

- Quartz heater
- Filter
- Feed tube (3/4")
- Support for water still
- Set of silicone tubes
- User manual

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT