**CENTRÍFUGA MINI RS-1008**  
**RS-1008 MINI CENTRIFUGE**  
**MINI-CENTRIFUGEUSE RS-1008**

**RS Lab**

Ref. GDC007

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**7.2.2. Cleaning PC tubes**

PC materials are low in chemical resistance against alkaline solutions. Avoid using neutral detergents with pH higher than 9. Note that pH of some neutral detergents are still higher than 9 even if diluted according to the instruction in the maker's catalog. Use detergent with its pH between 7 and 9.

**7.2.3. Autoclaving PA, PC and PP tubes**

PA begins softening at about 120°C; PC and PP at about 130°C. Autoclave PA tubes at 115°C (0.7kg/cm<sup>2</sup>) for 30 minutes and PC and PP tubes at 121°C (1.0kg/cm<sup>2</sup>) for 20 minutes. If these temperatures are exceeded, the tubes may be deformed.

When using a sterilizing chamber, please operate as follows:

- Place tubes in vertical position, mouths upward. If tubes are placed sideways, they may deform into an oval shape due to gravity.
- Remove tube caps to prevent from deformation or rupture.
- Wait until the sterilizing chamber cools down to the room temperature before the tubes are removed.

**7.2.4. Condition and life expectancy of tubes**

The life expectancy of plastic tubes depends on the characteristics of samples, speed of the rotor used, temperature, and so on. When the plastic tubes are used for centrifugation of aqueous samples (pH between 5 and 9), their life expectancies are defined as follows:

Operation at maximum speed:

High quality tubes (PA, PC, PP): 30-50 operations

Ordinary tubes(PA, PC, PP): around 10 operations

Life expectancy of tubes also depends on the pretreatment conditions such as cleaning and sterilization, lifetime can be cut down.

Note: Do not use damaged or cracked tubes.

**8. After-sales Service**

In order to ensure centrifuge operation safely and efficiently, it is necessary for regular maintenance. If centrifuge has problems do not attempt to repair it by yourself, contact your distributor and request to send it to the Technical Service.

**Environmental protection instructions**

Do not dispose of this equipment with normal household waste at the end of its life cycle; take it to a collection point for the recycling of electrical and electronic equipment. It does not contain elements that are hazardous or toxic to humans but improper disposal will harm the environment. The materials are recyclable as indicated by the marking. By recycling materials or other forms of reuse of old equipment, you are making an important contribution to the protection of the environment. Please contact your local authority for advice on the collection points.

**INDEX OF LANGUAGES**

Spanish .....	2-9
English .....	10-17
French .....	18-25

**CONTENT INDEX**

Safety reminder.....	10	5. Operation.....	14
1. Specifications .....	12	6. Maintenance.....	14
2. Required operational condition .....	12	7. Instructions for rotors and tubes .....	15
3. Installation .....	12	8. After-sales service .....	17
4. Rotor preparation.....	13		

**Safety reminder**

Common safety precautions

- Follow the instructions and procedures described in this manual to operate this centrifuge safely.
- Carefully read all safety messages in this manual and the safety instructions on the instrument.
- Safety messages are labeled as indicated below. They are in combination with signal words of "WARNING" and "CAUTION" with the safety alert symbol to call your attention to items or operations that could be dangerous to you or other persons using this instrument. The definitions of signal words are as follows:

**WARNING:** Personal Danger

Warning notes indicate any condition or practice, which if not strictly observed, could result in personal injury or possible death.

**CAUTION:** Possible damage to instrument

Caution notes indicate any condition or practice, which if not strictly observed or remedied, could result in damage or destruction of the instrument.

**NOTE:** Notes indicate an area or subject of special merit, emphasizing either the product's capability or common errors in operation or maintenance.

- Do not operate this centrifuge in any manner not described in this User Manual. When in doubt or have any troubles with this centrifuge, ASK FOR HELP.
- The precautions described in this User Manual are carefully developed in an attempt to cover all the possible risks. However, it is also important that you are alert for unexpected incidents. Be carefully operating this centrifuge.

**WARNING**

- This centrifuge is not explosion-proof. Never use explosive or flammable samples.
- Do not install the centrifuge in or near places where inflammable gases are generated, or chemicals are stored.
- Do not place dangerous material within 30cm around the centrifuge.
- Make sure to prepare necessary safety measures before using samples that are toxic, radioactive or contaminated with pathogenic micro-organisms at your own responsibility.
- If the instrument, rotor and/or accessories that has been contaminated by solutions with toxic, radioactive or pathogenic materials, clean it according to the decontamination procedure that you are specified.
- If you require services at site, please sterilize and decontaminate it in advance, and then notice the service center involved in the details of the particular materials.
- Do not handle the power cord or turn on or off the POWER switch with wet hands to void electrical shocks.
- For safety purposes, do not enter within 30cm around this centrifuge while it is in operation.
- While the rotor is rotating, never forcedly release the door lock.
- Unauthorized repairs, disassembly, and other services to the centrifuge except by our service center are strictly prohibited.

**CAUTION**

- This centrifuge must be located on one firm and level table.
- Make sure the centrifuge is horizontal before running.
- Make sure the angle between the door and cover is greater than 70 degrees when open the door.
- Be careful not put your fingers or hands between the door and cover when the door off.
- Do not move or relocate this centrifuge while it is running.
- If fluid spills in the rotor chamber, please promptly clean and dry with a dry cloth to avoid sample contamination.
- Ensure to remove any objects and fragments of the tubes dropped inside the rotor chamber before running this centrifuge.
- Cautions on rotors
  - (1) Always check for damages on the rotor surface before using it. Do not use the rotor if an abnormality is found.
  - (2) Do not exceed the allowable imbalance.
  - (3) Use the rotor and tubes within their actual capacities.
- If any abnormal condition occurs during operation, please stop it immediately and contact our service center through your distributor.
- Vibrations are likely to damage the centrifuge, contact our service center if abnormality observed.

**7.1. Rotors**

- The rotors and adapters supplied with this centrifuge are the following:

Angle rotor

PCR strips rotor

0.2mL adapter

0.5mL adapter

**7.2. Tubes**

**7.2.1. Cleaning and sterilizing tubes**

Table - Cleaning and sterilizing conditions for tubes

O: Applicable  
X: Inpplicable

			PA	PC	PP
Cleaning	Cleaning fluids	Acidic (pH 5 or lower)	X	X	X
		Acidic (higher than pH 5)	O	O	O
		Alkaline (higher than pH 9)	O	X	O
		Alkaline (pH 9 or lower)	O	O	O
		Neutral (pH 7)	O	O	O
		Warm water (up to 70°C)	O	O	O
	Ultrasonic cleaning	Neutral detergent (pH 7)	O	O	O
Sterilization	Autoclaving	115°C (0.7kg/cm <sup>2</sup> ) 30 minutes	O	O	O
		121°C (1.0kg/cm <sup>2</sup> ) 20 minutes	X	O	O
		126°C (1.4kg/cm <sup>2</sup> ) 15 minutes	X	X	X
	Boiling	15 to 30 minutes	O	O	O
	Ultraviolet sterilization	200-300nm	X	X	X
	Gas sterilization	Ethylene oxide	O	X	O
Formaldehyde		O	O	O	

PA: Polyallomer

PC: Polycarbonate

PP: Polypropylene

## 6.2. Rotor chamber

### CAUTION

- Do not directly pour water, neutral detergent or disinfectant solution into the rotor chamber. Otherwise fluids may leak into the drive units and cause corrosion or deterioration to the bearings.

If the rotor chamber needs cleaning, clean with cloth or sponge moistened with a neutral detergent solution. Sterilize the chamber by wiping with a cloth moistened with 70% ethanol solution.

## 6.3. Drive shaft

We recommend regular maintenance for drive shaft. You can wipe the drive shaft with soft cloth, and then apply a thin coat of silicon grease.

## 6.4. Door

Clean and sterilize the door using the same method as the step (1) above.

## 6.5. Rotor

If not in use for a long term, then detach the rotor and store it.  
For sample leaks in the rotor, rinse the rotor with water. Apply a thin coat of silicon grease to the rotor when it is completely dry.  
The rotor should be regular maintenance, recommend to cleaning it each 3 months to ensure tube and rotor holes keep clean, and then apply a thin coat of silicon grease.

## 7. Instructions for rotors and tubes

### CAUTION

- Read the instructions thoroughly for the correct use of rotor.
- Do not exceed the allowable maximum speed of tubes; be aware that the allowable maximum speed of some tubes may be lower than the rotor's maximum speed.

## 1. Technical specifications

Model	RS-1008
Code	GDC007
Fixed speed	7000rpm
Maximum RCF	2680xg
Rotors capacity	Angle rotor: 8 x 0.2/0.5/1.5/2mL Rotor for PCR strips: 4 strips PCR de 8 x 0.2mL
Run time	Continuous operation
Motor	DC brush motor
Noise level	≤45dB
Power (adapter AC/DC)	Input: AC 100-240V, 50/60Hz Output: DC 24V
Dimensions (mm)	(L)160×(A)170×(H)122
Weight	0,5kg

## 2. Required operational condition

### 2.1. Basic operational conditions

- Power: 100V-240V, 50Hz/60Hz
- Ambient temperature: 10°C-40°C
- Relative humidity: ≤80%
- No vibration and airflow around.
- No electric dust, explosive and corrosive gases around.

### 2.2. Transport and storage condition

- Storage temperature: -40°C-55°C
- Relative humidity: ≤93%

## 3. Installation

This section describes the instructions that you should abide when install the centrifuge to ensure your safety and the optimum performance. Before moving the centrifuge, the rotor must be removed.

**WARNING**

- Improper power supply may damage centrifuge.
- Make sure the power source conforms to the required power supply before connecting.

1. Place the centrifuge on a firm, flat and level table, ensure the four feet of this centrifuge stand on the table firmly. Avoid installing on the slippery surface or surface prone to vibration.
2. Ideal ambient temperature is  $20^{\circ}\text{C}\pm 5^{\circ}\text{C}$ , avoid placing the centrifuge in direct sunlight if temperature exceeds  $30^{\circ}\text{C}$ .
3. Keep clear of the centrifuge at least 10cm on both sides and at least 30cm behind it to guarantee the cooling efficiency.
4. Keep away from heat or water to avoid sample temperature issues or centrifuge failures.

**4. Rotor preparation****4.1. Prepare the samples****4.2. Inject the samples into tubes****CAUTION**

- Do not overload samples into the centrifuge which will cause leaking

**4.3. Keep the tubes balance**

- Although the centrifuge can accept sample balancing by eye, we recommend that you keep this centrifuge in a well-balanced condition to extend its life expectancy.
- Never intentionally run the centrifuge under unbalanced condition even though the allowable imbalance is not exceeded.

**4.4. Inspect the rotor**


- Check the rotor for damages before using.

**CAUTION**

- Any abnormality such as cracks are found, stop using the rotor.
- Only manufacturer's rotors must be used with the unit.

**4.5. Symmetrically load centrifuge tubes in rotor****CAUTION**

- Make sure the rotor and shaft are tightened, otherwise the rotor may be moved off while rotating and cause damage of the centrifuge and rotor.

**5. Operation**

1. Switch power ON located at the back of the unit.
2. Open the lid and load the centrifuge tubes symmetrically on to the rotor.
3. Close the centrifuge lid to start centrifugation process.
4. Switch power 'OFF' or open the lid and the rotor will stop gradually.
5. Open the lid and collect the samples.

**6. Maintenance****CAUTION**

- If do not follow the recommended instructions for cleaning or disinfecting may damage the centrifuge.

**6.1. Centrifuge**

If the centrifuge is exposed to ultraviolet rays for a long time, the color of the doors may be changed. After using, cover the centrifuge with a piece of cloth to protect it from direct exposure.

If the centrifuge needs cleaning, clean it with a cloth or sponge moistened with a neutral detergent solution.

Sterilize the centrifuge by wiping with a cloth moistened with 70% ethanol solution.