AGITADOR MAGNÉTICO CON CALEFACCIÓN ANALÓGICO Y DIGITAL
ANALOGUE AND DIGITAL MAGNETIC STIRRER WITH HEATING
AGITATEUR MAGNÉTIQUE ANALOGIQUE ET DIGITAL AVEC CHAUFFAGE

Ref. | Code | Réf.
LBB007 - Mod. 670/1
LBD022 - Mod. 670/2

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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TABLE OF CONTENTS

1. INTRODUCTION	4
2. TECHNICAL SPECIFICATIONS	4
3. OPERATION	5

1. INTRODUCTION

The Nahita-Blue magnetic stirrer with heating combines a modern design with a compact and robust structure, offering great versatility of use and easy handling. It is ideal for stirring and heating low-viscosity mixtures, both liquid and solid-liquid. It is available in two versions: analogue and digital.

2. TECHNICAL SPECIFICATIONS

Code	LBB007	LBD022
Model	670/1	670/2
Motor	Brushless motor	
Stand	Not included	Included
Maximum stirring capacity	2 l (H ₂ O)	
Control mode	Knob	Pantalla digital y perilla
Speed	200-1500 rpm	
Stirring bar size	30 mm	
Plate size	Ø 135 mm	
Plate temperature	RT-240 °C	RT-280 °C
Temperature sensor	NO	PT1000
Power	110 W	200 W
Voltage	110-220V 50/60 Hz	
Dimensions	255x155x90 mm	

3. OPERATION

- Application: Mainly used for stirring and heating liquid or solid-liquid mixture with low viscosity.
- Usage: It is widely used as an experimental tool in higher education institutions, environmental protection, health, epidemic prevention, the chemical industry, tap water facilities, and other applications.
- Using steps:

REFERENCE LBB007

- Turn on the device.
- Turn the left knob to adjust the heating power. The red light will turn on, and the device will start heating. Turn the knob again to stop the heating.
- Turn the right knob to adjust the speed. The green light will turn on when stirring begins. Turn the knob again to stop the operation.

REFERENCE LBD022

- Turn on the device.
- Turn the left knob to adjust the temperature and the right knob to adjust the speed. The values will appear on the display.
- Press any knob to start the operation.

■ Instrument use environment:

- Indoor use only
- Altitude: ≤ 2000 m
- Working temperature range: $+5$ °C a $+40$ °C
- Humidity range: $\leq 80\%$
- No vibrations or airflow that could affect performance.
- No conductive dust, explosive gases, or corrosive gases in the surrounding air.
- Maintain a safety distance of 10 cm around the instrument while it is in operation.

■ Safety precautions:

- Do not unplug or plug in the power connector or toggle the power button with wet hands.
- Do not unplug the power plug while the instrument is powered on.
- Do not perform maintenance or cleaning while the instrument is powered on.
- Do not install the instrument on an uneven, swaying, or vibrating work surface.