## 9. STOCKAGE ET TRANSPORT

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#### INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

- ◆ The equipment cannot reach the set temperature:
  - Check the room temperature is not higher than the set temperature.
 Stir speed cannot reach set point:
  - Excessive medium viscosity may cause abnormal speed reduction of the motor
- ◆ Heating cannot be started after set the temperature , or stirring cannot be started when adjust the control knob.
  - Check the control panel has damages which may have arisen during transport.
 If these faults are not solved, please contact your distributor to send the equipment to Nahita Technical Service Department.

## 8. MAINTENANCE AND CLEANING

- ◆ Proper maintenance helps to keep instruments working in a good state and lengthen its lifetime.
- ◆ Be careful not spray the detergent or water into the instrument when cleaning.
- ◆ Unplug the power line when cleaning.
- ◆ Only use the following products to clean the equipment:
  - For dyes: Isopropyl alcohol
  - For construction materials: Water containing tenside / Isopropyl alcohol
  - For cosmetics: Water containing tenside / Isopropyl alcohol
  - For food: Water containing tenside
  - For fuels: Water containing tenside
- ◆ Wear the proper protective gloves during cleaning of the instrument.
- ◆ Before using other method for cleaning or decontamination, please contact your distributor to avoid damages on the equipment.
- ◆ The enamel makes the hotplate easier to care for and more resistant to acids and bases. Because of it, however, the heating plate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.
- ◆ The instrument must be cleaned and put it into the initial packaging carton before sending to service for repair, avoiding the contamination of hazardous.
- ◆ Use the instrument in a dry clean room and stable temperature environment.

## 9. STORAGE AND TRANSPORTATION

- ◆ Keep it in dry and clean room with good ventilation and no corrosive gas
- ◆ Prevent it from wetting by the rain and avoid violent collision in transportation.

## INDEX OF LANGUAGES

Spanish .....	2 - 10
English .....	11 - 19
French .....	20 - 28

## INDEX OF CONTENTS

Preface and Warranty.....	11
1. Safety Instructions .....	12
2. Proper Use.....	14
3. Inspection.....	14
Reception	
Listing of Items	
4. Trial run.....	14
5. Control panel.....	15
6. Operating Modes.....	15
Analogical models	
Digital model	
7. Problems and solutions.....	17
8. Maintenance and Cleaning.....	17
9. Storage and transportation.....	18
10. Main technical features.....	18
11. Working condition.....	18

## PREFACE AND WARRANTY

Thank you for purchasing a Nahita Blue magnetic stirrer. Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

If it is necessary, contact your dealer to get the necessary help and advises for a proper use of the equipments. For this, user must give information about the serial number of the equipment, the problems found, the procedures followed to try to solve the problem and contact information.

Nahita Blue stirrers are warranted to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty please contact your distributor in order to send the equipment to Nahita Technical Service Department.

## 1. SAFETY INSTRUCTIONS

Read the user's manual before operating with the equipment. This equipment must be used only by qualified staff.

Be careful when touching the equipment: heating plate and metal case. Heating plate can reach temperatures up to 380°C.

Once the equipment has been turned off it can remain hot for several minutes

Connect the device to an earthed power supply to ensure safety of machine and experiment; the voltage indicated on the equipment must correspond to that of the power supply used.

This equipment is forbidden to use with inflammable, explosive, poisonous and strong corrosive solutions.

Pay attention to the safety and resistance to high temperatures of the container used. If the container is damaged, liquid could fall into the equipment with the risk of fire.

Make sure the equipment is placed and installed in a plain, horizontal and stable surface.

Pay attention to the set temperature; never leave unattended the equipment during operation and with heating function on.

Pay attention to overloading of material in relation to the capacity of the equipment. An overloading could produce overheating and damages in the equipment.

Non-authorized personnel is not allowed to disassemble and repair this machine. If the equipment presents a problem, contact to your distributor

faults. The heating system will be automatically turned off, the alarm will sound and the indicator "ALM" will be displayed. Please check if there is any visible damages on temperature sensor and if it is well connected to the equipment.

If "SEr" is displayed, it means the temperature sensor falls off and "ALM" indicator will be displayed, Please check if there is any visible damages on temperature sensor.

◆ Press any key to make it silence when the alarm sounds.

If no key is pressed for 1 min, display will automatically enter normal display mode.

### 4 Internal temperature parameter setting

◆ The model 692/1 has two temperature sensor; one in the heating plate and the external sensor Pt 100.

◆ Default parameters operate with the inner temperature sensor, to change it follow the next steps:

◆ Press [Ω] button until screen display Lc 0000, choose number 9 and press [Ω] again to confirm the selection.

◆ LCD displays EnS (default is 0) and choose the operating mode you want:

0: Readings with heating plate sensor

1: Readings with both sensors (Press [▼] to display the two temperatures)

2: Readings with external sensor

To exit the parameter setting press [Ω] key

### 5 Auto-tuning function

◆ The auto-tuning function, allows to establish automatically those inner parameters that are most suitable for a certain working temperature; this can be very useful in those applications in which the same solution and the same temperature are going to be always used.


◆ For this, place the container with the liquid to be stirred, put the temperature sensor inside and enter the auto-tuning mode by pressing  button for 5 seconds. The indicator "AT" will flash on display, and will stop flashing when auto-tuning finishes. The auto-tuning process can be cancelled by pressing  for 5 seconds again. When the auto-tuning process is finished, the most adequate inner parameters for the selected temperature will be automatically established.

## 7. PROBLEMS AND SOLUTIONS

The equipment does not turn on:

- Check whether the power cable is plugged
- Check whether the fuse is broken or loose

## 2 Control panel

### Key definition:

- Set the temperature value and inner parameters.
- In setting mode, select the digit to be modified. Out of setting mode, if the key is pressed for 5 s, it enters the auto-tuning program.
- Decrease the value.
- In normal display mode, press this key to turn off or on the backlit display. In setting mode, press this key to increase the value.

### 3 Operation mode

When turning the equipment on, all signs on display are lighted. After 3 s, display enters into normal display mode in which the value of real temperature is shown at the upper part and stirring speed at the lower part of display.

#### Speed setting

- ◆ Turn the speed adjustment knob clockwise until reaching the desired speed; the stirring speed value in rpm will be displayed.

#### Temperature setting

- ◆ Press to enter temperature setting mode. The indicator "SP" will be displayed at the upper part of display and below, the temperature set value; Use keys to change the set value and press to exit setting mode and save temperature set value.
- ◆ When real temperature is lower than set temperature, the heating system will start working and "HEAT" indicator will be displayed.

#### Over temperature alarm

- ◆ When real temperature is higher than 10°C the set temperature, then the alarm sounds and indicator "ALM" is displayed. If this happens because set temperature has been changed into a lower value, the "ALM" indicator will be displayed but no alarm will sound.
- ◆ In any case, the heating system of the equipment will be automatically turned off.

#### Abnormal temperature measurement alarm

- ◆ At the upper part of display indicator "Er-2" will be shown; this means that the temperature sensor faults, temperature exceeds the measuring range of sensor or controller itself

- ◆ During operation wear the necessary personal protections to avoid the risk of possible damages as:
  - Burns caused by splashing and evaporation of liquids
  - Intoxication caused by release of toxic or combustible gases.
- ◆ Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface; do not operate with the equipment in explosive atmospheres, with hazardous substances or under water.
  - Gradually increase or reduce the speed if:
    - The stirring bar breakaway because of too high speed
    - The instrument is not running smoothly, or container moves on the base plate.
- ◆ Temperature must always be set to at least 50A lower than the ignition point of the media used.
- ◆ Beware of hazards due to:
  - Flammable material or media with a low boiling temperature
  - Overloading of media
  - Unsafe container
- ◆ Process pathogenic materials only in closed vessels.
- ◆ In case the stirrer bar used is made of PTFE, please note:
  - ◆ Elemental fluorine, three fluoride and alkali metals will corrode the PTFE, and halogen alkenes make it expansion at room temperature. Molten alkali, alkaline earth metals or their solution, as well as the powder in second and third ethnic of the Periodic Table of elements will have chemical reaction with PTFE when temperature reaches 300 ~400 A
  - ◆ Before operation, always check the instrument and accessories are not damaged.
  - ◆ Do not use damaged components. Safe operation is only guaranteed with the accessories described in the "Listing of items" chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.
  - ◆ Ensure that the external temperature sensor is inserted in the media to a depth of at least 20 mm (only model 692/1).
  - ◆ When using metal vessels, do not place the temperature sensors on the bottom of the vessel.
  - ◆ Placing sensors on the vessel bottom can cause excessively high temperature readings, especially in media which have poor conductivity. The tip of the measuring sensor must be at least 5 mm from the vessel bottom, a distance of 10 mm is ideal.
  - ◆ The instrument can only be disconnected from the main power supply by pulling out the connector plug, not the cable.
  - ◆ The voltage stated on the label must correspond to the main power supply.
  - ◆ Ensure that the main power supply cable does not touch the heating base plate. Do not cover the device.
  - ◆ Keep away from high magnetic field.
  - ◆ Respect the minimum safety distances between devices, between the device and the wall and above the assembly (min.100 mm).

## 2. PROPER USE

The instrument is designed for mixing and / or heating liquids in schools, laboratories or factories. This device is not suitable for using in residential areas or environments that could be dangerous for user or instrument, as stated in Chapter 1.

## 3. INSPECTION

### 3.1 Reception

Unpack the equipment carefully and check for any damages which may have arisen during transport. If it happens, please contact your distributor.

**NOTE: If there is any apparent damage on the equipment, do not plug it into the power line.**

### 3.2 Listing of Items

Nahita Blue stirrers are supplied with the following accessories:

- Main unit	1 pc
- Power Cable	1 pc
- Magnetic bar	1 pc
- Rod support	1 pc
- Temp. sensor (model 692/1)	1 pc
- User's manual	1 pc

## 4. TRIAL RUN


- ◆ Make sure the required operating voltage and power supply voltage match.
- ◆ Ensure the socket is earthed reliably.
- ◆ Ensure the power switch of the equipment is off
- ◆ Plug in the power cable and turn the equipment on.
- ◆ Pour the liquid to be stirred in a vessel with a stirring bar.
- ◆ Put the vessel on the stirring plate.
- ◆ Set the stirring speed, the equipment starts stirring.
- ◆ Observe the stirring bar and LCD display (model 691/10).
- ◆ In those model with heating function (models 692 and 692/1), set the temperature; the equipment starts heating.
- ◆ Observe the real temperature on LCD display (model 692/1).
- ◆ Stop the heating and stirring functions, and turn the equipment off.

If the operations above are normal, the device is ready to operate. If these operations are not normal, the device may be damaged during transportation, please contact your distributor.

## 5. CONTROL PANEL


### Analogical models Models 682 y 692

- 1.1 Temperature adjustment knob (only 692 model)
- 1.2 Speed adjustment knob
- 1.3 Indicator pilot for speed and heat operation

### Digital model Model 692/1

- 2.1 Temperature adjustment knob
- 2.2 Digit selection
- 2.3 Increase button
- 2.4 Decrease button
- 2.5 Speed adjustment knob
- 2.6 LCD display

## 6. OPERATING MODES

### 6.1 Analogical models:

- ◆ Place the equipment on a plain, stable surface and then put the container with liquid on the hotplate.
- ◆ Switch the equipment on, the power light will light on.
- ◆ Turn the speed knob clockwise until reaching the desired stirring speed. The corresponding pilot will light on.

If the equipment is provided with heating function (model 692), turn the temperature knob clockwise until reaching the desired temperature. The corresponding pilot will light on.

### 6.2 Digital model:

#### 1 Main technical features

- Temperature sensor: Pt100 Thermal Resistive
- Hall type speed sensor
- Temperature setting range: 0-400°C; Temperature Measuring Range: -10/410°C
- Error of Temperature Measurement: < 0.5%