AUXILAB S.L.

Material de laboratorio
Laboratory supplies

DECLARACIÓN DE CONFORMIDAD CE
CE DECLARATION OF CONFORMITY
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi polarimeter model 412. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the Zuzi polarimeter model 412, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the polarimeter avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the polarimeter to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the polarimeter.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler.
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

7. TROUBLESHOOTING

<i>TROUBLE</i>	<i>CAUSE</i>	<i>SOLUTION</i>
Problem with light source	Light source is dirty Light source is damaged	Open the lid and clean the lamp Replace the lamp
After pressing the on/off switch, the equipment does not turn on	Energy does not reach the equipment The fuse is blown	Check the cable and the socket Check the 3 A fuse and replace if necessary
Too high noise	Mechanical friction	Gears are not well lubricated

8. RECOMMENDED ACCESSORIES

<i>Code</i>	<i>Description</i>
90404910	Test tube of 100 mm
90404920	Test tube of 200 mm
90404901	Glass joints for test tubes, 15 mm Ø
90404691	Circular rubber washer for test tube
90404692	Test tube stopper
90400001	Quartz piece for polarimeter calibration, value +5°
90400002	Quartz piece for polarimeter calibration, value -5°
90400003	Quartz piece for polarimeter calibration, value +17°
90400004	Quartz piece for polarimeter calibration, value -17°
90400005	Quartz piece for polarimeter calibration, value +35°
90400006	Quartz piece for polarimeter calibration, value -35°

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this polarimeter by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

6. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The polarimeter is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the polarimeter by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the polarimeter (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Zuzi Technical Assistance Department.

IMPORTANT: BEFORE CHANGING THE LAMP OR THE FUSE, MAKE SURE THAT THE POLARIMETER IS DISCONNECTED FROM THE NET. CHECK YOU ARE USING THE PROPER LAMPS, AS OTHER TYPES CAN PROVOKE MALFUNCTION.

Replacement of the fuse

- ◆ Fuses (2.1) are located at the back side of the polarimeter.
- ◆ With the help of a screw driver take the fuse compartment out and replace the fuse.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the polarimeter and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

¡ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

¡ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish	2-11
English	12-21
French	22-32

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	13
2. DESCRIPTION	14
3. TECHNICAL SPECIFICATIONS	15
4. THEORETICAL FOUNDATIONS	15
5. INSTALLATION / SETTING UP	17
6. MAINTENANCE AND CLEANING	20
7. TROUBLESHOOTING	21
8. RECOMMENDED ACCESSORIES	21
9. ANNEX I: CE CERTIFICATE	32

1. USES OF THE INSTRUMENT

The Zuzi polarimeter model 412 is a complete equipment that allows determination of the optical rotation of optically active substances. Thus, it is very useful in several fields of industrial and quality control processes as well as in research. Among the main application of the polarimeter there are:

- ◆ Research: isolation of crystallized substances, evaluation and characterization of optically active compounds, kinetics and changes in concentration.
- ◆ Pharmaceutical industry: determination of the optical rotation of medicines (antibiotics, codeine, analgesics) and drugs.
- ◆ Food industry: determination of the sugar and condiment content.
- ◆ Chemical industry: study of biopolymers, synthetic and natural polymers, etc.

2. DESCRIPTION

Figure 1-Front side

- 1.1 Control panel
- 1.2 Sample compartment


Figure 2-Back side

- 2.1 Fuses 1 and 3.15 A
- 2.2 RS232 output
- 2.3 Socket
- 2.4 Grille
- 2.5 Power on/off



Figure 3-Control panel

- 3.1 Display
- 3.2 Clear: zero
- 3.3 Measure
- 3.4 Re-M: repetition of measurement
- 3.5 Average: average of repeated measurements
- 3.6 Central button (1, 2, 3): results of repeated measurements Re-M

- ◆ The software has five different functions:

- Read: to perform sample measurements and register the last three measurements and their average value.
- Open: to open a previously saved file.
- Print: to print measured data (the computer must be connected to a printer).
- Save: to save the measured data.
- Exit: to close the program.

Concentration and content measurements

- ◆ Prepare several standard samples of different concentration and measure their optical rotation value. Represent the curve with the concentrations and their corresponding optical rotation values.
- ◆ Measure the optical rotation of the sample to be analysed and use the curve to get the concentration value of that sample.

Security

- ◆ The polarimeter must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ Put the polarimeter on top of a horizontal, plane, stable table, having a safety area of at least 30 cm per side.
- ◆ Do not place the polarimeter near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the working area.
- ◆ When you are not using the polarimeter for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement a fuse).
- ◆ Never try to repair the polarimeter by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the polarimeter (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

unplugged and protected from dust in order to avoid possible accidents and prolong its working life.

Setting up

- ◆ Turn the equipment on by pressing power on/off (1.3); "0.000" will be displayed.
- ◆ The equipment is ready to operate under the normal procedure.

Measurements

- ◆ To obtain the optical rotation (α) value of samples proceed as follows:
- ◆ Put the tube with distilled water or the blank solution into the sample compartment (1.2), let the displayed reading stabilize and press Clear (3.2).

Note: When placing the tubes make sure they are in the correct position and direction. Make sure the sample does not have air bubbles since measurements could be affected.

- ◆ Remove the blank solution or distilled water from the sample compartment and put the sample to be measured; the equipment will automatically measure and display the result. When result stabilizes #1 appears on display.
- ◆ To repeat measurement press Re-M (3.4). Up to 2 repetitions can be performed and registered in the equipment (#2 and #3 will be displayed). To display the result of the different measurements press the Central button (1,2,3) (3.6). By pressing the Average button (3.5) the average value of the different measurements will be displayed (AV will be displayed).
- ◆ Before measuring a new sample, it is necessary to make the zero. For this:
 - Remove the sample with the already analyzed sample from the sample compartment
 - Put the blank sample or distilled water, let the displayed reading stabilize and press Clear (3.2)
 - Put the new sample and the equipment will automatically measure and display the result

Note: before turning the equipment off, remove the last analyzed sample from the sample compartment, let the displayed reading stabilize and press Clear (3.2); then turn the equipment off.

Software installation

- ◆ Insert the CD-Rom and open the installation program by making double click on the Setup file. Accept all steps until the installation process finishes.
- ◆ Connect the polarimeter to the computer and turn it on.
- ◆ To start working, look for the software WZZ-2 and execute the file WZZ-2.
- ◆ Use the compatible RS port (COM-1 or COM-2) for communication with the equipment.


3. TECHNICAL SPECIFICATIONS

Code	50410012
Model	412
Range	$\pm 45^\circ$
Accuracy	$\pm (0.001^\circ + \text{valor medido} \times 0.05\%)$
Min. displayed value	0.005°
Speed	1.30°/Segundo
Monochromatic light source	LED light and filter (λ : 589.44 nm)
Output	RS232
Tubes	100 y 200 mm (included)
Power	220V 50Hz $\pm 10\%$
External dimensions	60x32x22 cm
Weight (approx.)	18 Kg

Note: RS232 parameters: Baudios 9600, Stop 1 bit, Data 8 bits.

4. THEORETICAL FOUNDATIONS

Figure 4. Internal diagram of the polarimeter

- 4.1 Light source
- 4.2 Aperture
- 4.3 Lens
- 4.4 Coloured filter
- 4.5 Polarizer A
- 4.6 Magnetic rotation coil
- 4.7 Sample
- 4.8 Polarizer B

- 4.9 Photo-multiplier
- 4.10 Pre-amplifier
- 4.11 Automatic high voltage
- 4.12 Frequency amplifier
- 4.13 Amplifier power
- 4.14 Servo motor
- 4.15 Worm-gear
- 4.16 Counter

Zuzi polarimeter 412 is based on the optical null principle to obtain accurate readings.

For this purpose, it uses a LED diode as the light source. Firstly, the aperture (4.2) and the lens (4.3) form a parallel light beam, that, after passing through the polarizer A (4.4) becomes a polarized light beam whose vibrating plane is $00'$ (Figure 5a). When the polarized light beam passes through the Faraday coil (4.5), its vibrating plane will swing in a $\pm\beta$ angle at a frequency of 50 Hz. After passing through the polarizer B (4.8) the light beam is projected into the photo-multiplier (4.9) so as an alternating current electric signal is generated.

According to the optical null principle, the equipment takes as the optical null position ($\alpha = 0^\circ$) the point in which the axes of the polarizers A and B are orthogonally arranged ($00' \perp PP'$) and therefore the light beam is completely blocked and does not reach the photo-multiplier.

When an optically active sample is placed, the vibrating plane will rotate into an α angle (Figure 5b) (optical rotation of the sample) at a certain frequency F that finally will become an electric signal of the same frequency F. This signal makes the servo motor spin (4.14), which at the same time makes the polarizer A turn in an α angle until the optical null position is restored and the frequency F signal disappears (Figure 5c).

The α angle that the polarizer A has rotated will be displayed as the reading of the optical rotation of the analyzed sample.

Figure 5a


Figure 5b



Figure 5c

Note: $00'$ is the polarization axe of polarizer A and PP' is the polarization axe of polarizer B

Figure 5a: The polarized light generated by the polarizer A vibrates in the plane $00'$.

Figure 5b: When light passes through an optically active sample its vibrating plane swings in a α angle.

Figure 5c: The polarizer A rotates in a α angle to compensate the optical rotation of the sample and to restore the optical null position; the α value will be displayed with the correct sign.

5. INSTALLATION / SETTING UP

Preliminary inspection

◆ Unwrap the polarimeter, take off the involving plastic and take off the protection in which it comes fitted. Take off all the protective items and, without connecting the polarimeter to the net, make sure that it does not present any damage because of the shipment. In case the polarimeter presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

◆ Please check that all the accessories are enclosed with the equipment:

- | | |
|-------------------------|-----------------|
| - 1 tube 100 mm-long | - Fuses |
| - 1 tube 200 mm-long | - Screw driver |
| - Set of rubber washers | - Schuko wire |
| - Set of glass joints | - User's manual |

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

◆ Please put the polarimeter on top of a horizontal, plane and stable table making a free space of at least 30 cm per side.

◆ Do not put the polarimeter near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc. Avoid vibrations in the working area.

◆ Avoid the presence of inflammable or toxic substances in the working area.

◆ The polarimeter is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz $\pm 10\%$ provided with earth wire and to the other end to the polarimeter connector (2.3).

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz $\pm 10\%$.

When you are not using the polarimeter for a long period of time please make sure it is