

Test de compatibilidad de pipetas

	10 µl		20 µl	100 µl	200 µl			1.000 µl													
	P. sin filtro		P. con filtro	P. con filtro	P. con filtro	P. sin filtro		P. con filtro	P. sin filtro												
	MJB032	MJB034	MJB035	MJB033	MJB036	MJD010	MJD011	MJD012	MJG005	MJG006	MJJ040	MJJ042	MJJ043	MJU041	MJU044	MJN038	MJN040	MJN041	MJN039	MJN042	
<b>BioHit</b>																					
BM-3 Biohit mLine pipeta manual, 0.1-3ul	●	●	●																		
BM-10 Biohit mLine pipeta manual, 1-10ul	●	●	●																		
BM-20 Biohit mLine pipeta manual, 2-20ul						●	●	●			●	●	●								
BM-100 Biohit mLine pipeta manual 10-100ul									●	●	●	●	●								
BM-200 Biohit mLine pipeta manual 20-200ul											●	●	●	●	●						
BM-1000 Biohit mLine pipeta manual 100-1000ul																●	●	●	●	●	●
BU-20 Pipeta Proline Plus 2-20ul						●	●	●			●	●	●								
BU-100 Pipeta Proline Plus 10-100ul									●	●	●	●	●								
BU-200 Pipeta Proline Plus 20-200ul											●	●	●	●	●						
BU-1000 Pipeta Proline Plus 100-1000ul																●	●	●	●	●	●
BH-20 Biohit Proline pipeta manual, 2-20ul											●	●	●								
BH-50 Biohit Proline pipeta manual, 5-50ul											●	●	●								
BH-100 Biohit Proline pipeta manual, 10-100ul									●	●	●	●	●								
BH-200 Biohit Proline pipeta manual, 20-200ul											●	●	●	●	●						
BH-1000 Biohit Proline pipeta manual, 100-1000ul																●	●	●	●	●	●
BH8-50 Biohit Proline pipeta manual 8 canales, 5-50ul in 0.5ul											●	●	●								
BL-120A 5-120ul con adaptador											●	●	●								
BP-100 5-100ul en incrementos de 1ul									●	●	●	●	●								
BP-500 10-500ul en incrementos 5ul																●	●	●	●	●	●
BP-1000 50-1000ul en incrementos 10ul																●	●	●	●	●	●
BP12-100 doce canales, 5-100ul en incrementos 1ul											●	●	●	●	●						
<b>Capp</b>																					
CU-2 Pipeta CappAero, 0.2-2 ul Gris estándar	●	●	●								●	●	●								
CU-10 Pipeta CappAero, 1-10 ul Gris estándar	●	●	●								●	●	●								
CU-20 Pipeta CappAero, 2-20 ul Gris estándar						●	●	●			●	●	●								
CU-50 Pipeta CappAero, 5-50 ul, Ecoline, Rojo											●	●	●								
CU-100 Pipeta CappAero, 10-100 ul, Ecoline, Naranja									●	●	●	●	●								
CU-200 Pipeta CappAero, 20-200 ul, Ecoline, Amarillo											●	●	●	●	●						
CU-1000 Pipeta CappAero, 100-1000 ul, Ecoline, Azul																●	●	●	●	●	●
<b>Eppendorf</b>																					
E2000-10 desde 1-10ul Air-Displacement Micropipeta																					
E2000-20 desde 2-20ul Air-Displacement Micropipeta																					
E2000-100 desde 10-100ul Air-Displacement Micropipeta									●	●	●	●	●								
E2000-200 desde50-200ul Air-Displacement Micropipeta											●	●	●	●	●						
E2000-1000 desde 100-1000ul Air-Displacement Micropipeta																●	●	●	●	●	●
E2100-2.5 1 desde 0.1-2.5ul Air-Displacement Micropipeta																					
E2100-10 desde 1-10ul Air-Displacement Micropipeta																					
E2100-20 desde 2-20ul Air-Displacement Micropipeta						●	●	●			●	●	●								
E2100-100 desde10-100ul Air-Displacement Micropipeta									●	●	●	●	●								
E2100-200 desde 20-200ul Air-Displacement Micropipeta											●	●	●	●	●						
E2100-1000 desde 100-1000ul Air-Displacement Micropipeta																●	●	●	●	●	●
ES-1000 Research Plus																●	●	●	●	●	●
ES-100 Research Plus 10-100ul									●	●	●	●	●								
ES-200 Research Plus 20-200ul											●	●	●								
EE-100 desde 10-100ul Pipeta electrónica multicanal									●	●	●	●	●								
EE-1000 desde 100-1000ul Pipeta electrónica																●	●	●	●	●	●
EE-10																					
EE12-10																					
<b>Fisher Scientific</b>																					
FFI-50 FisherBrand Finnpiptete I 5-50 ul											●	●	●								
FFI-200 FisherBrand Finnpiptete I 20-200 ul											●	●	●	●	●						
FFI-1000 FisherBrand Finnpiptete I 100-1000 ul																●	●	●	●	●	●
<b>Gilson</b>																					
P-10 Micropipeta Pipetman Air-Displacement, 1-10 ul	●	●	●	●	●																
P-20 Micropipeta Pipetman Air-Displacement, 2-20 ul						●	●	●			●	●	●								
P-20 Micropipeta Pipetman multicanal Air-Displacement																					
P-100 Micropipeta Pipetman Air-Displacement, 10-100 ul									●	●	●	●	●								
P-200 Micropipeta Pipetman Air-Displacement, 50-200 ul											●	●	●	●	●						
P-200 Micropipeta Pipetman multicanal Air-Displacement																					
<b>Labnet/HTL</b>																					
P-1000 Micropipeta Pipetman Air-Displacement, 100-1000 ul																●	●	●	●	●	●
LN-2 Labnet Labpette Discovery 0.2 hasta 2ul autoclavable	●	●	●	●	●																
LN-20 Labnet Labpette Discovery 2 hasta 20ul autoclavable						●	●	●			●	●	●								
LN-100 Labnet Labpette Discovery 10 hasta 100ul autoclavable									●	●	●	●	●								
LN-200 Labnet Labpette Discovery 20 hasta 200ul autoclavable											●	●	●	●	●						
LN-1000 Labnet Labpette Discovery 100 hasta 1000ul autoclavable																●	●	●	●	●	●
<b>Nichiryo</b>																					
NPX-2 Nichipet EX Autoclavable/UV digital estable 0.2-2ul	●	●	●	●	●																
NPX-10 Nichipet EX Autoclavable/UV digital estable 1-10ul	●	●	●	●	●																
NPX-20 Nichipet EX Autoclavable/UV digital estable 2-20ul						●	●	●			●	●	●								
NPX-100 Nichipet EX Autoclavable/UV digital estable 10-100ul									●	●	●	●	●								
NPX-200 Nichipet EX Autoclavable/UV digital estable 20-200ul											●	●	●	●	●						
NPX-1000 Nichipet EX Autoclavable/UV digital estable 100-1000ul																●	●	●	●	●	●
<b>Rainin</b>																					
ED1-200R Rainin EDP Pipeta electrónica, 20 - 200 ul											●	●	●	●	●						
ED2-25R Rainin EDP 2 Pipeta electrónica, 2.5 - 25 ul											●	●	●								
ED2-100R Rainin EDP 2 Pipeta electrónica, 10 - 100 ul									●	●	●	●	●								
ED2-1000R Rainin EDP 2 Pipeta electrónica, 100 - 1000 ul																●	●	●	●	●	●
ED-100R Rainin EDP Plus Pipeta electrónica, 10 - 100 ul									●	●	●	●	●								
<b>Socorex</b>																					
ED-1000R Rainin EDP Plus Pipeta electrónica, 100 - 1000 ul																●	●	●	●	●	●
S825-20 Acura® 825 Volumen variable 2-20						●	●	●			●	●	●								
S825-50 Acura® 825 Volumen variable 5-50											●	●	●								
S825-100 Acura® 825 Volumen variable10-100									●	●	●	●	●								
S825-200 Acura® 825 Volumen variable 20-200											●	●	●	●	●						
S825-1000 Acura® 825 Volumen variable 100-1000																●	●	●	●	●	●