

Test de compatibilidad de pipetas

| | Modelo | Capacidad | 10 µl | | | | | 20 µl | | | 100 µl | | | 200 µl | | 1.000 µl | | | | | |
|--------------------------------------|--------|-------------|---------------|--------|---------------|--------|--------|---------------|--------|---------------|---------------|--------|---------------|---------------|--------|---------------|--------|--------|--------|--------|--------|
| | | | P. sin filtro | | P. con filtro | | | P. con filtro | | P. con filtro | P. sin filtro | | P. con filtro | P. sin filtro | | P. con filtro | | | | | |
| | | | MJB032 | MJB034 | MJB035 | MJB033 | MJB036 | MJD010 | MJD011 | MJD012 | MJC005 | MJC006 | MJU040 | MJU042 | MJU043 | MJU041 | MJU044 | MJN038 | MJN040 | MJN041 | MJN039 |
| Digipette <i>ECO</i> | MGB019 | 0,2-2 ul | | | | | | | | | | | | | | | | | | | |
| | MGB020 | 1,0-10 ul | | | | | | | | | | | | | | | | | | | |
| | MGB021 | 2,0-20 ul | | | | | | | | | | | | | | | | | | | |
| | MGB022 | 20-100 ul | | | | | | | | | | | | | | | | | | | |
| | MGB023 | 50-200 ul | | | | | | | | | | | | | | | | | | | |
| | MGB024 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| <i>Autoclavable vol variable</i> | MGB001 | 0,1-2,5 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB002 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB003 | 2,0-20 ul | | | | | | | | | | | | | | | | | | | |
| | MGB004 | 5,0-50 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB005 | 10-100 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB006 | 20-200 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGB007 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| <i>Autoclavable volumen fijo</i> | MGD001 | 10 ul | • | • | • | • | • | | | | | | | | | | • | • | • | • | • |
| | MGD002 | 20 ul | | | | | | | | | | | | | | | | | | | |
| | MGD003 | 25 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD004 | 50 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGD005 | 100 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGD006 | 200 ul | | | | | | | | | | • | • | • | • | • | | | | | |
| | MGD007 | 500 ul | | | | | | | | | | | | | | | | • | • | • | • |
| | MGD008 | 1000 ul | | | | | | | | | | | | | | | | • | • | • | • |
| <i>Autoclavable 8 canales</i> | MGG001 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGG002 | 5,0-50 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGG003 | 10-100 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGG004 | 20-200 ul | | | | | | | | | | • | • | • | • | • | | | | | |
| | MGG005 | 30-300 ul | | | | | | | | | | | • | • | • | • | • | | | | |
| <i>Autoclavable 12 canales</i> | MGG006 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGG007 | 5,0-50 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGG008 | 10-100 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGG009 | 20-200 ul | | | | | | | | | | • | • | • | • | • | | | | | |
| Dosypette | MGG010 | 30-300 ul | | | | | | | | | | | | | | | | | | | |
| | MGM039 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGM040 | 5,0-50 ul | | | | | | | • | • | • | • | • | • | • | | | | | | |
| | MGM041 | 30-300 ul | | | | | | | | | | | | | | | | | | | |
| Accumax <i>Smart</i> | MGM042 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| | MGB010 | 0,1-2,5 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB011 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB012 | 2,0-20 ul | | | | | | | | | | | | | | | | | | | |
| Socorex | MGB013 | 5,0-50 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB014 | 10-100 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB015 | 20-200 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGB016 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| | MGB028 | 1,0-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB030 | 2,0-20 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| <i>Calibra digital 822</i> | MGB031 | 10-100 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB032 | 20-200 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGB026 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| | MGB034 | 0,1-2 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB035 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB036 | 1,0-10 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB037 | 2,0-20 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB038 | 5,0-50 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| <i>Acura manual 825</i> | MGB039 | 10-100 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGB040 | 20-200 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGB033 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| | MGB050 | 0,1-2 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB051 | 0,5-10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGB052 | 1,0-10 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGB053 | 2,0-20 ul | | | | | | | • | • | • | • | • | | | | | | | | |
| <i>Acura manual 826 XS</i> | MGB054 | 5,0-50 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGB055 | 10-100 ul | | | | | | | | • | • | • | • | • | | | | | | | |
| | MGB056 | 20-200 ul | | | | | | | | | | | | | | | | | | | |
| | MGB057 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| | MGB054 | 5,0-50 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGB055 | 10-100 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGB056 | 20-200 ul | | | | | | | | | | | | | | | | | | | |
| <i>Acura manual 815</i> | MGB057 | 100-1000 ul | | | | | | | | | | | | | | | | | | | |
| | MGD014 | 1 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGD015 | 5 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGD016 | 10 ul | • | • | • | • | • | | | | | | | | | | | | | | |
| | MGD017 | 10Y ul | | | | | | | • | • | • | • | • | | | | | | | | |
| | MGD018 | 15 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD019 | 20 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD020 | 25 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD021 | 30 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD022 | 32 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD023 | 40 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD024 | 50 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| | MGD025 | 60 ul | | | | | | | | | • | • | • | • | • | | | | | | |
| MGD026 | 70 ul | | | | | | | | | • | • | • | • | • | | | | | | | |
| MGD027 | 75 ul | | | | | | | | | • | • | • | • | • | | | | | | | |
| MGD028 | 80 ul | | | | | | | | | • | • | • | • | • | | | | | | | |

