

R

PROTECTION AND SECURITY

RB - Protection clothes	410
RBB - Coats.....	410
RBC - Protection clothes.....	410
RBD - Head cap.....	411
RBG - Glasses.....	411
RBJ - Masks	415
RBK - Aprons.....	415
RBL - Gloves.....	415
RBN - Shoe cover.....	416
RD - Showers and eyewashes	417
RDB - Stations.....	417
RDD - Eyewash bottles	417
RDZ - Others	419
RZ - Others	419
RZB - Table protectors	419
RZD - Wound treatment	419





RBB - Coats

Lab-coats, 100% cotton



- 1 | Snap buttons, 3 pockets and lapel collar

Code	For	Size
RBB001	Man	S
RBB002	Man	M
RBB003	Man	L
RBB004	Man	XL
RBB009	Woman	S
RBB010	Woman	M
RBB011	Woman	L
RBB012	Woman	XL

Lab-coats, cotton/polyester



- 1 | 65% polyester - 35% cotton
- 2 | Optical white colour, snap buttons, 3 pockets and lapel collar.

Code	For	Size
RBB005	Man	S
RBB006	Man	M
RBB007	Man	L
RBB008	Man	XL
RBB013	Woman	S
RBB014	Woman	M
RBB015	Woman	L
RBB016	Woman	XL

RBB - Coats

Lab-coats disposables



- 1 | Long-sleeved, single-use, polypropylene
- 2 | Velcro fastening (5 velcro fasteners), shirt collar, elasticated cuffs.
- 3 | 100% Latex-free product.

Code	Dimensions	Weight	Colour
RBB019	112x138 cm	35gr, PP	White
RBB020	112x138 cm	35gr, PP	Green
RBB021	112x138 cm	35gr, PP	Blue
RBB022	112x138 cm	25gr, PP	White
RBB023	112x138 cm	25gr, PP	Green
RBB024	112x138 cm	25gr, PP	Blue

RBC - Protection clothes

Overalls with hood, PP



- 1 | Single use long sleeve coverall made of Polypropylene,
- 2 | Zip fastening.
- 3 | Elasticated hood and cuffs.
- 4 | Hood for head protection.
- 5 | 100% Latex free product.

Code	Size	Colour
RBC001	M	White
RBC002	L	White
RBC003	XL	White
RBC004	XXL	White
RBC005	XXXL	White



RBD - Head cap

▶ Disposable head cap



1 | Disposable round polypropylene cap, 15gr.

Code	Ø approx.	Colour	Box
RBD001	30 cm	White	100 units
RBD002	30 cm	Green	100 units

RBD - Head cap

▶ Disposable head cap, fold



1 | Single-use polypropylene accordion cap

Code	Ø approx.	Colour	Box
RBD003	30 cm	White	100 units
RBD004	30 cm	Blue	100 units
RBD005	30 cm	Green	100 units
RBD006	30 cm	Rouge	100 units

RBG - Glasses



Anti-fog



Anti-fog Plus

A thin film of a few μm over the entire surface of the lens that prevents the formation of condensation by absorbing moisture droplets. It prevents the formation of fog, making it impossible to work safely.

The Anti-Fog Plus coating is permanent.



Anti-scratch



Anti-scratch Plus

Treatment applied to the surface of lenses which forms a barrier against surface damage. The lenses are protected against scratches and abrasions caused by rubbing against solid particles.

The Anti-scratch Plus treatment is permanent.

Regulations respected

STANDARD EN 166

REGULATION (EU) 2016/425

Field of use

3

Liquids (drops)

4

Dust with particle thickness > 5 μm

5

Gas, vapour, spray, smoke and dust with particle thickness < 5 μm

9

Splashes of molten metal and penetration of hot solids




▶ RBG - Glasses

▶ Protection mask 601



- 1 | Wide field of vision panoramic glasses. RX insert for corrective lenses.
- 2 | PPE category III (polycarbonate lens). Adjustable elastic strap and flexible frame.
- 3 | With ventilation systems: indirect and in the joint.

601 mask features



1. Anti-fog Plus



2. Anti-scratch Plus

Code	RBG001
Frame	Transparent - Waterproof
Field of use	3, 4, 5, 9
Lens	Transparent PC

▶ Protection mask 602



- 1 | Polycarbonate lens. Indirect ventilation system.
- 2 | Overlap with prescription glasses. Adjustable elastic band and soft frame.

Code	RBG002
Frame	Transparent
Field of use	3, 4
Lens	Transparent UC (untreated)
Lot	10


▶ RBG - Glasses


▶ Safety eyewear 505



- 1 | Neck strap included in the package.
- 2 | Extremely lightweight and wrap-around, soft non-slip rubber nose bridge
- 3 | Elastic structure, very flexible and resistant.

Safety eyewear 505 features


 = 22 g





2. Anti-scratch

Code	RBG005
Frame	Transparent
Lens	Transparent
Lot	10

Safety eyewear 505 features

 = 22 g





2. Anti-scratch 3. Anti-fog


Code	RBG006
Frame	Transparent
Lens	Transparent
Lot	10

▶ Safety eyewear 511 for small faces



- 1 | Adjustment of the pins in length and inclination.
- 2 | Certified version for small faces, with H marking according to EN166
- 3 | Frame made of anti-allergic material, built-in lateral protection.

Safety eyewear 511 features

 = 34 g



2. Anti-scratch

Code	RBG007
Frame	Blue
Lens	Transparent
Lot	10



▶ RBG - Glasses

▶ Safety eyewear 511



=34 g



- 1| Adjustment of the pins in length and inclination.
- 2| Frame made of anti-allergic material, built-in lateral protection.

Code	RBG009
Frame	Blue
Lens	Transparent
Lot	10

▶ Safety eyewear 511



=34 g



- 1| Adjustment of the pins in length and inclination.
- 2| Frame made of anti-allergic material, built-in lateral protection.

Code	RBG008
Frame	Black
Lens	Transparent
Lot	10

▶ Safety eyewear 519



=43 g




Code	RBG011
Frame	Black-grey
Lens	Transparent
Lot	10

▶ RBG - Glasses

▶ Safety eyewear 519



=43 g



- 1| Ventilation system practiced in the side protections.
- 2| Adjustable temples in length and inclination. Superposable to prescription glasses.

Code	RBG010
Frame	Black-grey
Lens	Transparent
Lot	10

▶ Safety eyewear 520



- 1| Ventilation system in the side shields.
- 2| Superimposable on prescription glasses. Made entirely from polycarbonate, with no metal parts.

Code	RBG012
Frame	Transparent
Lens	Transparent UC
Lot	10

▶ Safety eyewear 568



=23 g



- 1| Inexpensive glasses ideal for essential protection
- 2| Light, practical, made of polycarbonate.

Code	RBG015
Frame	Transparent
Lens	Transparent
Lot	10



▶ RBG - Glasses

▶ Safety eyewear 546



- 1 | Wrap-around lens with excellent side coverage
- 2 | Sideburns that guarantee an excellent fit and comfort.
- 3 | Non-slip rubber nose bridge. Very lightweight
- 4 | Protects against blue light up to 525 nm (UV525)

Safety eyewear 546 features

=25 g 1. Anti-fog plus 2. Anti-scratch Plus

Code	RBG014
Frame	Black-orange
Lens	Orange

Safety eyewear 546 features

=25 g 1. Anti-fog

Code	RBG022
Frame	Rouge-grey
Lens	Transparent

Safety eyewear 546 features

=25 g 1. Anti-scratch

Code	RBG033
Frame	Black-green
Lens	Transparent

▶ Safety eyewear 529



- 1 | Glasses with adjustable temples

=23 g 1. Anti-scratch

Code	RBG024
Frame	Black-green
Lens	Transparent
Lot	10

▶ RBG - Glasses

▶ Safety eyewear 505



- 1 | Extremely light and enveloping. Elastic, highly flexible and resistant structure. UV400 protection. New temples for improved ergonomics.
- 2 | Perfect fit and comfort thanks to flexible, adjustable nose bridge.
- 3 | Neck strap included in packaging.

=24 g 1. Anti-scratch 2. Anti-fog

Code	RBG029
Frame	Brown
Lens	Brown
Lot	10

▶ Face shield 699 with chin guard



- 1 | Interchangeable foam front headband. Top adjustment. 1mm thick polycarbonate face shield, foldable at 90°. Can be fit over protective glasses and face masks.
- 2 | Certified to EN166 and EN170 standards.

Code	Field of use
RBG017	3



▶ RBG - Glasses

▶ Medical lens model 701



1 | Class 1 optical quality lens with anti-fog and anti-scratch coating (marked U 1 S N). Adjustable elastic strap.

Code	RBG030
-------------	--------

▶ RBJ - Masks

▶ Surgical mask Type IIR



1 | Polypropylene mask, 3 layers and 3 folds. Single use.
 2 | Elastics to adjust to the ears and a mouldable strip for a better nasal support.
 3 | 100% latex and fibreglass free product.

Code	RBJ007
Size	U
Dimensions	17,5 x 9,5
Colour	Blue
Properties	
First layer	Polypropylene density 25 gr/m ³
Second layer	MeltBlown density 25 g/m ³
Third layer	Polypropylene density 20 g/m ³
Box	50 units
Lot	40

▶ RBK - Aprons

▶ Apron PS white



1 | Polyethylene (LDPE) apron, 0.02 mm thick, with head opening and fastening straps on the back.
 2 | For general purpose protection against minor hazards. Single use.
 3 | 100% Latex free product.

Code	RBK001
Size	U
Dimensions	120 X 78
Thickness	0,02
Colour	White
Box	100 units
Lot	10

▶ RBL - Gloves

▶ Silicon semi-glove



Code	Temperature max.	Size
RBL002	250 °C	Unique



Code	Temperature max.	Size
RBL001	260 °C	Unique



Code	Temperature max.	Size
RBL004	250 °C	Unique



RBL - Gloves

▶ Heat-resistant glove



1 | Thickness cotton: 850 g/m². For short term manipulations.

Code	Temperature max.	Size
RBL003	100 °C	Unique

▶ Gloves Thermal Kevlar-400



Code	Temperature max.	Size
RBL505	400 °C	Unique

▶ Nitrile gloves



1 | 100% Nitrile glove, powder-free. Highly elastic. Ambidextrous.

Code	Size	Colour	Box	Lot
RBL005	S 6/7	Cobalt Blue	100 units	10
RBL006	M 7/8	Cobalt Blue	100 units	10
RBL007	L 8/9	Cobalt Blue	100 units	10
RBL008	XL 9/10	Cobalt Blue	100 units	10

RBL - Gloves

▶ Neoprene gloves



1 | Chemical protection glove. BKLOPT - TYPE A.
2 | Flocked cotton. Dotted palm and fingers for better grip. Presented in a personalized individual pouch for improved storage.

Code	Length	Thickness	Size
RBL009	330 mm	0,68 mm	8
RBL010	330 mm	0,68 mm	9
RBL011	330 mm	0,68 mm	10

▶ P.V.C gloves. Smooth finish



1 | Chemical protection glove. AKLMPST - TYPE A.
2 | Sanitized® treatment for increased glove longevity and better hygiene for your hands. Interlock cotton support.

Code	Length	Size
RBL012	400 mm	9
RBL013	400 mm	10

▶ Special ICE gloves



1 | One-piece knitted assembly. 2-layer support: outer one-piece knitted polyamide and warm fleece lining.
2 | The unique formula of this coating allows for flexibility and dexterity at temperatures down to -50°C.

Code	Size
RBL014	8
RBL015	9
RBL016	10



RBN - Shoe cover

Disposables shoe cover



1 | PE plastic shoe covers. Ambidextrous. 100% Latex free product.

Code	Dimensions	Colour	Box	Lot
RBN001	410 x 155 mm	White	100 units	10
RBN002	410 x 155 mm	Green	100 units	10
RBN003	410 x 155 mm	Blue	100 units	10

RDB - Stations

Eye wash station



- 1 | Station made of high impact polystyrene and polyethylene bottle.
- 2 | Can be filled with distilled water or other liquid indicated for ophthalmological use.
- 3 | 500 mL bottle is not included (Ref. RDD001) Dimensions 300x300 mm

Code	Description
RDB001	French version
RDB002	English version
RDB003	German version

RDD - Eyewash bottles

Eye wash bottle, PE




- 1 | Polyethylene bottle. Dimensions 63x73x213 mm
- 2 | Eyecup with lid and head with automatic valve
- 3 | It has a discharge tube.

Code	Description	Lot
RDD001	Eye wash bottle M45	10

RDB - Stations


Combined wall station



- 1 | Supplied with mirror and pictogram with instructions for use.
- 2 | Wall station with:
 - 1 pH neutral solution: for eye injuries caused by acids and alkalis
 - 1 saline solution: for eye injuries caused by dust, dirt, splinters...

Code	Solutions	Use
RDB007	1 x 200 mL pH neutral 1 x 500 mL saline	1 eye

Combined wall station



- 1 | Supplied with mirror and pictogram with instructions for use.
- 2 | Wall station with:
 - 1 pH neutral solution: for eye injuries caused by acids and alkalis
 - 1 saline solution: for eye injuries caused by dust, dirt, splinters...

Code	Solutions	Use
RDB009	1 x 500 mL pH neutral Duo 1 x 1000 mL saline Duo	2 eyes

Combined wall station, dust-proof iBOX



- 1 | Supplied with mirror and pictogram with instructions for use.
- 2 | Wall station with:
 - 1 pH neutral solution: for eye injuries caused by acids and alkalis
 - 1 saline solution: for eye injuries caused by dust, dirt, splinters...

Code	Solutions	Use
RDB011	1 x 200 mL pH neutral 1 x 500 mL saline	1 eye
RDB012	1 x 500 mL pH neutral Duo 1 x 1000 mL saline Duo	2 eyes



▶ RDD - Eyewash bottles

▶ pH neutral solution



- 1 | 4.9% sterile phosphate solution
- 2 | Indicated to prevent eye injuries and minor skin injuries caused by acids and alkalis.
- 3 | 3 year shelf life, have expiration date.
- 4 | With clear and precise application instructions and pictograms. With CE marking in accordance with medical equipment regulations.

Code	Capacity	Washing T. approx.	Use
RDD006	200 mL	2 min	1 eye
RDD007	500 mL (Duo)	2 min	2 eyes
RDD008	1000 mL	10 min	1 eye

▶ Saline rinsing solution



- 1 | 0.9 % sterile sodium chloride solution.
- 2 | Indicated to prevent eye injuries caused by dust, dirt, metal or wood splinters, solvents or oils.
- 3 | 3 year shelf life, they have an expiry date.
- 4 | With clear and precise application instructions and pictograms. With CE marking in accordance with medical equipment regulations.

Code	Capacity	Washing T. approx.	Use
RDD002	200 mL	2 min	1 eye
RDD003	500 mL	2 min	1 eye
RDD004	1000 mL	10 min	1 eye
RDD005	1000 mL (duo)	10 min	2 eyes

▶ RDB - Stations

▶ Shower and eye washer



- 1 | Large-diameter artichoke, operated by a rigid handle.
- 2 | With pedestal type eye wash, equipped with two sprayers of high flow of low pressure aerated water and manual side opening.

Code	Material
RDB004	stainless steel 304

▶ Eye washer



- 1 | Pedestal type and ground-anchored base. Automatic shut-off valve.
- 2 | With two low-pressure aerated water sprinklers.

Code	Material
RDB006	Corrosion-resistant plastic tray



RDZ - Others

Belt bag for bottles



1 | Adaptable to the belt or work clothes. Easy to open in case of accident

Code	Use
RDZ001	For bottles of 200 mL

RZB - Table protectors

Safety shield



1 | Made of polymethyl methacrylate. High stability with non-slip feet.

Code	Dimensions	Protection	Thickness
RZB001	242x365 mm	Hazardous chemicals, aerosols coarse sprays	5 mm
RZB002	305x480 mm	Beta-emitting isotopes, including P32	9 mm

Protection paper



1 | For protection of work surfaces and the user in case of spillage of aggressive liquids.

2 | Weight: 140 g/m²; Thickness: 250 μm.

Code	Dimensions	N° sheets
RZB013	460x570 mm	50
RZB014	460x570 mm	100

RZB - Table protectors

Protection film



1 | Protective film with self-adhesive layers covered with a 25 micron FEP (fluorinated copolymer).

Code	Material	Width	Length
RZB012	Aluminium	635 m	4,6 mm

RZD - Wound treatment

Plaster dispenser



1 | A small, modern plaster dispenser with a unique design.

2 | QuickFix® UNO is available in a glossy white or with a transparent cover. It can be screwed to the wall or fixed with adhesive plasters.

3 | Supplied without strips. This allows you to fill the dispenser with the strips of your choice, according to your needs.

Code	Colour	Dimensions
RZD001	Transparent	130x85x35 mm
RZD005	White	130x85x35 mm

NEW PRODUCT



Disposable sleeve



1 | Polyethylene. Adjustable elastic bands

Code	Colour	Pack
RBM001	Blue	C/100



NEW CENTRIFUGE RANGE NAHITA-BLUE:

INNOVATION AND
EFFICIENCY IN
THE LABORATORY



www.auxilab.es/en/



correo@auxilab.es



+34 693 941 198