


DRAINING RACKS

As essential elements in the daily lab work, table and wall draining racks are used to drain and dry all kind of laboratory material (glassware, plasticware, etc.) of small and medium size (beakers, measuring cylinders, flasks, tubes, etc.). Available both table and wall models.



ref. 63224730



ref. 63224750

Table models

- 1 | Made of plastified wire.
- 2 | With silicone rubber protectors and extractable tray in the base.
- 3 | They can stand on a table surface or can be hanged on the wall.

Reference	Pins	Dimensions (mm)	Type
63224730	32	470 x 360 x 160	Table / Wall
63224750	55	640 x 360 x 140	Table / Wall

Accessories

Reference	Description
63224759	Tray for draining rack



Wall model

- 1 | Made of polypropylene.
- 2 | Pins can be easily detached allowing adapting to material with special size or shape.
- 3 | Provided with drainpipe for emptying the lower tray.

Reference	Pins	Dimensions (mm)	Type
41000001	20	500x350x100	Wall