

**Escobillones | Brushes**

 Marca|Trademark:**Nahita**

Le ofrecemos una amplia gama de escobillones de laboratorio para la limpieza a fondo de tubos de ensayo, pipetas, probetas, buretas y todo tipo de frascos. Fabricados en alambre de acero zincado y con cerdas naturales, son ideales para mantener una correcta limpieza del material de vidrio de laboratorio asegurando así siempre un buen resultado en sus experimentos. Todos los escobillones presentan anilla para colgar en el extremo del mango.

We offer a wide range of laboratory brushes for cleaning assay tubes, pipettes, cylinders, burettes and any kind of flasks. Made of zinc-plated steel wire and natural bristle, they are ideal for properly cleaning glassware to assure a good performance of your experiments. All the brushes present a ring to be hanged at the end of the handle.



Referencia Code	Long. total Total length	Pelo (ØxL) Brush (ØxL)	Plumero Tuft	Para For
63506020	300 mm	10x70 mm	No	Tubos hemólisis / Hemolysis tubes
63506030	360 mm	15x70 mm	Sí / Yes	Tubos pequeños / Small tubes
63506040	350 mm	25x75 mm	Sí / Yes	Tubos medios / Medium tubes
63506050	350 mm	30x85 mm	Sí / Yes	Tubos grandes / Large tubes
63506070	380 mm	40x110 mm	No	Frascos y probetas, 250 mL / Bottles & cylinders, 250 mL
63506080	380 mm	50x130 mm	No	Frascos y probetas, 500 mL / Bottles & cylinders, 500 mL
63506090	440 mm	56x140 mm	No	Frascos y probetas, 1 L / Bottles & cylinders, 1 L
63506910	440 mm	80x160 mm	No	Frascos y probetas, 2 L / Bottles & cylinders, 2 L
63506150	500 mm	10x100 mm	No	Pipetas / Pipettes
63506100	850 mm	10x180 mm	Sí / Yes	Buretas 10 mL / Burettes 10 mL
63505110	870 mm	15x260 mm	Sí / Yes	Buretas 25 mL / Burettes 25 mL
63506130	450 mm	80x150 mm	No	Vasos / Glasses
63506590	210 mm	15x70 mm	No	Traqueal / Tracheal
63506960	270 mm	5/10x80 mm	No	Butirómetros / Butyrometers