

## EMBUDOS BÜCHNER

Büchner funnels are used in vacuum filtration processes and present a perforated plate ideal for placing a filter membrane. Made of high chemical resistant porcelain, Nahita Büchner funnels are glazed inside and outside and are available in different sizes to adapt to any user necessities.



Reference	Capacity	Total height	Inner Ø	Stem mean Ø
20237040	30 mL	70 mm	40 mm	13 mm
20237050	50 mL	80 mm	50 mm	15 mm
20237070	120 mL	125 mm	70 mm	20 mm
20237090	200 mL	145 mm	95 mm	21 mm
20237110	400 mL	165 mm	110 mm	26 mm
20237120	600 mL	170 mm	125 mm	30 mm
20237150	1000 mL	195 mm	135 mm	33 mm

\*Approx. dimensions.