

Como su propio nombre indica, los matraces de destilación se utilizan en aquellas reacciones que requieren la evaporación de líquidos y su posterior condensación. Para ello, los matraces de destilación están fabricados en vidrio borosilicato de gran resistencia al choque térmico y presentan una base esférica para el calentamiento uniforme de las soluciones a destilar. Disponen también de un cuello largo y estrecho y de un brazo lateral para conectar con un refrigerante y permitir así la condensación del líquido evaporado.

Distilling flasks are used in those reactions requiring evaporation and a later condensation of liquids. Because of this, distilling flasks are made of thermal shock resistant borosilicate glass and present a spherical base for a uniform heating of the samples to be distilled. They are also provided with a long and narrow neck and a side arm to be connected to a condenser, thus, allowing the condensation of the previously evaporated liquid.



Referencia Code	Capacidad Capacity	Ø Boca Neck Ø	Altura Height	Ø Brazo Arm Ø	Long. Brazo Arm length
10041210	100 ml	15 mm	174 mm	7 mm	100 mm
10041220	250 ml	19 mm	216 mm	8 mm	120 mm
10041250	500 ml	24 mm	275 mm	8 mm	150 mm
10041310	1000 ml	28 mm	320 mm	9 mm	180 mm