



Endo Glassware offers this series of short stem funnels which are ideal for transferring liquids into narrow neck flasks or for filtering any type of solution. Made of soda glass or borosilicate glass, they are available in different sizes.

► Borosilicate glass funnels



- 1 | Highly resistant to acids or bases and thermal shock. .
- 2 | For transferring liquids or filtering solution.
- 3 | 60° cone angle and borosilicate glass. .

Code	Cap.	Funnel Ø	Stem Ø	Height	Filter Ø
PNB001	10 ml	40 mm	4 mm	70 mm	50-60 mm
PNB002	30 ml	50 mm	7 mm	95 mm	50-60 mm
PNB003	100 ml	75 mm	9 mm	140 mm	90-110 mm
PNB004	150 ml	90 mm	10 mm	170 mm	110-125 mm
PNB005	200 ml	100 mm	11 mm	185 mm	110-125 mm

► Soda glass funnels



- 1 | Presentan una gran resistencia a la acción de ácidos y bases y al choque térmico.
- 2 | For transferring liquids or filtering solution.
- 3 | 60° cone angle and soda glass. .

Code	Cap.	Funnel Ø	Stem Ø	Height	Filter Ø
PNB006	400 ml	120 mm	12 mm	220 mm	125-150 mm
PNB007	750 ml	150 mm	15 mm	300 mm	150-180 mm