

NARROW BOTTLES PE WITHOUT CAP



- 1 | Round bottle with wide neck without cap
- 2 | Made in polyethylene (PE)
- 3 | Natural color

Code	Capacity	Ø ext. (total)	Ø int. (neck)	Heigh	GL
45315801	50 mL	40 mm	13 mm	84 mm	18
45315821	100 mL	45 mm	13 mm	102 mm	18
45315841	250 mL	60 mm	18 mm	138 mm	25
45315861	500 mL	75 mm	18 mm	169 mm	25
45315881	1000 mL	95 mm	21 mm	212 mm	28
45315901	2000 mL	120 mm	21 mm	260 mm	28

STANDARDIZED SCREW CAPS (PE)

■ Screw cap




Reference	GL
45315800	18
45315840	25
45315880	28



NARROW BOTTLES PE WITHOUT CAP

■ Wash bottle screw cap



Reference	GL
45303300	18
45303301	25
45303302	28



■ Short spout dropping cap



Reference	GL
45303303	18
45303304	25



■ Long spout dropping cap



Reference	GL
45303306	18
45303307	25
45303308	28

