

Termómetros de cristal líquido | Liquid crystal thermometers

Marca|Trademark: Nahita

Este tipo de termómetros basa su funcionamiento en la propiedad de los cristales líquidos de cambiar secuencialmente de color con la variación de la temperatura; así, pasan del rojo/amarillo al verde y posteriormente al azul, indicando el color verde la temperatura real. El formato de estos termómetros es el de una tira flexible auto-adhesiva que se adapta incluso a superficies curvas. De fácil y rápida lectura, presentan la ventaja de una duración de uso ilimitada en condiciones normales y de ser totalmente inocuos y carentes de peligro. Están disponibles en diversos rangos de lectura para adaptarse a cualquier necesidad del usuario.

This type of thermometers is based on the property of liquid crystals to change colour sequentially with temperature variations. They change from red/yellow to green and then to blue colour; real temperature is indicated in green. The format of these thermometers is a self-adhesive flexible strip that adapts even to curved surfaces. Very easy to read and use, they present the advantage of a limitless use in normal conditions and that of being totally innocuous. They are available in different ranges to adapt to any user's necessities.



Referencia / Code	Rango / Range	Medidas / Size
72501013	-3/14 °C; 27/57 °F	130x18 mm
72501031	14/31 °C; 57/88 °F	130x18 mm
72501049	32/49 °C; 90/120 °F	130x18 mm
72502004	-3/14 °C	58x18 mm
72502013	6/14 °C	58x18 mm
72502022	14/22 °C	58x18 mm
72502031	24/32 °C	58x18 mm
72502040	32/40 °C	58x18 mm
72502049	41/49 °C	58x18 mm