


Kitasatos para vacío | **Filtering**

Los kitasatos son unos matraces robustos de paredes gruesas y fabricados en vidrio borosilicato de excelente resistencia física y química. De uso extendido en laboratorios de química y biología, poseen una oliva lateral para ser conectados a una bomba de vacío o trompa de agua y pueden ser utilizados con diversos tipos de embudos para la filtración rápida de soluciones. Presentan además graduación aproximada y amplia superficie de rotulación en blanco.

Marca|**Trademark:** Nahita

Filtering flasks are sturdily constructed with thick walls and are made of borosilicate glass of an excellent physical and chemical resistance. Widely used in chemistry and biology laboratories, they present a lateral hose connection to be connected with some sort of vacuum assembly and they can be used with any kind of funnel for fast solution filtration. They also feature a large marking area and approximate graduation.



Referencia	Capacidad	Ø Base
Code	Capacity	Base Ø
12418221	250 mL	85 mm
12418251	500 mL	105 mm
12418311	1000 mL	140 mm