



1| Endo Glassware offers a wide range of multipurpose bottles of different capacities and ideal for storing any type of solution.

2| All bottles are available in transparent (white) glass or amber glass (for those photosensitive solutions or reagents).

#### ▶ Wide neck bottles



With standard ground joint and polyethylene stopper.

Code	Capacity	Base Ø	Neck Ø	Height	Glass
CHB001	100 ml	54 mm	95 mm	29/22	White
CHB005	100 ml	54 mm	95 mm	29/22	Amber
CHB002	250 ml	74 mm	125 mm	34/23	White
CHB006	250 ml	74 mm	125 mm	34/23	Amber
CHB003	500 ml	88 mm	165 mm	45/27	White
CHB007	500 ml	88 mm	165 mm	45/27	Amber
CHB004	1000 ml	112 mm	200 mm	60/46	White
CHB008	1000 ml	112 mm	200 mm	60/46	Amber

#### ▶ Narrow neck bottles



With standard ground joint and polyethylene stopper.

Code	Capacity	Base Ø	Neck Ø	Height	Glass
CHG001	100 ml	54 mm	95 mm	14/15	White
CHG005	100 ml	54 mm	95 mm	14/15	Amber
CHG002	250 ml	74 mm	125 mm	19/26	White
CHG006	250 ml	74 mm	125 mm	19/26	Amber
CHG003	500 ml	88 mm	165 mm	24/29	White
CHG007	500 ml	88 mm	165 mm	24/29	Amber
CHG004	1000 ml	112 mm	200 mm	29/32	White
CHG008	1000 ml	112 mm	200 mm	29/32	Amber