


▶ BPG - Slide storage

▶ Storage box for slides



- 1 | Made of red plastic, stackable.
- 2 | With 6 supports for 75 slides (total capacity of 450 slides).

Code	Description	Dimensions
BPG017	Storage box for 450 slides	34x23x10 cm

▶ Microscope slides storing system

- 1 | It is designed for classification and storage of microscope slides.
- 2 | It is composed of 4 sliding drawers, a closing cover and four fixed support feet or, as alternative, four swivel castors to guarantee mobility.



- 1 | Each single drawer can hold up to 13 supports for slides (ref. BPG026) and transparent cover plates (ref. BPG025), to be ordered separately.
- 2 | Each slides support BPG026 can contain 100 microscope slides, easily extractable. Each cavity is numbered. The complete Storing System can therefore contain up to 5,200 slides.
- 3 | Transparent lids BPG025 guarantee an easy check of the contents, protecting the slides from dust and other foreign bodies.

Code	Description
BPG025	Lid
BPG026	Slide support



- 1 | The metal structure makes it stable, while ensuring the value of lightness, functionality, strength and durability.
- 2 | It is made of impact-resistant PMMA and steel.

Code	Wheels	Dimensions
BPG035	Yes	490x475x780 mm
BPG036	No	490x475x700 mm



- 1 | Mobile version could be transformed into fixed and vice versa, using four fixed foot kit, or four wheel kit.

Code	Description
BPG037	Kit of four swivel caster for slide storing system
BPG038	Kit of four feet for slide storing system