



▶ Injection vials 2 mL, crimp closure



- 1 | Made of type 1 clear glass according to European Pharmacopeia.
- 2 | Dimensions 12x32 mm.
- 3 | Volume 2 mL.

Code	Glass	Closure	Mouth	Graduation
CMF001	Clear	Crimp	DN11	No
CMF002	Amber	Crimp	DN11	No
CMF003	Clear	Crimp	DN11	Yes

▶ Aluminum caps for DN11 vials



Code	Mouth	Specifications
CML002	DN11	Joint PTFE/silicone pre-slit


▶ Injection vials 2 mL, screw cap



- 1 | Made of type 1 clear glass according to European Pharmacopeia.
- 2 | Dimensions 12x32 mm.
- 3 | Volume 2 mL.

Code	Glass	Closure	Mouth	Graduation
CMF009	Clear	Screw	DN8	No
CMF011	Clear	Screw	DN8	Yes
CMF012	Amber	Screw	DN8	Yes

▶ PP screw caps for vials



Code	Mouth	Specifications
CML003	DN9-425	Septum PTFE/rubber
CML006	DN8-425	Septum PTFE/rubber
CML007	DN8-425	Septum PTFE/silicone
CML008	DN8-425	Septum PTFE/silicone, pre-slit