



▶ Vortex stirrer, RSLAB 6PRO

- 1 | Suitable for touch (low speed) or continuous operation (high speed).
- 2 | Adapters securely clicked onto appliance in any position.
- 3 | Sturdy cast with aluminium base for excellent stability.
- 4 | Specially designed rubber feet for excellent vibration damping.
- 5 | Eccentric with oil-less ball bearings.
- 6 | Protection class according to DIN EN60529: IP21.
- 7 | Provided with a standard top for tubes with diameter <30 mm. (code LGG016)



Code	LGB003
Model	RSLAB-6PRO
Speed range	0-2500 rpm
Shaking movement	Orbital (diameter 4mm)
Dimensions	127x130x160 mm
Weight	3,5 Kg
Protection class	IP21

▶ Accessories (Not included)



LGG009 Universal top for tube adapters. Continuous operation.



LGG014 Adapter for 8 round bottom tubes 20 mm



LGG010 Adapter for 48 PCR microtubes 6 mm.



LGG015 Tube holding rod from code LGG010 to LGG014 Touch operation.



LGG011 Adapter for 15 round bottom tubes 10 mm



LGG016 Standard top for tubes diameter 30 mm



LGG012 Adapter for 12 round bottom tubes 12 mm.



LGG018 Antiskid platform for tubes max. diameter 99 mm



LGG013 Adapter for 8 round bottom tubes 16 mm



LGG017 Rubber feet