



Filtration assembly, 300 mL - Endo Glassware filtration assembly is ideal for filtering any type of solution to particle contamination analysis and also for filtration of HPLC solvents.

- 1 | Glass pieces are made of borosilicate glass and present 40/35 ground joints while the clamp is made of anodized aluminium.
- 2 | Vacuum connection, integrated into the fritted glass support, and the flask mouth are over the filtered liquid exit, thus avoiding leaking of the filtered liquid into the vacuum tube.

3 | The fritted glass disk presents a diameter of 40 mm and the funnel is graduated

4 | The different components of the assembly are available individually for replacement of damaged pieces.



Code	Description
PJB001	Complete filtration assembly, 300 mL

Composed of:

Code	Description
PBH030	Erlenmeyer flask ground joint 40/35, 1000 mL
PJD003	Body with fritted glass disk
PJD002	Filtration funnel, 300 mL
PJD006	Anodized aluminium clamp