

► Graduated cone-shaped glass

1 | Endo Glassware measuring cups are conical in shape with amber graduation.

2 | They are used for measuring and transferring approximate volumes of liquid.

3 | The base is wide and stable, with a spout for convenient pouring.

4 | They are available in various capacities to suit the needs of each user.



CODE: CBD001



CODE: CBD002



CODE: CBD003

Code	Capacity	Diameter	Height
CBD001	50 ml	45 mm	145 mm
CBD002	1000 ml	130 mm	330 mm
CBD003	2000 ml	160 mm	370 mm

BEAKERS

[C] BEAKERS, FLASKS AND OTHER CONTAINERS

Note: the manufacturer reserves the right to modify the characteristics and photographs of the products