



CGG - Narrow mouth, glass

Bottle icons



Capacity



Total capacity



Height



Base




Thread



Weight

DIN 28, features



- 1 | Available in transparent and amber glass
- 2 | Complies with specifications of European and United States Pharmacopeia
- 3 | Guaranteed product Quality control system with physical and chemical analysis: verification of the dimensions, thermal shock resistance, hydrolytic resistance and pressure resistance.
- 4 | Made in soda glass type III.
- 5 | Thread DIN 28 (see available stoppers)
- 6 | Bottles and stoppers with batch number and analysis certificate

Bottles with DIN-28 threaded mouth, transparent glass



Code	Colour						
CGG001	Transp.	30 mL	35 mL	67 mm	35,0 mm	DIN 28	40
CGG002	Transp.	60 mL	67,5 mL	94 mm	38,6 mm	DIN 28	60
CGG003	Transp.	125 mL	133 mL	114,4 mm	49,0 mm	DIN 28	90
CGG004	Transp.	250 mL	260 mL	142,9 mm	61,7 mm	DIN 28	140
CGG005	Transp.	500 mL	512 mL	176,8 mm	76,8 mm	DIN 28	240
CGG006	Transp.	1000 mL	1020 mL	215,6 mm	97,5 mm	DIN 28	450

Bottles with DIN-28 threaded mouth, amber glass



Code	Colour						
CGG007	Amber	30 mL	35 mL	67 mm	35,0 mm	DIN 28	40
CGG008	Amber	60 mL	67,5 mL	94 mm	38,6 mm	DIN 28	60
CGG009	Amber	100 mL	108 mL	104,8 mm	46,1 mm	DIN 28	72
CGG010	Amber	125 mL	133 mL	114,4 mm	49,0 mm	DIN 28	90
CGG011	Amber	250 mL	260 mL	142,9 mm	61,7 mm	DIN 28	140
CGG012	Amber	500 mL	512 mL	176,8 mm	76,8 mm	DIN 28	240
CGG013	Amber	1000 mL	1020 mL	215,6 mm	97,5 mm	DIN 28	450



▶ CGL - Stoppers for DIN-28 bottles

▶ Safety thread stopper, for DIN-28 bottles



1 | Safety stopper with seal, made of polyethylene.

Code	Color tapón	Color precinto
CGL009	Blanco	Blanco

▶ Tapón roscado, for DIN-28 bottles



1 | Stoppers with seal, made of polyethylene

Code	Stopper colour	Seal colour
CGL010	Black	Black
CGL011	White	Red
CGL012	Blue	White
CGL013	Red	White
CGL014	Green	White
CGL015	Yellow	White

▶ Dropper stopper with plastic tube, for DIN-28 bottles



1 | Polyethylene stopper with seal, plastic teat and tube. The tube can be cut to adapt to tubes from 1000 ml to 30 ml.

Code	Description
CGL019	Stopper, teat and tube

▶ Dropper stopper with glass tube, for DIN-28 bottles



1 | Polyethylene stopper with seal, latex teat and glass tube

Code	Description
CGL020	Stopper, teat and tube (65 mm) for 30 mL bottles
CGL021	Stopper, teat and tube (90 mm) for 60 mL bottles
CGL022	Stopper, teat and tube (110 mm) for 125 mL bottles