



## Sterile stoppers



### Features

Designed for microbiology and cell culture applications.  
 Made from very flexible cellulose, tapered to fit most sizes of: Erlenmeyer, matrix flasks and tubes.  
 Sterile to 200 °C.  
 Permeable to air.  
 Transpirable.  
 Hermetically adjustable by pressing on the top of the glass vessel.


### RANGE OF SIZES

For Ø of 9 to 11,5 mm. Part No. **1001478**  
 For Ø of 12,5 to 14,5 mm. Part No. **1001479**  
 For Ø of 13,5 to 15,5 mm. Part No. **1001480**  
 For Ø of 17,5 to 18,5 mm. Part No. **1001481**  
 For Ø of 18 to 22 mm. Part No. **1001482**  
 For Ø of 25 to 27 mm. Part No. **1001483**  
 For Ø of 28 to 31 mm. Part No. **1001484**  
 For Ø of 31 to 35 mm. Part No. **1001485**  
 For Ø of 36 to 42 mm. Part No. **1001486**  
 For Ø of 44 to 48 mm. Part No. **1001487**



## Caps for biological culture tubes "Sero-Tap"

Made of anodised aluminium with internal AISI 304 stainless steel spring inside the cap.  
 Not hermetically sealed. Autoclavable at 1 bar.



Made of anodised aluminium with internal AISI 304 stainless steel spring inside the cap. Not hermetically sealed. Autoclavable at 1 bar.

Part No.	For tubes Ø 12 - 13 mm	Part No.	For tubes Ø 19 - 20 mm
<b>1012131</b>	Colour Blue	<b>1019201</b>	Colour Blue
<b>1012132</b>	Colour Silver	<b>1019202</b>	Colour Silver
<b>1012133</b>	Colour Gold	<b>1019203</b>	Colour Gold
<b>1012134</b>	Colour Garnet	<b>1019204</b>	Colour Garnet
<b>1012135</b>	Colour Black	<b>1019205</b>	Colour Black
Part No.	For tubes Ø 15 - 16 mm	Part No.	For tubes Ø 24 - 26 mm
<b>1015161</b>	Colour Blue	<b>1024261</b>	Colour Blue
<b>1015162</b>	Colour Silver	<b>1024262</b>	Colour Silver
<b>1015163</b>	Colour Gold	<b>1024263</b>	Colour Gold
<b>1015164</b>	Colour Garnet	<b>1024264</b>	Colour Garnet
<b>1015165</b>	Colour Black	<b>1024265</b>	Colour Black
Part No.	For tubes Ø 17 - 18 mm	Part No.	For tubes Ø 28 - 30 mm
<b>1017181</b>	Colour Blue	<b>1028301</b>	Colour Blue
<b>1017182</b>	Colour Silver	<b>1028302</b>	Colour Silver
<b>1017183</b>	Colour Gold	<b>1028303</b>	Colour Gold
<b>1017184</b>	Colour Garnet	<b>1028304</b>	Colour Garnet
<b>1017185</b>	Colour Black	<b>1028305</b>	Colour Black



INTERIOR HOLDING CLIP.

