



Digital, compact and functional centrifuges that have adopted a new vibration absorber contributing to a noise and vibration reduction



- 1 | Fixed- angle rotor.
- 2 | Digital speed adjustment.
- 3 | Operation during time periods up to 99 min.
- 4 | Safety switch that avoids centrifugation if the lid is

open and that automatically stops operation in case of opening the lid.

- 5 | Maintenance-free induction motor.

Code	GBF006	GBF007
Model	2698	2698/5
Capacity	8x15 mL	16x5 mL
Máx. RCF		2000 g
Máx. speed		4000 rpm
Timer		0-99 min
Dimensions		32x32x28 cm
Weight		15 Kg
Power supply		220 V / 50 Hz

SINGLE ROTOR [G] CENTRIFUGES

Note: manufacturer reserves the right to modify features and photographs of the goods