



▶ Narrow mouth, glass

Bottle icons



Capacity



Total capacity



Height



Base




Thread



Weight

DIN 18, features



- 1 | They meet the specifications of the European and US Pharmacopeia
- 2 | Made of soda-lime glass type III
- 1 | Available in transparent, amber glass, cobalt blue and green
- 2 | Guaranteed quality control system of the product with physical and chemical tests: verification of the dimensions and the thermal shock, hydrolytic and pressure resistances.
- 3 | Thread DIN 18

▶ Bottles with DIN18 threaded mouth, transparent glass



Code	Colour						
CGG014	Transp.	5 mL	6,7 mL	50,2 mm	22,0 mm	DIN 18	22,0
CGG015	Transp.	10 mL	12,5 mL	58,0 mm	24,8 mm	DIN 18	27,5
CGG016	Transp.	15 mL	20,0 mL	65,0 mm	28,5 mm	DIN 18	35,0
CGG017	Transp.	30 mL	35,0 mL	78,8 mm	33,0 mm	DIN 18	48,0
CGG018	Transp.	50 mL	57,0 mL	92,2 mm	37,2 mm	DIN 18	55,0
CGG019	Transp.	100 mL	105,0 mL	112,0 mm	44,5 mm	DIN 18	88,0

▶ Bottles with DIN18 threaded mouth, amber glass



Code	Colour						
CGG020	Amber	5 mL	6,7 mL	50,2 mm	22,0 mm	DIN 18	22,0
CGG021	Amber	10 mL	12,5 mL	58,0 mm	24,8 mm	DIN 18	27,5
CGG022	Amber	15 mL	20,0 mL	65,0 mm	28,5 mm	DIN 18	35,0
CGG023	Amber	30 mL	35,0 mL	78,8 mm	33,0 mm	DIN 18	48,0
CGG024	Amber	50 mL	57,0 mL	92,2 mm	37,2 mm	DIN 18	55,0
CGG025	Amber	100 mL	105,0 mL	112,0 mm	44,5 mm	DIN 18	88,0



Bottle icons



Capacity



Total capacity



Height



Base



Thread



Weight

▶ Bottles with DIN18 threaded mouth, cobalt blue glass



Code	Colour						
CGG026	Cobalt blue	5 mL	6,7 mL	50,2 mm	22,0 mm	DIN 18	22,0
CGG027	Cobalt blue	10 mL	12,5 mL	58,0 mm	24,8 mm	DIN 18	27,5
CGG028	Cobalt blue	15 mL	20,0 mL	65,0 mm	28,5 mm	DIN 18	35,0
CGG029	Cobalt blue	30 mL	35,0 mL	78,8 mm	33,0 mm	DIN 18	48,0
CGG053	Cobalt blue	50 mL	57,0 mL	92,2 mm	37,0 mm	DIN 18	55,0
CGG054	Cobalt blue	100 mL	105,0 mL	112,0 mm	44,0 mm	DIN 18	88,0

▶ Bottles with DIN18 threaded mouth, green glass



Code	Colour						
CGG030	Green	5 mL	6,7 mL	50,2 mm	22,0 mm	DIN 18	22,0
CGG031	Green	10 mL	12,5 mL	58,0 mm	24,8 mm	DIN 18	27,5
CGG032	Green	15 mL	20,0 mL	65,0 mm	28,5 mm	DIN 18	35,0
CGG033	Green	30 mL	35,0 mL	78,8 mm	33,0 mm	DIN 18	48,0



CGL - Stoppers for glass bottles

▶ Dropper stopper with glass tube, DIN 18



1 | Polyethylene stopper with seal, teat (latex free) and glass tube.

2 | Black colour



Code	Colour	Description
CGL024	Black	Dropper stopper for 5 mL bottles
CGL025	Black	Dropper stopper for 10 mL bottles
CGL026	Black	Dropper stopper for 15 mL bottles
CGL027	Black	Dropper stopper for 30 mL bottles
CGL028	Black	Dropper stopper for 50 mL bottles
CGL029	Black	Dropper stopper for 100 mL bottles



▶ CGL - Stoppers for glass bottles

▶ Dropper stopper with glass tube, DIN 18



- 1 | Polyethylene stopper with seal, teat (latex free) and glass tube.
- 2 | White colour

Code	Colour	Description
CGL003	White	Dropper stopper for 5 mL bottles
CGL004	White	Dropper stopper for 10 mL bottles
CGL005	White	Dropper stopper for 15 mL bottles
CGL006	White	Dropper stopper for 30 mL bottles
CGL007	White	Dropper stopper for 50 mL bottles
CGL008	White	Dropper stopper for 100 mL bottles

▶ Stopper DIN 18



- 1 | Kid safety stopper with with seal and inner dropper. (code CGL002)
- 2 | Made of polyethylene and available in yellow or black.

Code	Description	Colour
CGL002	Safety stopper with inner dropper	Yellow
CGL030	Stopper DIN18	Black