



Automatic pipettes, eco series

- 1 | The automatic pipettes Eco series from Digipette, especially useful for educational purposes, present a high ease of use and durability at a really economic price.
- 2 | The ergonomic design of the pipette adapts perfectly to the hand to avoid fatigue after long periods of use and the two button operation is comfortable and requires only one hand. They present a button for aspirating and dispensing the volume and

another button for tip ejection avoiding the manipulation of the tip with the hand and possible contaminations.

- 3 | Each model is colour coded to easily identify the pipette and to indicate the type of tip that must be used with each model.

- 4 | The volume adjustment is performed with only one hand and is shown in a 3-digit display.



Code	Capacity	Accuracy (%)	Precision (%)	Tips
MGB020	1.0 µL	± 4.0	≤ 2.0	MJB011
	10.0 µL	± 1.2	≤ 7.0	
MGB022	20.0 µL	± 3.7	< 1.2	MJJ016
	100.0 µL	± 1.9	< 0.4	
MGB025	1000.0 µL	± 1.5	≤ 0.6	MJN015
	5000.0 µL	± 0.8	≤ 0.3	

Accessories - Tips for micropipettes



- 1 | Made of polypropylene and fully autoclavable.

Code	Capacity	Colour	Units/bag
MJB011	0.1-10 µL	White	1000
MJN016	2-200 µL	Yellow	1000
MJN015	100-1000 µL	Blue	500
MJQ009	1000-5000 µL	White	300