



## Roll and tilt stirrer 683



### Features

- 1| Is ideal for mixing blood samples, viscous substances and liquid solid suspensions where minimum aeration is required. It also constitute a powerful and silent unit which can run continuously for many hours unattended.
- 2| It presents 7 rollers with simultaneous rolling and rocking motion

getting sinusoidal wave effect that ensures a thorough mixing motion.

- 3| Moreover adjustable speed between 10 and 80 rpm and the possibility of accommodate various types of tubes make the Nahita roll and tilt stirrer perfectly adaptable to user necessities.

### Technical features

Code	LJD001
Speed range	10-80 rpm
Number of rollers	7
Roller length	240 mm
Power	220 V / 50 Hz
Dimensions	400x250x100 mm
Weight	4 Kg