


▶ FCB - Bossheads

- 1 | Laboratory bossheads are used to other lab pieces such as retort rings, clamps, rods, etc. to support frames.
- 2 | Made of metal and very durable, Nahita laboratory bossheads are available in different models and sizes to adapt to any necessity.

▶ Double bossheads



- 1 | Metallic and available in different sizes and colours the allow assembling pieces in right angle.
- 2 | With adjusting wheels made of plastic.

Code	Description	Length	Rods up to
FCP001	Painted small	53 mm	13 mm
FCP002	Painted medium	70 mm	16 mm
FCP003	Painted large	84 mm	21 mm
FCP004	Painted extra large	97 mm	25 mm
FCP005	Zinc plated small	53 mm	13 mm
FCP006	Zinc plated medium	70 mm	16 mm
FCP007	Zinc plated large	84 mm	21 mm
FCP008	Zinc plated extra large	97 mm	25 mm
FCP009	Zinc plated rounded	100 mm	16 mm

CODE:FCP001- FCP004  
Painted

CODE:FCP005- FCP008  
Zinc



CODE:FCP009  
Zinc rounded

▶ Double, resin coated



Code	Rods up to
FCP010	16 mm

▶ Double rotary




Code	Rods up to
FCP011	20/13 mm
FCP012	18 mm

▶ FCB - Bossheads

▶ Bosshead model Phywe

- 1 | Metallic and large, its strong format is ideal for heavy loads.
- 2 | To hold two rods up to 24 mm diameter at a right angle.

Code	Colour	Length	Rods up to
FCP013	Gris	170 mm	24 mm



▶ Universal bossheads

- 1 | Metallic and very robust they allow cross, parallel and "T" arrangements.
- 2 | The wheels to adjust the screw are made of plastic.

Code	Colour	Length	Rods up to
FCP014	Chrome plated	55 mm	16 mm
FCP015	Painted	55 mm	16 mm



▶ Bossheads for laboratory frameworks

- 1 | Bossheads type P allow the assembly of two rods in right angle while the Trevin bosshead can be used to mount multiple supports in three different angles.
- 2 | Both types present hexagonal hollow screw for Allen wrenches.

Code	Modelo	Rods up to
FCP017	Typo P	12 mm
FCP016	Typo P	10 mm
FCP018	Typo P	14 mm
FCP019	Typo P	16 mm
FCP020	Trevin	9-13 mm

