



► Spherical flasks


- 1 | We offer this wide range of spherical flasks made of thermal shock resistant borosilicate glass.
- 2 | With narrow neck and a white marking area, these flasks are ideal for distillation processes or those reactions in which a uniform sample heating is necessary

- 3 | They are available in different capacities to cover all users' necessities.

► Flat bottom



Flat bottom FLASK, glass




Code	Cap.	d1.	d2.	h.
CDG001	50 mL	50 mm	22 mm	105 mm
CDG002	100 mL	65 mm	22 mm	108 mm
CDG003	250 mL	85 mm	34 mm	140 mm
CDG004	500 mL	105 mm	34 mm	180 mm
CDG005	1000 mL	135 mm	43 mm	210 mm
CDG006	2000 mL	160 mm	52 mm	250 mm
CDG007	3000 mL	188 mm	59 mm	282 mm

BORO
3.3

► Round bottom



Round bottom FLASK, glass



Code	Cap.	d1.	d2.	h.
CDJ001	50 mL	50 mm	22 mm	105 mm
CDJ002	100 mL	65 mm	22 mm	108 mm
CDJ003	250 mL	85 mm	34 mm	140 mm
CDJ004	500 mL	105 mm	34 mm	180 mm
CDJ005	1000 mL	135 mm	43 mm	210 mm
CDJ006	2000 mL	175 mm	50 mm	270 mm
CDJ007	5000 mL	230 mm	56 mm	345 mm

BORO
3.3