



▶ PBC - Soxhlet extractor

1| Soxhlet extraction systems are glass equipments used to extract soluble components from a solid sample into an organic solvent.  
2| Within its offer of glass equipment, Endo Glassware includes this series of Soxhlet extractors of different capacities and sizes made of

chemical and thermal resistant borosilicate glass.  
3| Each equipment is made up of a flat bottom flask, an extractor body and a Dimroth condenser

▶ 50 mL



Code	Description
PBC001	Soxhlet extractor, 50 mL

Compuesto por:

Code	Description	Ground joint
PBN026	Dimroth condenser	29/32
PBR001	Extractor body	19/26 y 29/32
PBK003	Flat bottom flask, 100 mL	19/26

▶ 125 mL



Code	Description
PBC002	Soxhlet extractor, 125 mL

Compuesto por:

Code	Description	Ground joint
PBN027	Dimroth condenser	45/40
PBR002	Extractor body	29/32 y 45/40
PBK008	Flat bottom flask, 250 mL	29/32

▶ 250 mL



Code	Description
PBC003	Extractor Soxhlet, 250 mL

Compuesto por:

Code	Description	Ground joint
PBN028	Dimroth condenser	55/44
PBR004	Extractor body	29/32 y 55/44
PBK010	Flat bottom flask, 500 mL	29/32