


## COMECTA Tubes with screw cap



### 1. 10 ml round bottom tubes Ø16 x 100mm.

Polystyrene, screenprinted and with screw cap.  
Packed 50 units.

Part No. **5004000** Non-sterile.

Part No. **5004001** Sterile.

### 2. 10 ml round bottom tubes Ø16 x 100mm.

Polystyrene and with screw cap.

Part No. **5004002** Non-sterile. Packed 1000 units.

Part No. **5004003** Sterile. Packed 200 units.

### 3. 13.5 ml conical bottom tubes. Ø16 x 110mm.

Polystyrene and with screw cap.

Rcf: 8,000 g.

Part No. **5004004** Non-sterile. Packed 1000 units.

Part No. **5004005** Sterile. Packed 200 units.

### 4. 15 ml round bottom tubes Ø16 x 120mm.

Polystyrene and with screw cap.

Part No. **5004006** Non-sterile. Packed 750 units.

Part No. **5004007** Sterile. Packed 150 units.

### 5. 15 ml conical bottom tubes. Ø17 x 120mm.

Polystyrene, screenprinted and with screw cap.

Rcf: 3,500 g.

Packed 150 units.

Part No. **5004008** Non-sterile.

Part No. **5004009** Sterile.

### 6. 15 ml conical bottom tubes. Ø17 x 120mm.

Polystyrene, screenprinted and with screw cap.

Rcf: 5,000g, Autoclavables.

Packed 150 units.

Part No. **5004010** Non-sterile..

Part No. **5004011** Sterile.

### 7. 50 ml flat bottom conical tubes. Ø30 x 115mm.

Polystyrene, screenprinted and with screw cap.

Autoclavables.

Packed 100 units.

Part no. **5004012** Non-sterile.

Part no. **5004013** Sterile.

### 8. 50 ml conical bottom tubes. Ø30 x 115mm.

Polystyrene, screenprinted and with screw cap.

Rcf: 5,000g, Autoclavables.

Packed 100 units.

Part no. **5004014** Non-sterile.

Part no. **5004015** Sterile.

### 9. 50 ml conical bottom tubes. Ø30 x 115mm.

**For photosensitive materials.**

Polystyrene, screenprinted and with screw cap.

Rcf: 5,000g, Autoclavables at 121 °C.

Resistant to temperatures -80°C.

Packed 25 units.

Part no. **5004016** Sterile..

### 10. Racks for cylindrical and conical tubes

| Part no.       | Made of        | Tubes number | For tubes up to Ø |
|----------------|----------------|--------------|-------------------|
| <b>1001715</b> | inox. AISI 304 | 15 (5 x 3)   | 17 mm             |
| <b>1025524</b> | inox. AISI 304 | 24 (8 x 3)   | 17 mm             |
| <b>1001716</b> | inox. AISI 304 | 8 (4 x 2)    | 30 mm             |



## COMECTA Petri dishes In polystyrene.

### 1. Plates without compartments. Ø 90 mm. 3 winds

Part No. **5004017** Aseptic. Packed 20 units.

Part No. **5004018** Sterile. Packed 20 units.

### 2. Two compartment plates. Ø 90 mm. 3 winds.

Part No. **5004019** Aseptic. Packed 20 units.

Part No. **5004020** Sterile. Packed 20 units.

### 3. Three compartment plates. Ø 90 mm. 3 winds.

Part No. **5004021** Aseptic. Packed 20 units.

Part No. **5004022** Sterile. Packed 20 units.

### 4. Four compartments plates. Ø 90 mm. 3 winds.

Part No. **5004023** Aseptic. Packed 20 units.

Part No. **5004024** Sterile. Packed 20 units.

### 5. Contact Petri dishes crosslinked for counting.

Part No. **5004025** Sterile. Ø 55 mm. Packed 10 units.

Part No. **5004026** Sterile. Ø 90 mm. Packed 20 units.

### 6. Square Petri dishes. 120x120 mm. 3 winds.

Part No. **5004027** Sterile. Packed 10 units.

### 7. Transport and storage support for Petri dishes of Ø 90 mm.

Capacity: 48 plates.

Part No. **5004028**

