



PBG - Distilling flasks



- 1 | Distilling flasks are used in those reactions requiring evaporation and a later condensation of liquids.
- 2 | Because of this, Endo Glassware distilling flasks are made of thermal shock resistant borosilicate glass and present a spherical base for a uniform heating of the samples to be distilled
- 3 | They are also provided with a long and narrow neck and a side arm to be connected to a condenser, thus, allowing the condensation of the previously evaporated liquid.

Code	Capacity	Ø Neck	Height	Ø Arm	Arm length
PBG001	100 ml	15 mm	174 mm	7 mm	100 mm
PBG002	250 ml	19 mm	216 mm	8 mm	120 mm
PBG003	500 ml	24 mm	275 mm	8 mm	150 mm
PBG004	1000 ml	28 mm	320 mm	9 mm	180 mm