





### ► Polypropylene spatula

With and end in a spoon form and the other end like a spatula. They have a big utility in any laboratory to prepare mixtures or take powder products. They are made of polypropylene; they are very soft and have a big resistance to chemical products and solvents at room temperature. They can be sterilized in

autoclave. They are available in red and blue and with different sizes:

- Spoon rough dimensions (AxB): 20x30 mm
- Spatula rough dimensions (CxD): 14x35 mm



Code	Colour	Length (L)
FKB030	Blue	16 cm
FKB033	Red	16 cm
FKB031	Blue	18 cm
FKB034	Red	18 cm
FKB032	Blue	20 cm
FKB035	Red	20 cm