





► Polypropylene racks




CODE BJB001

- 1| Autoclavable, rectangular shaped and available in different colours and models to accommodate tubes up to 30 mm in diameter.
- 2| Are easily stackable and disassembled to save space and also ideal to hold samples on water baths since they will not float.
- 3| Holes are labelled with alphanumeric marks.
- 4| Approximate dimensions: 150x247x70 mm.

Code	Holes	Diameter	Alignment	Colour
BJB001	90	13 mm	6x15	White
BJB002	90	13 mm	6x15	Blue
BJB003	90	13 mm	6x15	Red
BJB004	90	13 mm	6x15	Yellow
BJB005	60	16 mm	5x12	White
BJB006	60	16 mm	5x12	Blue
BJB007	60	16 mm	5x12	Red
BJB008	60	16 mm	5x12	Yellow
BJB009	40	20 mm	4x10	White
BJB010	40	20 mm	4x10	Blue
BJB011	40	20 mm	4x10	Red
BJB012	40	20 mm	4x10	Yellow
BJB013	24	25 mm	3x8	White
BJB014	24	25 mm	3x8	Blue
BJB015	24	25 mm	3x8	Red
BJB016	24	25 mm	3x8	Yellow
BJB017	21	30 mm	3x7	White
BJB018	21	30 mm	3x7	Blue
BJB019	21	30 mm	3x7	Red
BJB020	21	30 mm	3x7	Yellow

► Multiposition rack



- 1| Prepared to hold tubes of different sizes depending on the position and made of polypropylene, autoclavable.
- 2| It presents a solid and very resistant rectangular structure, easily stackable and ideal to keep samples in a bath since it does not float.
- 3| It is also provided with lateral hooks to join several racks.

Code	Capacity	Alignment	Dimensions
BJB026	32 tubes of 11 mm Ø	4x8	175x95x52 mm
	12 tubes of 16 mm Ø	2x6	
	4 tubes of 29 mm Ø	1x4	




racks and supports [B] GENERAL USE

▶ Plexiglas test tubes rack






Code	Capacity	Ø tubes	Alignment	Dimensions
BJB027	12 tubos	14 mm	2x6	105x55x50 mm
BJB028	24 tubos	14 mm	2x12	205x55x50 mm
BJB029	48 tubos	14 mm	4x12	220x90x50 mm
BJB030	6 tubos	22 mm	2x3	140x95x100 mm
BJB031	12 tubos	22 mm	2x6	175x95x100 mm
BJB032	6 tubos	32 mm	2x3	140x95x75 mm
BJB033	12 tubos	32 mm	2x6	225x95x75 mm
BJB034	24 tubos	32 mm	2x123	-

▶ Multi support racks



Code BJB023

1 | Made of PP, with a good resistance to chemical agents and autoclavable, these racks allow supporting tubes in a vertical position. Thanks to their design, they can also be used to dry tubes in inverted position and to support and dry electrophoresis and chromatography plates.

Code	Capacity	Ø tubes	Colour	Dimensions
BJB021	40 tubos	17-20 mm	 Red	23x10x7 cm
BJB022	50 tubos	14-17 mm	 White	18x10x7 cm
BJB023	80 tubos	10-13 mm	 White	23x8.5x6.5 cm


▶ Gradilla desmontable



1 | Fabricada en polipropileno, puede ser desmontada fácilmente para su limpieza o esterilización en autoclave.

Code	Capacity	Ø tubes	Colour	Dimensions
BJB035	31 tubos	13 mm	 White	27x6.5x12 cm

▶ Anodized aluminium racks



1 | Its manufacturing material provides a great durability and resistance to chemical and corrosive agents. Available in three different sizes, are "Z" shaped and its wide base ensures a high stability and allows one-hand operation. Are easily stackable.

Code	Nº holes	Diameter	Dim. (LxAxH)
BJD001	50	12 mm	170x995x55 mm
BJD002	50	14,5 mm	200x110x55 mm
BJD003	50	17,0 mm	225x120x68 mm

RACKS AND SUPPORTS [B] GENERAL USE

Note: the manufacturer reserves the right to modify the characteristics and photographs of the products



racks and supports [B] GENERAL USE

▶ Stainless steel racks (wire racks)



Autoclavable and very resistant, their special design makes them ideal for conical bottom tubes. All racks present square holes and reinforced edges.

Code	N° holes	Diameter	Dimensions	Alignment
BJD004	12	13 mm	80x55x85 mm	3x4
BJD005	24	13 mm	150x60x85 mm	3x8
BJD006	50	13 mm	175x90x85 mm	5x10
BJD007	100	13 mm	220x220x95 mm	10x10
BJD008	12	16 mm	135x85x50 mm	2x6
BJD009	24	16 mm	175x70x85 mm	3x8
BJD010	50	16 mm	210x105x85 mm	5x10
BJD011	72	16 mm	265x135x105 mm	6x12
BJD012	100	16 mm	240x240x105 mm	10x10
BJD013	24	18 mm	190x75x90 mm	3x8
BJD014	50	18 mm	230x120x90 mm	5x10
BJD015	72	18 mm	285x145x100 mm	6x12
BJD016	100	18 mm	240x240x105 mm	10x10
BJD017	12	25 mm	190x68x105 mm	2x6
BJD018	24	25 mm	250x100x105 mm	3x8
BJD019	50	25 mm	300x155x105 mm	5x10
BJD020	72	25 mm	105x365x180 mm	6x12

▶ Racks for vials




Code	N° holes	Diameter	Alignment
BJD036	50	12 mm	10x5
BJD037	24	12 mm	6x4

▶ Stainless steel racks, for type Falcon tubes



Code	N° holes	Diameter	Alignment
BJD034	24	18 mm	6x4
BJD035	24	30 mm	6x4

▶ Stainless steel racks



Code	N° holes	Diameter	Alignment
BJD022	24	10 mm	6x4
BJD023	50	12 mm	10x5
BJD026	24	16 mm	6x4
BJD027	50	16 mm	10x5
BJD028	24	18 mm	6x4
BJD029	50	18 mm	10x5
BJD030	24	20 mm	6x4
BJD031	50	20 mm	10x5
BJD032	24	22 mm	6x4
BJD033	24	30 mm	6x4

RACKS AND SUPPORTS [B] GENERAL USE

Note: the manufacturer reserves the right to modify the characteristics and photographs of the products