



Lactometers are glass instruments used to measure milk density and thus, determine if it has been mixed with water or cream has been partially removed.

Glass body with scale in lacto-densimetric degrees (15° equals 1015g/mL). Calibrated at 20°C and optionally available with thermometer for deviation correction.

Individually packaged in a case. Nahita lactometers are of the Quevenne type whose graduated stem presents a scale between 15 and 40 that correspond to density thousandths over the unit (the number 32 of the lactometer indicates a density of 1032 Kg/m<sup>3</sup>).



Code	Range	Divisions	Length	Diameter	Thermometer
KLG014	15-40	1	215 mm	19 mm	No
KLG015	15-40	1	290 mm	23 mm	Yes (0-40 °C)