



- 1 | Thanks to their weight they are used to stabilize recipients (bottles, flasks, etc.) avoiding them to float when submerged in fluids or water baths
- 2 | Made of steel covered with PVC they are available in different sizes

Code	Ø Inner	Weight (approx.)
FZB001	48 mm	270 g
FZB002	56 mm	520 g
FZB003	61 mm	780 g
FZB004	74 mm	960 g

## STABILIZING RINGS

### [F] LABORATORY ASSEMBLIES AND TOOLS

Note: the manufacturer reserves the right to modify the characteristics and photographs of the products