



All of them are filled with coloured spirit avoiding mercury contamination problems in case of breakage.

They are ideal to measure temperature of fridges, freezers or refrigerating chambers thanks to its temperature range reaching very low temperatures.



- 1 | External brown graduation
- 2 | Red spirit filled

Code	Description	Divi.	Range
KJP001	Solid-stem	1 °C	-40/40 °C



- 1 | Plastic body
- 2 | Length, 210 mm
- 3 | Blue spirit filled

Code	Description	Divi.	Range
KJP005	With hook	1 °C	-40/45°C