

1 | Nahita offers this group of liquid hydrometers with specific graduation and scale to measure the density of certain liquids (salts, syrup, milk, honey...).



Code	Model	Range	Divisions	Length	Diamètre
KLN001	Acid hydrometer	0-45 °Be	1	250 mm	11 mm
KLN002	Concentrated acid hyd	0-70 °Be	1	285 mm	11 mm
KLN003	Syrup hydrometer	0-50 °Be	1	250 mm	11 mm
KLN004	Milk hydrometer	0-25 °Be 1000-1040 Kg/m3	1 2	195 mm	19 mm
KLN005	Familiar milk hyd.	0-100 % agua	1/4	150 mm	12 mm
KLN006	Lye hydrometer	0-45 °Be	1/2	255 mm	12 mm
KLN007	Salt hydrometer	0-35 °Be	1/2	210 mm	13 mm
KLN008	Cider hydrometer	0-7 % alcohol	1/10	240 mm	18 mm
KLN009	Tannin hydrometer	0-60 °Be	1	260 mm	20 mm