



1 | Nahita offers this group of liquid hydrometers with specific graduation and scale to measure the density of certain liquids (salts, syrup, milk, honey...).



Code	Model	Range	Length	Diamètre
KLN003	Syrup hydrometer	0-50 °Be	250 mm	11 mm
KLN006	Lye hydrometer	0-45 °Be	255 mm	12 mm
KLN011	Heavy liquid	0- 70° Be/1,000-2,000 g/mL	290 mm	15 mm
KLN012	Light liquid	10- 70° Be/0,700-1,000 g/mL	305 mm	15 mm
KLN013	-	0- 50° Be/1,000-1,500 g/mL	305 mm	18 mm
KLN014	-	50 70° Be/1,500-2,000 g/mL	305 mm	15 mm
KLN015	-	70- 70° Be/0,700-2,000 g/mL	150 mm	10 mm