

► Clevenger apparatus, “heavy” essential oil determination



- 1| Essential oil determination for oil heavier than water
- 2| Heat resistant low expansion 3.3 borosilicate glass
- 3| Round bottom flask 1000 mL, 24/29
- 4| Liebig condenser 300 mm, 24/29

Code	Description	Capacity
PBB017	Clevenger apparatus “heavy” essential oil determination	1000 mL

► Clevenger apparatus, “light” essential oil determination



- 1| Essential oil determination for oil lighter than water
- 2| Heat resistant low expansion 3.3 borosilicate glass
- 3| Round bottom flask 1000 mL, 24/29
- 4| Liebig condenser 300 mm, 24/29

Code	Description	Capacity
PBB018	Clevenger apparatus “light” essential oil determination	1000 mL