



POLYPROPYLENE TOPS FOR CULTURE TUBES “BACTI-CAPS”

General features:
 Adaptable to the tubes by three flexible internal flaps that hold onto the edge of the tube, not a hermetic seal, for minimal sample evaporation.
 Highly resistant to culture media, acids, alkalis, alcohol and ether.
 Easy to clean, can be sterilized at 121 °C (1 atm, vapour pressure).
 Not suitable for dry sterilization.
 Come in 5 colours:
 YELLOW, BLUE, WHITE, RED and GREEN.

Part No.	For tube Ø 12 - 13 mm	Part No.	For tube Ø 15 - 16 mm	Part No.	For tube Ø 18 - 19 mm	Part No.	For tube Ø 24 - 25 mm	Part No.	For tube Ø 37 - 38 mm
5726013	Yellow Colour	5726016	Yellow Colour	5726019	Yellow Colour	5726025	Yellow Colour	5726238	White Colour
5726113	Blue Colour	5726116	Blue Colour	5726119	Blue Colour	5726125	Blue Colour	Boxes of 25 units	
5726213	White Colour	5726216	White Colour	5726219	White Colour	5726225	White Colour		
5726413	Red Colour	5726416	Red Colour	5726419	Red Colour	5726425	Red Colour		
5726513	Green Colour	5726516	Green Colour	5726519	Green Colour	5726525	Green Colour		
Packet of 100 units.		Packet of 100 units.		Packet of 100 units.		Packet of 100 units.			



WEIGHING BOATS “ACCU-WEIGH”

Anti-static, Suitable for weighing with reproducible tare weight.
 Free from contaminants, Wide range of sizes.



General features:
 Made from high impact polystyrene with a flat base, round edge, smooth surface, contamination free and hydrophilic.
 Flexible and easy to manipulate, suitable for liquids and solids.
 Supplied in two forms:
 - Diamond, black or white colours.
 - Rectangular, white colour.

Part No.	Shape	Capacity ml	External dimensions mm			Weight g	Colour	Quantity
			Height	Width	Depth			
5726004	Diamond	5	3.8	31	54	0.25	White	100 units
5726005	Diamond	5	3.8	31	54	0.25	Black	100 units
5726006	Diamond	30	13.5	55	85	1.00	White	100 units
5726007	Diamond	30	13.5	55	85	1.00	Black	100 units
5726008	Diamond	100	20	92	132	2.6	White	100 units
5726009	Diamond	100	20	92	132	2.6	Black	100 units
5726010	Rectangular	7	8.5	44	44	0.65	White	100 units
5726011	Rectangular	100	24	80	80	2.2	White	100 units
5726012	Rectangular	250	24	134	134	6.1	White	100 units



ALUMINIUM DISHES

Pressed with crenulated edge.
 No static.
 Withstands temperatures up to 450 °C.
 Can be used as a lid by mounting one on top of another.
 Ideal for sample testing of solids, humidity, mixtures of emulsions and moisture determinations by weight loss.
 Suitable use as a weighing boat, cooling or evaporation of samples.
 Single use and disposable.

Part No.	Thickness mm	Height / Ø (usable) mm		Quantity
5726000	0,2	5	42	packet of 112 units
5726001	0,2	10	57	packet of 110 units