



- 1 | Nahita offers this series of mortars made of porcelain and ideal for grinding, crushing and mixing substances of low hardness.
- 2 | With a wide and stable base and a spout for an easy pouring.

- 3 | Each mortar is supplied together with its porcelain pestle with the handle glazed and the grinding surface unglazed.
- 4 | Mortar pestles are also available individually.

▶ Porcelain mortars with pestle, unglazed



1 | Glazed outside except for the base and unglazed inside.

Code	Capacity	Diameter	Height	Pestle length	Ø Ball
BBC001	25 ml	60 mm	35 mm	65 mm	19 mm
BBC002	50 ml	80 mm	45 mm	90 mm	26 mm
BBC003	100 ml	100 mm	45 mm	110 mm	27 mm
BBC004	300 ml	130 mm	62 mm	125 mm	36 mm
BBC005	500 ml	160 mm	74 mm	145 mm	40 mm
BBC006	1000 ml	216 mm	90 mm	215 mm	46 mm
BBC007	2000 ml	254 mm	110 mm	245 mm	53 mm

▶ Porcelain mortars with pestle, glazed



1 | Glazed outside except for the base and unglazed inside.

Code	Capacity	Diameter	Height	Pestle length	Ø Ball
BBB001	100 ml	100 mm	45 mm	110 mm	27 mm
BBB002	300 ml	130 mm	62 mm	125 mm	36 mm
BBB003	500 ml	150 mm	74 mm	145 mm	40 mm
BBB004	1000 ml	210 mm	90 mm	195 mm	46 mm

▶ Pestles for porcelain mortars



Code	Length	Ø Handle	Ø Ball	For mortars
BBN001	65 mm	11 mm	19 mm	60 mm
BBN002	90 mm	15 mm	26 mm	80 mm
BBN003	110 mm	15 mm	27 mm	100 mm
BBN004	125 mm	18 mm	36 mm	130 mm
BBN005	145 mm	20 mm	40 mm	160 mm
BBN006	200 mm	27 mm	45 mm	216 mm