



**TALL BEAKERS AND JUGS**

Mirror finish stamped beakers and jugs with a curved lip. Easily cleaned and resistant to large numbers of chemical agents.

1 Standard	2 With Spout	3 With handle and spout		
Part No.	Part No.	Part No.	Ø x Height	Volume
8002902	8012902	8092902	72 x 84 mm.	250 ml.
8002905	8012905	8092905	86 x 129 mm.	500 ml.
8029010	8129010	8929010	102 x 150 mm.	1000 ml.
8029020	8129020	8929020	132 x 184 mm.	2000 ml.
8029030	8129030	8929030	160 x 180 mm.	3000 ml.
8029040	8129040	8929040	180 x 200 mm.	4000 ml.
8029050	8129050	8929050	200 x 220 mm.	5000 ml.



**LIDS FOR BEAKERS AND JUGS**

Part No.	Beaker size	Part No.	Beaker size
8000250	250 ml.	8003000	3000 ml.
8000500	500 ml.	8004000	4000 ml.
8001000	1000 ml.	8005000	5000 ml.
8002000	2000 ml.		



**BUCKET WITH HANDLE**

Buckets. Graduated from 1 to 12 litres. 320 Ø x 300 mm high. Volume 12 litres.

Part No. **8000252**

**Bucket lids**

Part No. **1002521**




**MEDICINE CUPS**

Part No.

**8085530** Graduated. Capacity 30 ml.

**8085531** Graduated. Capacity 50 ml.

**8085532** With handle, no graduations. Capacity 300 ml.



**SERILISATION JAR**

100 Ø x 155 mm high. Volume 1000 ml.

Part No. **8080904**



**ROUND BASIN**

370 Ø x 145 mm high. Volume 9000 ml.

Part No. **8000140**



**SUPPORT FOR ROUND BASIN**

Portable with rubber wheels and brake.

Basin Part No. 8000140.

Part No. **6010140**



**SUPPORT FOR BUCKETS**

Portable with rubber wheels and brake.

For bucket Part No. 8000252.

Part No. **6010252**